AVANOS*

BALLARD* MINI-BAL SAMPLING CATHETER

Instructions for Use


Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5


Fig. A

Instructions for Use

Rx Only: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Description

Mini-BAL is a catheter within a catheter. The 12 Fr inner catheter (Fig. 2a) is contained within a 16 Fr outer sheath (Fig. 2b). The inner catheter has a rounded mushroom shaped tip (Fig. 2c) that seals and covers the end of the outer catheter. This design permits the advancement of the inner catheter into the distal bronchial segments reducing exposure to upper airway microorganisms. The mushroom tip is radiopaque for radiographic confirmation of catheter location. There is a port for the delivery of low flow oxygen (Fig. 2d) when using this device on a non-intubated patient. A large diameter stopcock allows connection to suction source and specimen trap (Fig. 2e).

Designed to be compatible with BALLARD* Multi-Access Port (MAP) Closed Suction System (CSS).

Indication

Mini-BAL is used in the diagnosis of lung disease by allowing collection of bronchoalveolar lavage (BAL) specimens from deep within the lung (Fig. 5). The use of a bronchoscope is not necessary. This catheter may be used in both intubated and non-intubated adult patients.

Contraindication

There are no absolute contraindications for BAL (Bronchoalveolar Lavage). Relative contraindications to this procedure include:

- Lack of patient cooperation
- FVC (forced vital capacity) of less than one liter
- Asthma with moderate airway obstruction
- Hypercapnia
- Hypoxemia that cannot be corrected to a saturation greater than 90%
- Serious cardiac arrhythmia
- Myocardial infarction within the previous 6 weeks
- Uncorrected bleeding susceptibility
- Hemodynamic instability
- Hemoptysis of unknown origin

⚠Warning

Non-bronchoscopic bronchial alveolar lavage (BAL) may be expected to have complications similar to those seen in bronchoscopic bal procedures, which may include: pneumothorax, pneumonitis, post procedure fever, bronchospasm and bleeding.

Do not reuse, reprocess, or resterilize this medical device. Reuse, reprocessing, or resterilization may 1) adversely affect the known biocompatibility characteristics of the device, 2) compromise the structural integrity of the device, 3) lead to the device not performing as intended, or 4) create a risk of contamination and cause the transmission of infectious diseases resulting in patient injury, illness, or death.

⚠Caution

This medical device contains DEHP (Di(2-ethylhexyl) phthalate) that is presently classified in the European Union as a presumed human reproductive toxicant based on data from animal studies. There is no conclusive scientific evidence that exposure to DEHP contained in medical devices has caused harmful effects in humans. A risk assessment, which took into account DEHP exposure of all indicated patient populations, including those who are potentially at increased risk, was performed for this device and the conclusion is that the device is safe when used as directed.

Preparation

Prepare patient for procedure based on hospital protocol.

Suggested equipment:

- Mini-BAL (Ref 142)
- Sterile gloves
- Dressing pack
- Sterile saline
- 20/50 ml syringes
- Sputum trap
- Suction tubing

Adult Intubated Patients

1. Remove Mini-BAL from package.
2. Remove protective clip from tip of Mini-BAL catheter.
- 3A. Set-Up Protocol (For use directly with artificial airway)
 - a. Insert the Mini-BAL catheter into and through the access port elbow. (Fig. 2f)
 - b. Disconnect the ventilator circuit from the endotracheal or tracheostomy tube.
 - c. Attach access port elbow to the ventilator circuit. (Fig. 3)
 - d. Attach access port elbow to the endotracheal or tracheostomy tube while guiding the Mini-BAL catheter into the lumen (Fig. 4).
 - e. Position the curve of directional tip catheter for right or left lung (oxygen port (Fig. 2d) on same side).
 - f. Attach suction adapter (Fig. 1a) and 20/50 ml syringe with saline to 3-way stopcock. (Fig. 2e)
 - g. Attach sputum trap and suction apparatus to Mini-BAL and proceed with flushing and sampling protocols.
- 3B. Set-Up Protocol (For Use with BALLARD* MAP CSS)
 - a. Depress plunger button and rotate BALLARD* MAP Closed Suction Catheter (CSC) so the "Alternate Therapy Port" is in-line with the airway. (Fig. A)

⚠Caution: Do not rotate with instrument inserted through the "Alternate Therapy Port".
 - b. Open cap of Seal Cassette (Fig. Aa) on BALLARD* MAP CSC.
 - c. Remove Mini-BAL adapter (Fig. 1d) from package and attach onto Seal Cassette. (Fig. Aa)
 - d. Insert Mini-BAL through the Mini-BAL adapter into Seal Cassette. (Fig. A)
 - e. Follow steps 3e – g above.
4. Flushing Protocol
 - a. Advance Mini-BAL 5 cm beyond the end of the endotracheal or tracheostomy tube.
 - b. Clear tip by flushing with 2 ml saline.
5. Sampling Protocol
 - a. Hold outer sheath in place and advance the inner catheter until resistance is met. (Fig. 5)
(Inner catheter should now be in a wedge position.)

- b. Lock outer sheath in position by sliding blue slidelock (Fig. 2g) on Elbow to locked position.
- c. Perform BAL according to protocol.
6. After obtaining the final sample, pull the inner catheter back until the black strip is visible behind the control hub seal.
- 7A. Detachment Protocol (For use directly with artificial airway)
 - a. Remove Mini-BAL and elbow from the circuit and reconnect the ventilator.
- 7B. Detachment Protocol (For Use with BALLARD® MAP CSS)
 - a. Open the blue slidelock of Mini-BAL adapter and retract Mini-BAL completely.
 - b. Depress plunger button and rotate BALLARD® MAP Closed Suction Manifold to the "Off" or "Suction" position.
 - c. Remove the Seal Cassette and Mini-BAL Adapter by rotating counter clockwise.
 - d. Replace with new Seal Cassette (Fig. 1c) found in Mini-BAL package after each use or instrument insertion.
10. Low flow oxygen, up to 5 LPM, may be administered by oxygen line to the oxygen port (Fig. 2d). Oxygen should not be initiated until catheter placement within the trachea has been confirmed.
11. To reach the wedge position, hold the outer sheath in place and advance the inner catheter until resistance is met. (Fig. 5) (Achievement of a proper wedge position will ensure adequate lavage solution returns. If lavage return is inadequate, slight withdrawal of the inner catheter may be necessary. If the patient coughs fluid from the airway, the wedge is inadequate and the catheter needs to be advanced).
12. Perform BAL according to protocol.
13. After obtaining the final sample, pull the inner catheter back until the black strip is visible behind the control hub seal. Remove the Mini-BAL from the patient airway.

Adult Non-Intubated Patients


1. Prepare the upper airway according to standard anesthetic protocol for transnasal bronchoscopy.
2. Attach suction adapter (Fig. 1a) and 20/50 ml syringe with saline to 3-way stopcock. (Fig. 2e)
3. Attach sputum trap and suction apparatus to Mini-BAL.
4. Guide Mini-BAL through the nostril until the tip of the catheter can be seen in the midline of the pharynx. The directional tip should be on the midsagittal plane.
5. Place the head in a neutral position or with the chin tipped slightly toward the chest.
6. Instruct the patient to take a slow deep breath, advance the catheter during this inspiration.
7. Position the curve of directional tip catheter for right or left lung (oxygen port (Fig. 2d) on same side). To guide the directional tip, rotate the oxygen port right or left toward the desired lung.
8. Continue to advance the Mini-BAL until the catheter tip is approximately 15 cm beyond the vocal cords.
9. Flush the tip by instilling 2 ml saline. Check for up and down movement of fluid within the catheter with respiration to confirm airway placement.

⚠ Warning

Oxygen flow through the Mini-BAL device should not be initiated until catheter placement within the trachea or bronchi has been confirmed by evidence of expiratory vapor inside the catheter, to and fro movement of flush solution, or x-ray confirmation.

⚠ Caution

Low flow oxygen, up to 5 LPM, prior to, and during the procedure may be desirable. Oxygen should be administered by nasal cannula prior to the procedure and then delivered through the Mini-BAL oxygen port during the procedure. Pulse oximetry should be used to monitor oxygen saturation. Appropriate monitoring and emergency resuscitation equipment should be readily available.

 Diameter	 Length	 Single Use Only	 Sterilized using ethylene oxide
 Do not use if package is damaged	 Do not resterilize	 Not made with natural rubber latex	Contains  PHT DEHP
Rx Only	 Caution	 Consult instructions for use	

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*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.
© 2018 AVNS. All rights reserved. 2019-07-01
Distributed in the USA by Avanos Medical Sales, LLC,
Avanos Medical, Inc., 5405 Windward Parkway, Alpharetta, GA 30004 USA.
Sponsored in Australia by Avanos Medical Australia Pty Ltd,