**REF** 002600, 002602, 002611, 002600-14U, FM2625, SFT2625, 1320

**EN**

**AirLife® Cushion Nasal Cannula**

**Device Description:** An over-the-ear style cannula with nasal prongs used to deliver supplemental oxygen to the patient's nostrils (Figure A). Includes crushed resistant tubing with either a trumpet or U/Connect-IT Adapter. Device is non-sterile, disposable and for single patient use.

Fig. A

**Intended Purpose:** Intended for oxygen to be delivered to the patient at a prescribed concentration and flow rate. Oxygen tubing is connected to an oxygen source on one end and a cannula on the other end to deliver oxygen to the patient.

**Indications for Use:** For patients prescribed supplemental oxygen by a nasal cannula. To treat hypoxia, relieve shortness of breath and pre-oxygenate for surgical procedures in spontaneously breathing patients.

**Environment:** Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

**MR Safety:** MR Safe

**Patient Target Group:** The intended population includes neonatal, pediatric and adult patients.

Spontaneously breathing individuals with a variety of breathing conditions requiring supplemental oxygen.

**Expected Clinical Benefits:**

- Nasal cannulas are a conduit that deliver supplemental oxygen to a patient's nasopharynx to treat or prevent hypoxemia in acute and chronic values, vital signs, and/or work of breathing.
- Supplemental oxygen helps reduce COPD complications by stabilizing pulmonary hypertension, reducing secondary polycythemia, and decreasing arrhythmias.
- Supplemental oxygen can decrease anxiety, improve sleep, improve mental alertness, and/or improve health-related quality of life in people with chronic respiratory disease. Using oxygen during exercise can improve endurance, heighten performance, and/or decrease the sensation of breathlessness. This helps patients complete physical activities, improve ability to walk and/or perform other physical activities for an extended period.
- The crush resistant tubing helps reduce kinking and allows for greater than 75% of initial flow to be delivered should the tubing become kinked.

**Contraindications:** None known.

**Warnings:**

- Patient may become hypoxic if oxygen flow is interrupted.
- If oxygen is in use, do not use near flame or heat source.
- Position tubing to prevent patient from becoming tangled in tubing or strangulation.
- Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
- To reduce the risk of misconnections and patient injury, always trace tubing and lines back to their point of origin to verify that correct connections are made.
- High levels of oxygen may cause oxygen-induced hyperventilation, oxygen-induced hypercapnia, or oxygen toxicity. In premature infants, it may also cause retrolental fibroplasia.
- Do not reuse on more than one patient as this may increase risk of contamination and infection.

**Cautions:**

- Consult physician if the patient develops an infection, skin irritation, or material sensitization.
- Flow rates above 6 LPM may cause increased resistance or back pressure in the tubing.
- To prevent interruption of therapy or injury, do not let children play with the oxygen tubing.
- To prevent excess moisture accumulation in the tubing, do not use with heated humidification.
- To prevent equipment damage and/or interruption of therapy, do not let pets play with the oxygen tubing.
- Do not place tubing under rugs or other objects as it may obstruct flow.
- To prevent from becoming tangled in tubing, position the tubing in front of you before turning to sit or stand.

**Residual Risks:** Refer to Warnings and Cautions.

**Necessary User Qualifications:** This device is restricted to sale by or on the order of a physician. The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings and cautions as stated in the IFU. Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

**Instructions for Use:**

- Attach the oxygen tubing to oxygen source.
- Adjust the flow control knob to the prescribed liter flow.
- Check connections for leaks and verify flow from the prongs.
- Position the nasal cannula with the prongs facing upward and curved toward the face. Insert the two prongs into

**EN**

- the nostrils.
- Wrap the headset loop up and over the ears and tuck the loop under the chin, or adjust and secure the headset loop behind the head. If applicable, position foam covers over top of ears.
- Squeeze the sides of the bolo and glide the bolo up under the chin, or until snug around the head.
- Leave enough space to fit at least the width of two fingers between the bolo and chin.
- Discard and replace cannula if soiled or damaged. Do not sterilize.

**Compatibility to Other Devices:** AirLife tubing connectors are designed to work with gas outlets that comply to ISO 5359:2014. This includes, but is not limited to, stand-alone and portable oxygen concentrators, oxygen tanks, oxygen flowmeters, and bubble humidifiers. The female threaded connector, also referred to as U/Connect-it, is designed to be used with a male DISS 1240 oxygen fitting, 5 mm - 7 mm tapered oxygen fitting, or a ¼ inch (6.35 mm) tapered oxygen fitting, e.g., nipple & nut adapter (Figure B).

- Place the threaded connector between your thumb and forefinger below the threaded section (Figure C).
- Position the connector's end with either the threaded male fitting or the tapered fitting.
  - For attachment to male fitting, firmly push the connector onto the oxygen fitting and turn clockwise until a secure fit is achieved. (Figure C).
  - For attachment to a tapered fitting, with slight pressure, slide the connector onto the gas fitting and turn clockwise until a secure fit is achieved.


Fig. B


Fig. C

- Gently tug on tubing to ensure secure fit.

**Safe Disposal:** Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

**Storage:** Recommended normal storage temperature 0°C/+32°F to +50°C/+122°F.

**Incident Reporting:** Notice to the user and/or patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. If your physical symptoms worsen or you experience a sudden change in your condition (e.g., increased shortness of breath, fever, dizziness) call your doctor.

**Performance Characteristics:**

	002611	002602	002600
Patient	Neonatal	Pediatric	Adult
Prong Size	Curved	Curved	Curved
Tubing Length	7 ft (2.1m)	7 ft (2.1m)	7 ft (2.1m)
Tubing End Connector	Trumpet	Trumpet	Trumpet
Flow Rate	3 l/min	6 l/min	6 l/min
Special Feature	N/A	N/A	N/A
	002600-14U	FM-26225	SFT2625
Patient	Adult	Adult	Adult
Prong Size	Curved	Curved	Curved
Tubing Length	14 ft (4.3m)	25 ft (7.6m)	25 ft (7.6m)
Tubing End Connector	U-Connect IT	Trumpet	Trumpet
Flow Rate	6 l/min	6 l/min	6 l/min
Special Feature	N/A	Foam Ear Wraps	Soft Headset Tubing
	001320	001225	
Patient	Adult	Adult	
Prong Size	Curved, Flared Tip	Curved, Non-Flared Tip	
Tubing Length	7 ft (2.1m)	7 ft (2.1m)	
Tubing End Connector	Trumpet	Trumpet	
Flow Rate	6 l/min	l/min	
Special Feature	N/A	N/A	

**EN**

**REF** Catalogue Number  Manufacturer

**LOT** Lot Number  Expiration Date

**X** Quantity  Authorized Representative in the European Community

 Not Made with Natural Rubber Latex  Medical Device

 Not Made with DEHP  Cautions

 Consult Instructions for Use

**Rx ONLY** U.S. Federal Law restricts this device to sales by or on order of a physician

**CLEAN, READY TO USE** Clean, Ready to Use

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© 2019 Vyaire. Vyaire, the Vyaire Logo and AirLife are trademarks or registered trademarks of Vyaire Medical, Inc., or of its affiliates.
AirLife Australia Holding Pty. Ltd.
Vyaire Medical, Inc.