REF 001200, 001201,
001262, 001260

EN Oxygen Mask

Device Description: A medium concentration, under-the-chin, oxygen mask that is placed over a patient's nose or mouth to administer oxygen. Includes a crush-resistant 7 ft (2.1 m) tubing with trumpet end connector. The device is non-sterile, disposable and for single-patient use.


Fig. A

Intended Purpose: Intended for oxygen to be delivered to the patient at a prescribed concentration and/or flow rate. Oxygen tubing is connected to an oxygen source on one end and a mask on the other end to deliver oxygen to the patient.

Indications for Use: For patients requiring supplemental oxygen. To treat hypoxia, relieve shortness of breath and pre-oxygenate for surgical procedures in spontaneously breathing patients.

Environment: Hospitals, sub-acute, medical clinics, pre-hospital-, home, surgical centers, skilled nursing facilities.

MR Safety: MR Safe

Patient Target Group: The target population is adult patients.

Expected Clinical Benefits: The oxygen mask delivers up to 15 LPM of oxygen to reverse or prevent hypoxemia, and provide relief of shortness of breath in oral and nasal breathers.

Contraindications: None known.

Warnings:

- The oxygen connector in mask should not be connected directly to any endotracheal or tracheostomy tube inlet.
- If oxygen is in use, do not use near a flame or heat source.
- Position tubing to avoid strangulation.
- Device contains components that may present choking hazard.
- Patient may become hypoxic if oxygen flow is interrupted.
- Use a flow rate \geq 4 LPM to prevent CO₂ accumulation and re-breathing.
- Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
- To reduce the risk of misconnections and patient injury, always trace tubing from gas source to the medical device before connecting.
- Do not reuse on more than one patient as this may increase risk of contamination and infection.

Cautions:

- Consult physician if patient develops infection, skin irritation or material sensitization.
- Do not place anything on supply tubing that may obstruct flow.

Material(s) Disclosure: Risks and Precautionary Measures related to Phthalates: This instruction pertains to the phthalate symbol marked on the device or its packaging. If this device is used for the treatment of children, or treatment of pregnant or nursing women; please note that the following types of procedures may increase the risk of exposure to phthalates: Exchange transfusion in neonates, total parenteral nutrition in neonates, multiple procedures in sick neonates, hemodialysis in peripubertal males, male fetus and male infant of pregnant women, and lactating women; and massive blood infusion into trauma patients. Although these procedures have the potential for increased risk of exposure, conclusive evidence of human health risks has not been established. As a precautionary measure, to reduce the potential for unnecessary exposures to phthalates, the product must be used in accordance with the instructions for use, and practitioners should refrain from using this product beyond the period of time the product is medically necessary or needed.

Residual Risks: Refer to warnings and cautions.

Necessary User Qualifications: This device is restricted to sale by or on the order of a physician. The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings and cautions as stated in the IFU. Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

Instructions for Use:

- Attach one end of oxygen tubing to mask and the other end to oxygen source.
- Adjust liter flow as prescribed by physician.
- Place mask over patient's face covering both nose and mouth.
- Slip strap over patient's head and adjust for comfort and fit.
- Pinch nose clip to provide a seal around the nose.
- Check for leaks at connections and for proper placement of the mask.

Compatibility to Other Devices: AirLife tubing connectors are designed to work with gas outlets that comply to ISO 5359:2014. This includes, but is not limited to oxygen tanks, oxygen flowmeters, and oxygen blenders. The tubing trumpet end connector is compatible with a 5 mm - 7 mm tapered oxygen fitting or 1/4 inch (6.35 mm) tapered oxygen fitting (Figure B). The style of the tapered fitting will vary depending on the oxygen gas source or oxygen device.

- Place the ribbed or trumpet end connector between your

- EN**
- thumb and forefinger (Figure C).
 - Align the end of the connector to tapered fitting, then push the end connector over the tapered fitting until snug (Figure C).
 - Gently tug on the tubing to ensure a secure fit.


Fig. B


Fig. C

Safe Disposal: Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

Storage: Recommended normal storage temperature 0°C/+32°F to +50°C/+122°F.

Incident Reporting: Notice to the user and/or patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. If your physical symptoms worsen or you experience a sudden change in your condition (e.g., increased shortness of breath, fever, dizziness) call your doctor.

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© 2019 Vyair. Vyair, the Vyair Logo and AirLife are trademarks or registered trademarks of Vyair Medical, Inc., or of its affiliates.
AirLife Australia Holding Pty. Ltd.
Vyair Medical, Inc.