
EN 3-in-1 Oxygen Mask

Device Description: Under-the-chin style, 3 in 1 oxygen mask that is placed over a patient's nose or mouth to administer oxygen. Includes crush resistant oxygen tubing with U-Connect-It end connector. Disposable, single patient use. Mask can be configured as a nonrebreather mask, rebreather mask, or medium concentration mask.

Oxygen Mask Consists of:

1. Elastic Strap
2. Nose Clip
3. Adult Oxygen Mask
4. Oxygen Bag Connector
5. Green Tape
6. Oxygen Bag
7. 7 feet (2.1m) Crush Resistant Oxygen Tubing
8. U/Connect-It Adapter (Threaded Adapter)
9. Oxygen Connector
10. Rebreather Gasket Synthetic
11. Vinyl Tip Adapter
12. Gasket Synthetic Polyisoprene

Intended Purpose: Intended for oxygen to be delivered to the patient at a prescribed concentration and/or flow rate. Oxygen tubing is connected to an oxygen source on one end and a mask on the other end to deliver oxygen to the patient.

Indications for Use: For patients requiring supplemental oxygen. To treat hypoxia, relieve shortness of breath and pre-oxygenate for surgical procedures in spontaneously breathing patients.

Environment: Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

MR Safety: MR Safe

Patient Target Group: The target population is adult and pediatric patients.

Expected Clinical Benefits: The oxygen mask delivers up to 15 LPM of oxygen to reverse or prevent hypoxemia, and provide relief of shortness of breath in oral and nasal breathers.

Contraindications: None known.

Warnings:

- The oxygen connector in mask should not be connected directly to any endotracheal or tracheostomy tube inlet.
- If oxygen is in use, do not use near a flame or heat source.
- Position tubing to avoid strangulation.
- Device contains components that may present choking hazard.
- Patient may become hypoxic if oxygen flow is interrupted.
- Use a flow rate \geq 4 LPM to prevent CO₂ accumulation and re-breathing.
- Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
- To reduce the risk of misconnections and patient injury, always trace tubing from gas source to the medical device before connecting.
- Do not reuse on more than one patient as this may increase risk of contamination and infection.
- When used as a nonrebreather mask, assure adequate and uninterrupted oxygen flow to patient. If oxygen flow is inadequate or interrupted, patient may not be able to inhale, or may not be able to inhale adequate volume.

Cautions: When using mask with reservoir bag, make sure bag is inflated before attaching to patient.

- Consult physician if patient develops infection, skin irritation or material sensitization.
- Do not place anything on supply tubing that may obstruct flow.

Material(s) Disclosure: Risks and Precautionary Measures related to Phthalates: This instruction pertains to the phthalate symbol marked on the device or its packaging. If this device is used for the treatment of children, or treatment of pregnant or nursing women; please note that the following types of procedures may increase the risk of exposure to phthalates: Exchange transfusion in neonates, total parenteral nutrition in neonates, multiple procedures in sick neonates, hemodialysis in peripubertal males, male fetus and male infant of pregnant women, and lactating women; and massive blood infusion into trauma patients. Although these procedures have the potential for increased risk of exposure, conclusive evidence of human health risks has not been established. As a precautionary measure, to reduce the potential for unnecessary exposures to phthalates, the product must be used in accordance with the instructions for use, and practitioners should refrain from using this product beyond the period of time the product is medically necessary or needed.

Residual Risks: Refer to warnings and cautions.

Necessary User Qualifications: This device is restricted to sale by or on the order of a physician. The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings and cautions as stated

EN in the IFU. Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

Instructions for Use: Mask is assembled as a high concentration, nonbreather mask with safety vent.

1. Attach one end of oxygen tubing to mask and the other end to oxygen source.
2. Set flow rate to 10 LPM. Occlude valve between mask and reservoir bag until bag fills up.
3. Place mask over patient's face covering both nose and mouth.
4. Slip strap over patient's head and adjust for comfort and fit. Pinch nose clip to provide a seal around the nose.
5. Adjust oxygen flow to ensure that the reservoir bag does not fully collapse during inspiration.
6. Check for leaks at connections and for proper placement of the mask.

To use as a rebreather, remove flexible disks located on side port and inside of mask on top of reservoir bag.

1. Attach tubing to oxygen source and set flow rate 8 L/min to 10 L/min.
2. Follow steps 3-6 above.



To use as a medium concentration mask, remove reservoir bag, connector with baffle, and flexible disc on side port. Then, insert oxygen connector (supplied loose in bag) and reconnect tubing.

1. Attach tubing to oxygen source and set to prescribed flow rate.
2. Place mask over patient's face covering both nose and mouth.

3. Slip strap over patient's head and adjust for comfort and fit. Pinch nose clip to provide a seal around the nose.
4. Check for leaks at connections and for proper placement of the mask.

Compatibility to Other Devices: AirLife tubing connectors are designed to work with gas outlets that comply to ISO 5359:2014. This includes, but is not limited to, stand-alone and portable oxygen concentrators, oxygen tanks, oxygen flowmeters, and bubble humidifiers. The female threaded connector, also referred to as U/Connect-it, is designed to be used with a male DISS 1240 oxygen fitting, 5 mm - 7 mm tapered oxygen fitting, or a 1/4 inch (6.35 mm) tapered oxygen fitting, e.g., nipple & nut adapter (Figure B).

1. Place the threaded connector between your thumb and forefinger below the threaded section (Figure C).
2. Position the connector's end with either the threaded male fitting or the tapered fitting.
 - a. For attachment to male fitting, firmly push the connector onto the oxygen fitting and turn clockwise until a secure fit is achieved. (Figure C)
 - b. For attachment to a tapered fitting, with slight pressure, slide the connector onto the gas fitting and turn clockwise until a secure fit is achieved.
3. Gently tug on tubing to ensure secure fit.

Safe Disposal: Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

Storage: Recommended normal storage temperature 0°C/+32°F to +50°C/+122°F.

Incident Reporting: Notice to the user and/or patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. If your physical symptoms worsen (e.g., increased shortness of breath, fever, dizziness) call your doctor.

Performance Characteristics:


	001203	001210	001268
Patient	Adult	Adult	Pediatric
Safety Vent	Yes	No	Yes
Oxygen Tubing	7 ft (2.1 m)	7 ft (2.1 m)	7 ft (2.1 m)

EN

REF Catalogue Number

LOT Lot Number


X Quantity


LATEX Not Made with Natural Rubber Latex

PHT Made with DEHP

i Consult Instructions for Use

Rx ONLY U.S. Federal Law restricts this device to sales by or on order of a physician

CLEAN, READY TO USE Clean, Ready to Use

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© 2019 Vyaire. Vyaire, the Vyaire Logo and AirLife are trademarks or registered trademarks of Vyaire Medical, Inc., or of its affiliates.
AirLife Australia Holding Pty. Ltd.
Vyaire Medical, Inc.