

**REF 0217**

UNITS - 1

**Rx Only**Not Made  
With Natural  
Rubber LatexNot Made  
With DEHP

CE 0482



VixOne™ Small Volume Nebulizer  
 Nebulizador para Pequeños Volúmenes VixOne™  
 Nebulizador de Pequeno Volume VixOne™  
 Nebuliseur à faible Volume VixOne™  
 Nebulizzatore per Volumi ridotti VixOne™  
 VixOne™ Vernebler für Kliene Volumina

**SUGGESTED PROCEDURE:**

1. These products are intended for use in delivering aerosol medications.
2. All products are to be used in accordance with the applicable prescription.
3. Remove the Nebulizer from the package and add the prescribed medication.
4. Assemble the Nebulizer Top to the Bowl per the Instruction Sheet.
5. Attach the Nebulizer to the components included in this kit and then connect the supply tube to the gas source 50 psi (413 kPa) max, compressed air and/or oxygen. Rotate the Threaded Nut (if included) 1/4 turn to ensure proper fit.
6. Ensure that all connections are secure, and then start the gas flow 4-10 LPM.
7. Confirm appropriate aerosol function.



Do not use oxygen near ignition source  
 not tested for anesthetic systems



Performance information may not apply to drugs supplied in suspension or high viscosity form, consult drug supplier.

**Westmed**


5580 S. Nogales Highway  
 Tucson, Arizona 85706 USA  
 Phone: 800-975-7987  
 Fax: 520-294-6061  
 www.westmedinc.com

ASSEMBLED IN MEXICO  
 LABEL PN 77731, Rev. 01

**LOT**

MT Promedt Consulting GmbH  
 Altenhofstr. 80  
 66386 St. Inbert, Germany

## Assembly Illustration



### Cleaning Instructions:

The VixOne™ Nebulizer must be disassembled and cleaned periodically in order to operate correctly. Westmed recommends disassembly and cleaning of all nebulizer components at the end of each treatment.

1. Wash and dry your hands.
2. Disassemble the Nebulizer .
3. **Acute Care Setting:** Rinse all nebulizer component parts per acute care setting cleaning protocol.
4. **Home Care Setting:** Rinse or wash all nebulizer component parts in warm soapy water. Rinse and allow to dry thoroughly.
5. Reassemble and store for next use.



### **CAUTION:**

Do Not Use Vinegar. Do not boil or microwave as high temperatures may melt some of the plastics parts.