

REF 100940

UNITS - 1
RX Only



Uni-HEART® Lo-Flo
Continuous Nebulizer
Adult Mask, 7' Tubing
w/Threaded Nut

Use with Flow Rates between 2-4 LPM

 **CAUTION:**

Use of sterilization or cleaning solutions
may leave this device nonfunctional.



MT Promedt Consulting GmbH
Altenhofstr. 80
66386 St. Ingbert, Germany
Patent Pending

Westmed 

5580 S. Nogales Highway
Tucson, Arizona 85706 USA
Phone: 800-975-7987
Fax: 520-294-6061
www.westmedinc.com

U.S. Patent 5,287,847
Made In U.S.A
LABEL PN 79539, Rev. 01
Page 1 of 2



LOT



(01)00709078000294



(17)YYMMDD



(10)LLLLLL

A. Instructions for Continuous Nebulization

1. Remove Uni-HEART® from the package.
2. Connect Uni-HEART® to a flow meter using oxygen supply tubing.
CAUTION: Ensure tubing connections are secured.
3. Pour prescribed medication into reservoir through the aerosol outlet opening or connect to your IV setup via the IV Port as shown in Section B diagram.
NOTE: The top of the Uni-HEART® is welded to the reservoir and cannot be opened.
4. Adjust flow meter to desired flow rate of 2 L/min. **Output is approximately 4.5 mL/hr.** For higher output, set flow meter to 4L/min. **Output is approximately 9 mL/hr.**
IMPORTANT: The listed flow rate and outlet values are only intended as general guideline ($\pm 20\%$). This is not a substitute for the clinical judgment of qualified personnel.
5. Connect Aerosol Mask directly to the outlet of Uni-HEART®. Place Aerosol Mask on patient. Instruct patient to breathe normally through Aerosol Mask.
NOTE: If additional oxygen is desired, use a Nasal Cannula under the Aerosol Mask.
6. Check liquid consumption q 1 hr. x 2 to determine output of setup.
7. Increase or decrease the nebulizer flow rate by small increments to obtain desired output.
8. Aerosol performance is based on using a 50 psi pressure source.
NOTE: The output of each Uni-HEART® can be affected by several factors. It is important for users to follow this "user correctable" approach when utilizing continuous nebulization therapy (CNT).

B. Continuous Nebulization Setup with Uni-HEART® Nebulizer

