
en

BiNAPS® Supplemental Package Insert

Operating Environment Conditions

Temperature: 5°C to 40°C

Relative Humidity: 15 to 95%

Atmospheric Pressure: 80-106 kPa

Handling & Storage Environment Conditions

Temperature: 5°C to 40°C

Relative Humidity: 15 to 95%

Atmospheric Pressure: 50-106 kPa

Salter Labs® Sleep Diagnostics Nasal Cannula used with BiNAPS® device are considered a Single-Use, Disposable, Type B Applied Part

Refer to Sleep Diagnostics Equipment for system connection instructions.

Cannula may be discarded in trash after each use.

BiNAPS® interface box is reusable and may be disposed of after 5 years of use. Dispose of product according to applicable local regulations.

NOTE: BiNAPS® is intended only for connection to a Medical Electronics system complying with IEC 60601-1.

BiNAPS® is classified as an internally powered, continuous operation medical electronics device with sleep diagnostics cannula that are Type B Applied Parts.

WARNING: No modification of this equipment is allowed.

CAUTION: Do not touch Airflow and Snore lead contacts and the patient simultaneously.

NOTE: Do not use cellular phones, tablets or other non-medical electronics devices that can emit wireless electromagnetic signals near the device. All Sleep Diagnostics Equipment should be positioned at an appropriate distance away from the BiNAPS® device to minimize interference that could affect proper operation.

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