
O2/EtCO2 Divided Nasal Cannula

REF 4700, 4701, 4702, 4703, 4704, 4705, 4706, 4707, 4708, 4787 and Associated SO PNs

EN O2/EtCO2 Divided Nasal Cannula

DEVICE DESCRIPTION: The device is a non-sterile, disposable, single patient use cannula that provides supplemental oxygen while sampling patients exhaled gases.

INTENDED PURPOSE: The O2/EtCO2 Divided Nasal Cannula is intended to be used where low flow oxygen is administered and exhaled gas is monitored.

INDICATIONS FOR USE: The O2/EtCO2 Divided Nasal Cannula is intended to be used with spontaneously breathing patients who require supplemental oxygen and exhaled gas monitored.

ENVIRONMENT: Hospital, sub-acute, pre-hospital.

PATIENT TARGET GROUP: Infant, pediatric, adult.

EXPECTED CLINICAL BENEFITS: Capnography provides early detection and warning of changes in a patient's ventilatory status and cardiac perfusion. This enables clinicians to adjust medication dosage or apply therapeutic intervention to prevent hypoventilation and apnea. Supplemental oxygen up to 5 LPM can be delivered to treat hypoxemia measured by changes in SpO2, blood gases and relief of shortness of breath.

CONTRAINDICATIONS: None Known.

- WARNINGS:**
- If oxygen is in use, do not use near flame or heat source.
 - Position tubing to avoid strangulation.
 - Patient may become hypoxic if oxygen flow is interrupted.
 - Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
 - To reduce the risk of misconnections and patient injury, always trace tubing from gas source to the medical device before connecting.
 - Contains Tris(nonylphenyl) phosphite (TNPP) above 0.1% w/w.

- CAUTIONS:**
- Consult physician if patient develops infection, skin irritation or material sensitization.
 - Do not place anything on supply tubing that may obstruct flow.
 - Not for use on patients with obstruction of both nostrils.
 - Accuracy of CO2 monitoring can be diminished during the time the patient is mouthbreathing or if the patient has a deviated septum.

RESIDUAL RISKS: Refer to Warnings and Cautions.

NECESSARY USER QUALIFICATIONS: This device is restricted to sale by or on the order of a physician. There are no additional requirements for special facilities, special training, or particular qualifications for the use of this device.

INSTRUCTIONS: Follow instructions to avoid risk of infection and contamination.

- Connect the gas sampling line to the capnography monitor.
- Connect the oxygen tubing to oxygen source.
- Set prescribe flow.
- Check for flow from prongs.
- Insert prongs into nares.
- Wrap headset loop up and over ears.
- Squeeze the sides of the bolo and glide up under the chin.
- Check monitor display for EtCO2 value.

This product is disposable and is not intended for prolonged use greater than 7 days.

SAFE DISPOSAL: Dispose of device in accordance with local, state or national regulations.

INCIDENT REPORTING: Contact your country's Competent Authority and Manufacturer to report any serious incident.

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