REF: 4MSF8

EN

DEVICE DESCRIPTION: ETCO2 Filtered Sampling Line
The device is non-sterile, disposable, single-patient use sampling line that enables sampling of exhaled gases. The sampling line available in different end connectors.

INTENDED PURPOSE: The Sample Lines are intended to be used where exhaled gas is monitored.

INDICATIONS: The Sample Lines are intended to be used to sample exhaled gas in intubated or non-intubated patients requiring expired gas monitoring.

ENVIRONMENT: Hospital, sub-acute, pre-hospital

PATIENT TARGET GROUP: Infant, pediatric, adult
Shelf life: 3 years

EXPECTED CLINICAL BENEFITS:
1) Notifies clinicians of changes in the patient’s ventilator status
2) Provides information on the effectiveness of cardiac resusciation
3) Provides verification of placement of ETT in trachea.

CONTRAINDICATIONS: None known
WARNINGS: Position tubing to avoid strangulation

CAUTIONS: Consult physician if patient develops infection, skin irritation or material sensitization
Do not place anything on supply tubing that may obstruct flow
Excessive patient secretion or a build-up of liquids in the tubing may occlude the sampling and require more frequent replacement.

INSTRUCTIONS: Contact your country’s Competent Authority to report any serious incident
Dispose of device in accordance with local, state or national regulations
Follow instructions for proper care of device to prevent risk of contamination and injury.
This product is disposable and is not intended for prolonged use
Check for leaks at connections and for proper placement prior to use
Check for proper position to ensure detection of exhaled gas.

Replace sampling line when blockage is indicated by the capnography monitor.

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