Oxygen Cannula- Soft Infant

840080

EN Oxygen Cannula- Soft Infant

DEVICE DESCRIPTION: The device is non-sterile, disposable and for single-patient use. An oxygen cannula is a device with two prongs used to deliver supplemental oxygen to the patient's nostrils.

INTENDED PURPOSE: A nasal oxygen cannula is a two-pronged device used to administer oxygen to a patient through both nostrils. Depending on the cannula style flow rates will vary from > 0 LPM to ≤ 3 LPM.

INDICATIONS: For patient's prescribed supplemental oxygen by a nasal cannula. To treat hypoxia, relieve shortness of breath and pre-oxygenate for surgical procedures in spontaneously breathing patients.

ENVIRONMENT: Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

PATIENT TARGET GROUP: Newborns, infant.

EXPECTED CLINICAL BENEFITS:

- Nasal cannulas are a conduit that deliver supplemental oxygen to a patient's nasopharynx to treat or prevent hypoxemia in acute and chronic respiratory conditions. It is monitored by improvement in SpO₂ values, PaO₂ values, vital signs, and/or work of breathing.
- Supplemental oxygen helps reduce COPD complications by stabilizing pulmonary hypertension, reducing secondary polycythemia, and decreasing arrhythmias.
- Supplemental oxygen can decrease anxiety, improve sleep, improve mental alertness, and/or improve health-related quality of life in people with chronic respiratory disease.
- Using oxygen during exercise can improve endurance, heighten performance, and/or decrease the sensation of breathlessness. This helps patients complete physical activities, improve ability to walk and/or perform other physical activities for an extended period.
- The 3-channel tubing helps reduce kinking and allows for greater than 75% of initial flow to be delivered should the tubing become kinked.

CONTRAINDICATIONS: None Known.

WARNINGS:

- Patient may become hypoxic of oxygen flow is interrupted.
- If oxygen is in use, do not use near flame or heat source.
- Position tubing to prevent patient from becoming tangled in tubing or strangulation.
- Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
- To reduce the risk of misconnections and patient injury, always trace tubing and lines back to their point of origin to verify that correct connections are made.
- High levels of oxygen may cause oxygen-induced hypoventilation, oxygen-induced hypercapnia, or oxygen toxicity. In premature infants, it may also cause retrolental fibroplasia.

- The facepiece contains Tris (nonylphenol) phosphite (TNPP) above 0.1% w/w.

CAUTIONS:

- Consult physician if the patient develops an infection, skin irritation, or material sensitization.
- Flow rates above 3 LPM may cause increased resistance or back pressure in the tubing.
- To prevent interruption of therapy or injury, do not let children play with the oxygen tubing.
- To prevent excess moisture accumulation in the tubing, do not use with heated humidification.
- To prevent equipment damage and/or interruption of therapy, do not let pets play with the oxygen tubing.
- Do not place tubing under rugs or other objects as it may obstruct flow.
- To prevent from becoming tangled in tubing, position the tubing in front of you before turning to sit or stand.

RESIDUAL RISKS: Refer to Warnings and Cautions.

NECESSARY USER QUALIFICATIONS: This device is restricted to sale by or on the order of a physician. There are no additional requirements for special facilities, special training, or particular qualifications for the use of this device.

INSTRUCTIONS:

- Select the correct size cannula for the patient.
- Connect to the oxygen source.
- Adjust the flow control knob to the prescribed liter flow. Check for flow from the prongs.
- Position the nasal cannula with the prongs facing upward and curved toward the face. Insert the two prongs into the nostrils.
- Wrap the headset loop up and over the ears and tuck the loop under the chin, or adjust and secure the headset loop behind the head.
- Squeeze the sides of the bolo and glide the bolo up under the chin, or until snug around the head.
- Leave enough space to fit at least the width of two fingers between the bolo and chin.

- Discard and replace cannula if soiled or damaged. Do not sterilize.

SAFE DISPOSAL: Dispose of device in accordance with local, state or national regulations.

INCIDENT REPORTING: Contact your country's Competent Authority and Manufacturer to report any serious incident.

PERFORMANCE CHARACTERISTICS:

Cannula	16SOFT-PRE	16SOFT-NEO	16SOFT-INF	16SOFT-TOD
Flow Rate	3 LPM	3 LPM	3 LPM	3 LPM
Prong Size	0.060"	0.060"	0.073"	0.073"
Patient Population	Infant (Premature)	Infant (Neonate)	Infant	Infant (Toddler)

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MT Promedt Consulting GmbH