## Oxygen Cannula- Soft Adult / Pediatric

## 840081

## EN Oxygen Cannula- Soft Adult / Pediatric

**DEVICE DESCRIPTION:** The device is non-sterile, disposable and for single-patient use. An oxygen cannula is a device with two prongs used to deliver supplemental oxygen to the patient’s nostrils.

**INTENDED PURPOSE:** A nasal oxygen cannula is a two-pronged device used to administer oxygen to a patient through both nostrils. Depending on the cannula style flow rates will vary from > 0 LPM to ≤ 6 LPM.

**INDICATIONS:** For patient’s prescribed supplemental oxygen by a nasal cannula. To treat hypoxia, relieve shortness of breath and pre-oxygenate for surgical procedures in spontaneously breathing patients.

**ENVIRONMENT:** Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

**PATIENT TARGET GROUP:** The intended population includes infant, pediatric and adult patients. Spontaneously breathing individuals with a variety of breathing conditions requiring supplemental oxygen.

**EXPECTED CLINICAL BENEFITS:**

- Nasal cannulas are a conduit that deliver supplemental oxygen to a patient’s nasopharynx to treat or prevent hypoxemia in acute and chronic respiratory conditions. It is monitored by improvement in SpO<sub>2</sub> values, PaO<sub>2</sub> values, vital signs, and/or work of breathing.
- Supplemental oxygen helps reduce COPD complications by stabilizing pulmonary hypertension, reducing secondary polycythemia, and decreasing arrhythmias.
- Supplemental oxygen can decrease anxiety, improve sleep, improve mental alertness, and/or improve health-related quality of life in people with chronic respiratory disease.
- Using oxygen during exercise can improve endurance, heighten performance, and/or decrease the sensation of breathlessness. This helps patients complete physical activities, improve ability to walk and/or perform other physical activities for an extended period.
- The 3-channel tubing helps reduce kinking and allows for greater than 75% of initial flow to be delivered should the tubing become kinked.

**CONTRAINDICATIONS:** None Known.

**WARNINGS:**

- Patient may become hypoxic of oxygen flow is interrupted.
- If oxygen is in use, do not use near flame or heat source.
- Position tubing to prevent patient from becoming tangled in tubing or strangulation.
- Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
- To reduce the risk of misconnections and patient injury, always trace tubing and lines back to their point of origin to verify that correct connections are made.
- High levels of oxygen may cause oxygen-induced hypoventilation, oxygen-induced hypercapnia, or oxygen toxicity. In premature infants, it may also cause retrolental fibroplasia.

- The facepiece contains Tris (nonylphenol) phosphite (TNPP) above 0.1% w/w.

**CAUTIONS:**

- Consult physician if the patient develops an infection, skin irritation, or material sensitization.
- Flow rates above 6 LPM may cause increased resistance or back pressure in the tubing.
- To prevent interruption of therapy or injury, do not let children play with the oxygen tubing.
- To prevent excess moisture accumulation in the tubing, do not use with heated humidification.
- To prevent equipment damage and/or interruption of therapy, do not let pets play with the oxygen tubing.
- Do not place tubing under rugs or other objects as it may obstruct flow.
- To prevent from becoming tangled in tubing, position the tubing in front of you before turning to sit or stand.

**RESIDUAL RISKS:** Refer to Warnings and Cautions.

**NECESSARY USER QUALIFICATIONS:** This device is restricted to sale by or on the order of a physician. There are no additional requirements for special facilities, special training, or particular qualifications for the use of this device.

**INSTRUCTIONS:**

- Select the correct size cannula for the patient.
- Connect to the oxygen source.
- Adjust the flow control knob to the prescribed liter flow. Check for flow from the prongs.
- Position the nasal cannula with the prongs facing upward and curved toward the face. Insert the two prongs into the nostrils.
- Wrap the headset loop up and over the ears and tuck the loop under the chin, or adjust and secure the headset loop behind the head.
- Squeeze the sides of the bolo and glide the bolo up under the chin, or until snug around the head.
- Leave enough space to fit at least the width of two fingers between the bolo and chin.

- Discard and replace cannula if soiled or damaged. Do not sterilize.

**SAFE DISPOSAL:** Dispose of device in accordance with local, state or national regulations.

**INCIDENT REPORTING:** Contact your country’s Competent Authority and Manufacturer to report any serious incident.

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
Salter Labs 30 Spur Drive El Paso, TX 79906 USA Made in Mexico www.Sun-Med.com	UK REP	MT Promedt Consulting Ltd. Beaver House, 23-38 Hythe Bridge Street Oxford, OX1 2EP United Kingdom
EC REP	MT Promedt Consulting GmbH Ernst-Heckel-Straße 7 66386 St. Ingbert Germany	MedEnvoy UK Limited 85, Great Portland Street – First Floor London, W1W 7LT United Kingdom
CH REP	MedEnvoy NL-IM-000000248 Prinses Margrietplantsoen 33 Suite 123 2595 AM, The Hague The Netherlands	MedEnvoy Switzerland Gotthardstrasse 28 6302 Zug Switzerland

**CE 2797**

840081 Rev. D 2023/06