EN

English

Salter Labs™ Venturi Oxygen Mask

DEVICE DESCRIPTION:

The device is non-sterile, disposable and for single-patient multi-use.

- Percent O2 Lock: A single unit variable venturi system designed to deliver seven fixed oxygen concentrations via an adult mask to a spontaneously breathing patient. The concentration of oxygen range from 24 to 50%.
- Accu-Flow: A single unit interchangeable venturi system designed to deliver five fixed oxygen concentration via an adult mask to a spontaneously breathing patient. The concentration of oxygen range from 24 to 40%.

INTENDED PURPOSE:

Percent O2 Lock Mask is a device placed over a patient's nose or mouth to administer oxygen.

Accu-Flow Entrainment Mask is a device placed over a patient's nose or mouth to administer oxygen.

INDICATIONS:

For patients requiring supplemental oxygen. To treat hypoxia, relieve shortness of breath and pre-oxygenate for surgical procedures in spontaneously breathing patients.

ENVIRONMENT:

Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

PATIENT TARGET GROUP:

Pediatric, adult.

INTENDED USER:

Trained medical professionals.

EXPECTED CLINICAL BENEFITS:

The oxygen mask delivers up to 8 LPM of oxygen to reverse or prevent hypoxemia, and provide relief of shortness of breath in oral and nasal breathers.

CONTRAINDICATIONS:

None Known.

WARNINGS:

- If oxygen is in use, do not use near flame or heat source.
- Position tubing to avoid strangulation.
- Device contains components that may present choking hazard.
- Patient may become hypoxic if oxygen flow is interrupted.
- Not intended to be used with oxygen concentrator.
- Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
- To reduce the risk of misconnections and patient injury, always trace tubing from gas source to the medical device before connecting.

CAUTIONS:

- Consult physician if patient develops infection, skin irritation or material sensitization.
- Do not place anything on supply tubing that may obstruct flow.
- This is a single-patient multi-use product, reuse on more than one patient may cause a risk of cross-contamination.

RESIDUAL RISKS:

Refer to Warnings and Cautions.

NECESSARY USER QUALIFICATIONS:

The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings, and cautions as stated in the IFU.

Rx Only: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

INSTRUCTIONS:

Follow instructions to avoid risk of infection and contamination.

Percent O2 Lock:

- 1) To set delivered oxygen %, pull the sleeve back until the spring is completely compressed, then turn the sleeve to the prescribed setting.
- 2) Release the spring to lock sleeve in place.

- 3) Set the recommended flow rate.
- 4) Place mask over nose and mouth.
- 5) Adjust elastic strap and nose bar for a snug fit.
- 6) Check for leaks at connection and check the medical device works properly.
- 7) To add humidity, slip the humidifier cup over the venturi/diluter jet and connect aerosol tubing to humidifier cup.

Accu-Flow:

- 1) Select the diluter jet for the desired FiO2.
- 2) Attach to corrugated tube on mask.
- 3) Set the recommended flow rate.
- 4) Place mask over nose and mouth.
- 5) Adjust elastic strap and nose bar for a snug fit.
- 6) Check for leaks at connection and check the medical device works properly.
- 7) To add humidity, slip the humidifier cup over the venturi/diluter jet and connect aerosol tubing to humidifier cup.

This product is disposable and is not intended for prolonged use greater than 14 days.

CLEANING/REUSE INSTRUCTIONS

A discontinued use of this product on the same patient, with periods of use separated by a period of non-use, is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.

SAFE DISPOSAL:

Decontaminate and dispose of all potentially biohazardous material. Dispose of device in accordance with local, state or national regulations.

INCIDENT REPORTING:

Notice to the user that any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myairlife.com, and the competent authority of the Member State in which the user and/or patient is established. If the patient's physical symptoms worsen or experience a sudden change in their condition (e.g., increased shortness of breath, fever, dizziness) call the physician.

PERFORMANCE CHARACTERISTICS:

- For use with Flow Rate 3 - 8 LPM
- For use with Bonded O2 Tubing adaptor

The unit should provide the following percentage of oxygen at these stated flows:

Percent O2 Lock			
3 LPM:	24 ± 1%	6 LPM:	35 ± 2%
	26 ± 1%		40 ± 2%
	28 ± 1%		50 ± 2%
	31 ± 1%		

Accu-Flow			
3 LPM:	24 ± 1%	8 LPM:	35 ± 2%
	26 ± 1%		40 ± 2%
6 LPM:	31 ± 1%		

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MT Promedt Consulting GmbH
MT Promedt Consulting Ltd.
MedEnvoy UK Limited