T-Piece Jet Nebulizer

T-Piece Jet Nebulizer

DEVICE DESCRIPTION: The device is non-sterile, disposable and for single-patient reuse. Hand-held pneumatic small volume nebulizer with reservoir which connects to a mouthpiece or aerosol mask. The one-way valve (mask) opens during exhalation. The one-way valve helps to enhance aerosol delivery. The recommended driving gases are air or oxygen. The device may be used with a compressor or a gas source capable of delivering a flow rate of 4 to 8 LPM at a nominal operating pressure between 25 to 52 psi.

INTENDED PURPOSE: The Nebulizer is intended for use where liquid (solutions, suspensions, and emulsions) medications are to be delivered to a patient in aerosol form

INDICATIONS FOR USE: Any spontaneously breathing patient requiring aerosolized respiratory medication.

ENVIRONMENT: Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

PATIENT TARGET GROUP: Infant, pediatric, adult.

EXPECTED CLINICAL BENEFITS: Small volume nebulizer aerosolizes medication allows for quicker and easier absorption of medication into the lungs. Aerosolized bronchodilators help to relieve bronchospasms, decrease the frequency of exacerbations and improve lung function; therefore, relieving shortness of breath, and improving health status and exercise tolerance. Other aerosolized medications help to reduce sputum viscosity, decrease airway inflammation and improve pulmonary function.

CONTRAINDICATIONS: None Known.

WARNINGS:

- Position tubing to avoid strangulation.
- Device contains components that may present choking hazard.
- Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
- To reduce the risk of misconnections and patient injury, always trace tubing from gas source to the medical device before connecting.

CAUTIONS:

- Clinician or caregiver should avoid inhaling aerosolized medication.
- Prior to cleaning, during dismantling and reassembling, and after each use, remove the nebulizer from the driving gas source.

RESIDUAL RISKS: Refer to Warnings and Cautions.

NECESSARY USER QUALIFICATIONS: This device is restricted to sale by or on the order of a physician. There are no additional requirements for special facilities, special training, or particular qualifications for the use of this device.

INSTRUCTIONS: Follow cleaning instructions to avoid risk of infection and contamination.

1) Insert mouthpiece into large port of tee-adpater.

2) If applicable, attach reservoir tube to other side.

3) Screw cap onto medication cup.

4) Attach the tubing to the bottom port of the medication chamber.

5) Ensure green cone is seated on bottom of medication cup.

6) Add medication.

7) Attach tee-piece assembly to nebulizer.

8) Set gas flow rate between 4 LPM and 8 LPM (25-52 PSI).

9) Ensure mist is visible.

10) Instruct patient to keep lips sealed around mouthpiece during treatment.

CLEANING INSTRUCTIONS: This product is disposable and is not intended for prolonged use greater than 7 days. After each use, thoroughly rinse the nebulizer with sterilized water. DO NOT submerge the supply tubing in water. Disassemble all nebulizer components; do not remove any of the one-way valves. Rinse with sterilized water, shake out any excess. Allow the device to dry, re-assemble and place in a dry, dust free location. Replace the nebulizer every 7 days, or after 21 cleaning cycles, whichever is first.

SAFE DISPOSAL: Dispose of device in accordance with local, state or national regulations.


INCIDENT REPORTING: Contact your country's Competent Authority and Manufacturer to report any serious incident.

PERFORMANCE CHARACTERISTICS:

Nebulizer performance characteristics stated in the instructions for use are based upon testing that utilizes adult ventilatory patterns and are likely to be different for pediatric or infant populations. The liquid medication for nebulization, in particular a suspension and/or high-viscosity solution, can alter the particle size distribution curve, the mass median aerodynamic diameter (MMAD), aerosol output, and/or aerosol output rate, which can then be different from what has been disclosed. All stated performance characteristics were derived from the nebulization of 4mL of albuterol 0.1% (M/V) concentration in 0.9% sodium chloride solution.

| T-Piece Jet Nebulizer | | |
|---|--|-------|
| Test Solution: 4ml of albuterol 0.1% (M/V) concentration in 0.9% sodium chloride solution | | |
| Performance Characteristic | Source Pressure: 50 psi | |
| Aerosol Output Rate (ml/min) | 4 LPM | 6 LPM |
| Aerosol Output Rate (ug/min) | 0.18 | 0.27 |
| Aerosol Output Rate (ug/min) | 181 | 270 |
| Aerosol Output (ug) | 0.81 | 0.92 |
| Aerosol Output (ug) | 807 | 916 |
| Percentage of Fill Volume Emitted | 20.2% | 22.9% |
| Percentage of Fill Volume Emitted per minute | 4.5% | 6.7% |
| Nebulization Time (min) | 14.0 | 7.83 |
| Residual Volume (ml) | 1.67 | 1.68 |
| Aerosol Output Fraction | 0.81 | 0.92 |
| Particle Size and Source Pressure of 50 psi at 8 LPM | Mass Median Aerodynamic Diameter (MMAD): 3.73 (um) | |
| Geometric Standard Deviation (GSD): 2.87 | | |

| Stage No | Effective Cut-Off Dia. (um) | Cumulative % particle of albuterol under size | T-Piece Jet Nebulizer | | |
|----------|-----------------------------|---|-----------------------|----------|----------|
| | | | Sample 1 | Sample 2 | Sample 3 |
| 8 | 0.00 | 0.0% | 0.0% | 0.0% | 0.0% |
| 7 | 0.68 | 8.5% | 8.1% | 8.8% | 9.1% |
| 6 | 0.90 | 13.1% | 13.5% | 11.9% | 13.0% |
| 5 | 1.23 | 22.0% | 23.3% | 24.5% | 24.0% |
| 4 | 2.87 | 40.7% | 44.2% | 36.9% | 40.9% |
| 3 | 4.85 | 60.7% | 64.0% | 57.5% | 60.4% |
| 2 | 7.50 | 74.0% | 76.9% | 70.3% | 74.8% |
| 1 | 8.27 | 79.0% | 83.4% | 74.2% | 79.5% |

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