Nebulizer In-Line SVN

EN Nebulizer In-Line SVN

DEVICE DESCRIPTION: The device is non-sterile, disposable and for single-patient use. A pneumatic small volume nebulizer that fits into the breathing circuit to deliver medicated aerosols to the patient.

INTENDED PURPOSE: Nebulizer is intended to generate aerosols that are delivered directly to the patient.

INDICATIONS FOR USE: Any spontaneously breathing patient requiring aerosolized respiratory medication.

ENVIRONMENT: Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

PATIENT TARGET GROUP: Infant, pediatric, adult.

EXPECTED CLINICAL BENEFITS: Small volume nebulizer aerosolize medication allows for quicker and easier absorption of medication into the lungs. Aerosolized bronchodilators help to relieve bronchospasms, decrease the frequency of exacerbations and improve lung function, therefore, relieving shortness of breath, and improving health status and exercise tolerance. Other aerosolized medications help to reduce sputum viscosity, decrease airway inflammation and improve pulmonary function.

CONTRAINDICATIONS: None Known.

WARNINGS:

- Position tubing to avoid strangulation.
- Device contains components that may present choking hazard.
- Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
- To reduce the risk of misconnections and patient injury, always trace tubing from gas source to the medical device before connecting.
- Contains Octamethylcyclotetrasiloxane above 0.1% w/w.

CAUTIONS:

• Clinician or caregiver should avoid inhaling aerosolized medication.

RESIDUAL RISKS: Refer to Warnings and Cautions.

NECESSARY USER QUALIFICATIONS: This device is restricted to sale by or on the order of a physician. There are no additional requirements for special facilities, special training, or particular qualifications for the use of this device.

INSTRUCTIONS: Follow cleaning instructions to avoid risk of infection and contamination.

- Screw cap onto medication cup.
- Attach the tubing to the bottom port of the medication chamber.
- Ensure green cone is seated on bottom of medication cup.
- Add medication.
- Insert mask nebulizer port into ventilator tee assembly.
- Set gas flow rate between 5 LPM and 8 LPM (25-52 PSI).
- Ensure mist is visible.
- Adjust ventilator setting as needed to compensate for additional flow.
- At end of treatment, turn off flow and readjust ventilator settings.

CLEANING INSTRUCTIONS: This product is disposable and is not intended for prolonged use greater than 7 days. After each use, thoroughly rinse the nebulizer with sterilized water. DO NOT submerge the supply tubing in water. Disassemble all nebulizer components; do not remove any of the one-way valves. Rinse with sterilized water, shake out any excess. Allow the device to dry, re-assemble and place in a dry, dust free location. Replace the nebulizer every 7 days, or after 21 cleaning cycles, whichever is first.

SAFE DISPOSAL: Dispose of device in accordance with local, state or national regulations.

INCIDENT REPORTING: Contact your country’s Competent Authority and Manufacturer to report any serious incident.

PERFORMANCE CHARACTERISTICS:

Nebulizer performance characteristics stated in the instructions for use are based upon testing that utilizes adult ventilatory patterns and are likely to be different for pediatric or infant populations. The liquid medication for nebulization, in particular a suspension and/or high-viscosity solution, can alter the particle size distribution curve, the mass median aerodynamic diameter (MMAD), aerosol output, and/or aerosol output rate, which can then be different from what has been disclosed. All stated performance characteristics were derived from the nebulization of 4mL of albuterol 0.1% (M/V) concentration in 0.9% sodium chloride solution.

Nebulizer In-Line SVN			
Test Solution: 4mL of albuterol 0.1% (M/V) concentration in 0.9% sodium chloride solution			
	Source Pressure: 50 psi		
Performance Characteristic	4 LPM	6 LPM	8 LPM
Aerosol Output Rate (ml/min)	0.18	0.27	0.31
Aerosol Output Rate (µg/min)	181	270	311
Aerosol Output (ml)	0.81	0.92	0.97
Aerosol Output (µg)	807	916	967
Percentage of Fill Volume Emitted	20.2%	22.9%	24.2%
Percentage of Fill Volume Emitted per minute	4.5%	6.7%	7.8%
Nebulization Time (min)	14.0	7.83	6.25
Residual Volume (ml)	1.67	1.68	1.51
Aerosol Output Fraction	0.81	0.92	0.97
Particle Size Data at Source Pressure of 50 psi at 8 LPM			
Mass Median Aerodynamic Diameter (MMAD)	3.73 (µm)		
Geometric Standard Deviation (GSD)	2.87		

Nebulizer In-Line SVN					
Efficacy (Cut-Off dia. (µm))					
Stage No	Effective	Nebulizer In-Line SVN			
		Cumulative % particle mass of albuterol under size			
		Mean	Sample 1	Sample 2	Sample 3
8	0.00	0.0%	0.0%	0.0%	0.0%
7	0.68	8.5%	8.1%	8.8%	9.1%
6	0.90	13.1%	13.5%	13.9%	13.9%
5	1.22	22.9%	23.3%	23.6%	25.0%
4	2.87	40.7%	44.2%	36.9%	40.9%
3	4.85	60.7%	65.0%	57.5%	60.4%
2	7.50	74.0%	76.9%	70.3%	74.8%
1	8.27	79.0%	83.4%	74.2%	79.5%

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MT Promedt Consulting GmbH
MT Promedt Consulting Ltd.
MedEnvoy UK Limited