Humidifier 7100

EN **Humidifier 7100**

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DEVICE DESCRIPTION: The device is non-sterile, disposable and for single patient use. The bubble humidifier consists of a cylindrical plastic bottle to hold water, a diffuser which gas passes through and a lid with safety pressure relief valve.

INTENDED PURPOSE: The bubble humidifier is intended to add moisture to breathing gases for administration to patients in homecare, hospital, extended care and hospice.

INDICATIONS: For patients who require humidificaton of supplemental breathing gases.

ENVIRONMENT: Hospitals, homecare, extended care facilities, and hospice environments.

PATIENT TARGET GROUP: Infant, pediatric, adult.

EXPECTED CLINICAL BENEFITS: The bubble humidifier adds ≥10 mg H₂O/L of moisture to the gas flow. The humidification helps to relieve symptoms caused by dry gas, such as sinus irritation and nasal congestion.

CONTRAINDICATIONS: None known.

⚠️ WARNINGS:

- If oxygen is in use, do not use near flame or heat source.
- Using flow rates outside the recommend flow rate range may affect the function of the safety pop-off valve.
- Patient may become hypoxic if oxygen flow is interrupted.
- Leaks may partially or fully reduce oxygen being delivered to the patient.
- The humidifier is not suitable for patients whose upper airways have been bypassed by an endotracheal tube or tracheotomy.
- Contains lead above 0.1% w/w in brass alloy component.

CAUTIONS:

- Follow cleaning instructions to avoid risk of infection and contamination.

RESIDUAL RISKS: Refer to Warnings and Cautions.

NECESSARY USER QUALIFICATIONS: This device is restricted to sale by or on the order of a physician. There are no additional requirements for special facilities, special training, or particular qualifications for the use of this device.

INSTRUCTIONS: Keep bubble humidifier upright during use.

- Add water to the max level (290 mL). Sterile or distilled water is recommended to reduce possibility of contamination or scaling. If using tap water, boil and cool before using. Replace and tighten the lid and check that it is sealed properly. Change water daily. Maintain water above min level (70 mL).
- Attach bubble humidifier to oxygen source.
- To test for leaks: Important test instructions refer to graphics. Set flow to 3 LPM and block outlet port. Humidifier should sound in 5 seconds or less. If not, check for loose connections.
- Attach the tubing. If the tubing has an obstruction to flow, it will cause the bubble humidifier to emit a sound. Pop-off will occur when pressure reaches 2.5-3.5 PSI.

CLEANING INSTRUCTIONS:

- Do not sterilize.

•The useful life of the Bubble Humidifiers is based on the environment of use. For home use, the device is to be replaced every 90 days (3 months) or sooner if it malfunctions or if the user is unable to remove mineral deposits by cleaning.

- For hospital use, the bubble humidifier is to be replaced every 14 days.

SAFE DISPOSAL: Dispose of device in accordance with local, state or national regulations.

INCIDENT REPORTING: Contact your country’s Competent Authority and Manufacturer to report any serious incident.

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