



AVANOS*




BALLARD* 500 PEDIATRIC T TRACHEOSTOMY HEAT AND MOISTURE EXCHANGER

SPECIFICATIONS:

Tidal volume range - up to: 500 ml
Breathing Resistance: .25 cm H₂O at 30 l/min
Deadspace: 12 ml

 Single Use Only	Rx Only	Storage Temperature: -30°C / -22°F +40°C / +104°F	1 Qty	 1157
--	---------	---	-------	--

 Pharma Systems AB, Rubanksgatan 9, SE-741 71, Knivsta, Sweden.
Distributed in the USA by Avanos Medical Sales, LLC, 5405 Windward Parkway,
Alpharetta, GA 30004 USA. In USA, 1-844-428-2667.

www.avanos.com

Product of Sweden

INTENDED USE: The BALLARD* 500 Pediatric T HME is a single-use heat and moisture exchanger which transfers heat and humidity from the gas expired by the patient to the inspired gases. The Pediatric 500 T HME is intended for use on patients who are breathing spontaneously through a tracheostomy tube.

CONTRAINDICATION: Do not use the 500 Pediatric T HME on patients producing fulminating secretions in their airways and lungs.

DIRECTIONS FOR USE:

- Change the HME every 24 hours, or sooner if necessary
- Not for use with nebulizers or active humidifiers

CAUTIONS:

- During the use of the 500 Pediatric T HME, the patient should be closely monitored and proper airway care administered if complications arise.
- Before use, verify that the 500 Pediatric T HME has no occlusions and that air will flow through it.

WARNING:

- The 500 Pediatric T HME is designed for single use only and must not be cleaned and reused.

CONNECTORS ISO STD:

Patient: 15mm F

*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.
© 2018 AVNS. All rights reserved.

06-M1-036-01

 LOT

Use By 

(01)00609038938806

