

BALLARD* 1000 HEAT AND MOISTURE EXCHANGE FILTER

Instructions for Use

Bacterial/ Viral Filtration

Indication

The BALLARD* 1000 HMEF can be used to humidify inspired gases and to protect the patient from various types of airborne bacteria.

Contraindication

The BALLARD* 1000 HMEF should not be used together with active humidifiers or during aerosolized medication delivered via the airways.

Instructions for Use

- Change the HMEF every 24 hours, or sooner if necessary.
- Dead space bacterial filter: 55 ml
- Not for use with nebulizers or active humidifiers.

⚠ Cautions

- Read these directions completely before using the product.
- All tubing and its connections to the BALLARD* 1000 HMEF should be properly attached and checked for leakage prior to use.
- It is recommended that the BALLARD* 1000 HMEF is changed every 24 hours or sooner if necessary.
- Before use, verify that the BALLARD* 1000 HMEF has no occlusions and that air will flow through it.

⚠ Warning

- The BALLARD* 1000 HMEF is designed for single use only and should not be cleaned and reused.
- Compensation of ventilation may be necessary when using the BALLARD* 1000 HMEF in applications close to the patient, since it will add 55 ml of dead space to the system.




Connection of the BALLARD* 1000 HMEF to Breathing Circuit

The BALLARD* 1000 HMEF should be placed between the proximal end of the artificial airway and the Y-piece of the breathing circuit.

Connectors ISO STD

Patient: 15mm F / 22mm M

Machine: 15mm M / 22mm F

 Single Use Only	Rx Only	Storage Temperature: -30°C/-22°F  +40°C/+104°F	 REF 157
---	---------	---	---