Instructions for Use: Tru-Flow™ Flowmeters

Not Made with Natural Rubber Latex

EN Indications for Use: *Used to measure Oxygen flow rate in liters per minute*

READ ALL INSTRUCTIONS PRIOR TO USING

Only trained individuals may install, use or maintain Flowmeters and oxygen fittings. Please read and understand the operating instructions before use of any oxygen equipment. Should you need assistance in understanding the function of this equipment, contact your authorized distributor for further information.

WARNINGS and CAUTIONS

-  DO NOT use oil, grease, organic lubricant or ANY combustible material around the flowmeter.
-  DO NOT autoclave.
-  DO NOT immerse in ANY liquid.
-  DO NOT gas sterilize (ETO).
-  DO NOT smoke in the area where oxygen is in use.
-  DO NOT use if inspection shows any visible damage or the ball float is not at the bottom of the flow tube pier to connecting to a gas source.

ENVIRONMENT

Oxygen is a non-flammable gas; it will, however, substantially increase the rate which materials burn. Oil grease and other hydrocarbons, when combined with oxygen, will become highly combustible and should never be used on any part of an oxygen flowmeter or other equipment used with oxygen.

Flowmeters are intended to be used with a 345 kPa (50 psi +/- 5 psi) gas source. NEVER attempt to attach a flowmeter directly to a cylinder.

For maximum accuracy the optimum operating temperature should be 70°F (21°C)

Flowmeters have 1/8" -18 NPT inlet port for connection to an adapter(s). Other flowmeters may be equipped with a DISS Oxygen nut or quick connect stem. Attaching accessories which are very restrictive to flow may affect accuracy.

OPERATION

1. Attach the flowmeter to a 345 kPa (50 psi +/- 5 psi) gas source. Ensure the flowmeter is in a vertical position.
2. Turn the control knob counter-clockwise to begin and increase flow.
3. Adjust to the desired flow setting. Take readings at the center of the ball.
4. Turn the control knob clockwise to decrease and stop flow. Do not over tighten the control knob.

NOTE: Always turn the control knob "OFF" (clockwise) when not using the flowmeter.

CONNECTIONS

Inlet: 1/8" NPT internal thread

Outlet: CGA DISS 1240 9/16-18 UNF-2A

MAINTENANCE AND REPAIR

 **Cleaning:** Flow tube and cover are made of polycarbonate material. Solvents and cleaning solutions may seriously affect the strength and transparency of the material. Use only cleaning solutions recommended for use with plastic. If cleaning of the unit is desired, the exterior surfaces can safely be wiped off with a clean, lint-free cloth and isopropyl alcohol.

Maintenance: This unit should be periodically checked and repaired, if necessary, or when malfunction is suspected or occurs.

Repairs: Repairs must be performed by a qualified service technician.

WARRANTY

Components of this system are warranted under use against all manufacturing defects from the date of purchase for 7 years. Any failure resulting from defective parts or faulty workmanship, as determined during evaluation by the manufacturer will be repaired or replaced under warranty. This warranty will be null and void for any unit that has been subjected to misuse, negligence, accident or repairs other than those performed by manufacturer.

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Manufactured for SunMed, LLC