# Instructions for Use: IGuard™ Eye Protector

**EN** **DEVICE DESCRIPTION:** The device is intended to cover and protect the patient's eyes against possible corneal abrasion, lacerations, and sclera hemorrhage.

Blue polyether / polyurethane foam, latex free synthetic rubber adhesive laminated to PE foam, release liner, PET eye cups. The product will be laminated together to encapsulate the eye cups between the foam and adhesive. The adhesive is a skin contact approved adhesive.

**INTENDED PURPOSE:** The I-Guard Eye Protector is used to cover and protect the patient's eyes.

**INDICATIONS:** The device is intended to cover and protect the patient's eyes against possible corneal abrasion, lacerations, and sclera hemorrhage during surgery.

**ENVIRONMENT:** Hospital, pre-hospital, surgical centers, medical clinics

**PATIENT TARGET GROUP:** Pediatric and Adult

**EXPECTED CLINICAL BENEFITS:** To shield patient's eyes during medical procedures.

**CONDRAINICATIONS:** None known

**WARNINGS:**

- External pressure on the eye protector may lead to adverse effects.
- Mild erythema may result in lengthy procedures.
- Do not resterilize. Resterilization could change the physical characteristics of the material and should not be attempted.
- Do not use on patient whose eyes do not fit within the openings of the lenses.
- This device is not designed to bear weight. Do not place objects on the eye protector.
- Head support suggested when patient is in prone or lateral position in order to prevent pressure around eyes.
- Frequently check to ensure there is no pressure on the patients' eyes during the use of the device and that the eyes remain within the opening of the lenses.
- Store at 40-80°F (4-27°C).
- Not made with natural rubber latex.

**INSTRUCTIONS:** Suggested procedure:

- Remove paper backing.
- Position over eyes.
- Ensure that patient's eyes are covered by lenses.
- Press edges to secure proper placement.
- Head support suggested when patient is in prone or lateral position in order to prevent pressure around eyes.
- Check to ensure there is no pressure on patient's eye during procedure.
- Remove slowly and gently.

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