
AVANOS*

MICROCUFF* ENDOTRACHEAL TUBE

Instructions for Use

AVANOS* MICROCUFF* Endotracheal Tube

Rx Only: Federal (USA) law restricts this device to sale by or on the order of a physician.

Description

AVANOS* MICROCUFF* Endotracheal Tubes are available cuffed in a variety of sizes and styles; refer to the package label for specific size and style. AVANOS* MICROCUFF* Endotracheal Tubes are designed with a Polyurethane High Volume Low Pressure Cuff and have an attached pilot balloon with a one-way Luer-slip adapter. AVANOS* MICROCUFF* Endotracheal Tubes are available in Magill and Oral Curved shapes, with (adult) and without (pediatric) a Murphy Eye. As a reference during intubation, AVANOS* MICROCUFF* Endotracheal Tubes have depth marks in centimeters, which indicate the distance to the distal tip. In addition, the endotracheal tube has a black intubation depth marking above the cuff for reference in determining the cuff and distal tip position after they have passed through the vocal cords and out of the field of vision. All AVANOS* MICROCUFF* Endotracheal Tubes have a radiopaque line.

Indications

AVANOS* MICROCUFF* Endotracheal Tubes are indicated for airway management by nasal or oral intubation of the trachea in pediatric and adult patients.

Contraindications

As with any PVC endotracheal tube, use of AVANOS* MICROCUFF* Endotracheal Tubes in procedures which will involve the use of a LASER beam or electrosurgical active electrode in the immediate area of the device is contraindicated. Contact of the endotracheal tube with a LASER beam or electrosurgical active electrode, especially in the presence of oxygen-enriched mixtures could result in rapid combustion of the endotracheal tube with harmful thermal effects and with emission of corrosive and toxic products, including hydrochloric acid (HCl).

Warnings

- Deflate cuff prior to repositioning the tube. Movement of the endotracheal tube with cuff inflated could result in patient injury or damage to the cuff, requiring a tube change.
- Do not over inflate cuff. Over inflation may result in airway damage.
- Do not reuse, reprocess, or resterilize this medical device. Reuse, reprocessing, or resterilization may 1) adversely affect the known biocompatibility characteristics of the device, 2) compromise the structural integrity of the device, 3) lead to the device not performing as intended, or 4) create a risk of contamination and cause the transmission of infectious diseases resulting in patient injury, illness, or death.

Cautions

General

- Do not use if package has been previously opened or damaged.
- Exposure to elevated temperatures and ultraviolet light should be avoided during storage.
- Intubation and extubation should be performed using currently accepted medical techniques.
- Each tube's cuff, pilot balloon, and valve should be tested by inflation before use.
- Avoid damaging the cuff during intubation. If the cuff is damaged, the tube should not be used.
- Inflation of cuff by "feel" or by using a measured amount of air is not recommended since resistance is an unreliable guide during inflation.
- In selecting the sealing pressure, an intracuff pressure measuring device should be used in conjunction with minimal Occluding Volume or Minimum Leak techniques. Cuff pressure should be monitored. Any deviation from the selected seal pressure should be investigated and corrected immediately.
- Should extreme chin-to-chest flexing of the head or movement of the patient (e.g. to a lateral or a prone position) be anticipated after intubation, use of a reinforced endotracheal tube should be considered.
- When the patient's position is altered after intubation, it is essential to verify that the tube position remains correct in the new patient position.
- Tubes should be securely anchored to avoid unnecessary tube movement.
- Seat the connector firmly in both the endotracheal tube and adapter on the ventilation equipment to prevent disconnection during use.
- A bite block should be used in cases where the patient may bite down and flatten the endotracheal tube.
- Non-standard dimensioning of some connectors on ventilators or anesthesia equipment may make secure mating with the endotracheal tube's 15 mm connector difficult.
- Good clinical judgment should be used when selecting complementary devices, such as bronchoscopes and tracheal suction catheters, which are intended to be inserted through the endotracheal tube. One clinical review recommends using a suction catheter with outer diameter no greater than 70% of the inner diameter of the endotracheal tube.³
- The inflation line check valve may interfere with Magnetic Resonance Imaging (MRI) picture clarity. Ensure the valve is positioned away from the area being scanned.
- Diffusion of nitrous oxide mixture, oxygen, or air may either increase or decrease cuff volume and pressure. To decrease such diffusion, inflating the cuff with the same gas mixture that will contact the cuff's external surface is recommended.
- Diffusion of water molecules through the polyurethane cuff may lead to the formation of condensate within the inflation system. Condensation in the cuff does not impact product performance; however, condensation in the pilot balloon may lead to inaccurate cuff pressure readings when measured with a cuff pressure manometer. Minimal Occluding Volume or Minimal Leak Techniques are not affected by condensation and may be used to maintain an adequate seal.
- Three-way stopcocks or other devices should not be left inserted in the inflation valve for extended periods of time. The resulting stress could crack the valve housing and allow the cuff to deflate.
- The use of lidocaine topical aerosol has been associated with the formation of pinholes in PVC cuffs (Jayasuriya, K.D., and Watson, W.F.: P.V.C. Cuffs and Lignocaine-based Aerosol, Brit. J. Ann. 53:1368, 1981). The same authors report that lidocaine hydrochloride solution does not have this effect. Lab tests with polyurethane cuffs have shown to produce similar negative interactions as PVC cuffs with lidocaine topical aerosol. Also similar to PVC cuffs, lab testing showed polyurethane cuffs did not have these negative interactions with lidocaine hydrochloride solution.
- Follow the manufacturer's application instructions when using lubricating jellies with endotracheal tubes. If excessive amounts of jelly dry onto the inner surface of the endotracheal tube, it could result in either a lubricant plug or a clear film that partially or totally blocks the airway.

- Do not resterilize.
- This medical device contains DEHP (diethylhexylphthalate) that is presently classified in the European Union as a presumed human reproductive toxicant based on data from animal studies. There is no conclusive scientific evidence that exposure to DEHP contained in medical devices has caused harmful effects in humans. A risk assessment, which took into account DEHP exposure of all indicated patient populations, including those who are potentially at increased risk, was performed for this device and the conclusion is that the device is safe when used as directed.

Adverse Reactions

Reported adverse reactions associated with the use of tracheal tubes are many and diverse. Standard textbooks and the scientific literatures should be consulted for specific adverse reaction information.

The order of listing is alphabetical and does not indicate frequency or severity. Reported adverse reactions include: abrasion of the arytenoid cartilage vocal process; cartilage necrosis; cicatrix formation; consequences of failure to ventilate including death; damage to the perichondrium; development of dense or diffuse fibrosis invading the entire glottic area; emphysema; endobronchial aspiration; endobronchial intubation (hypoxemia); endotracheobronchial aspiration; epistaxis; esophageal intubation (stomach distention); excoriated membranes of the pharynx; eye trauma: fibrin deposition; formation of subglottic web; fracture-luxation of cervical column (spinal injury); fragmentation of cartilage; glottic edema (supraglottic, subglottic, retroarytenoid); granuloma of the inner arytenoid area; infections (laryngitis, sinusitis, abscess, respiratory tract infection); inflammation; intermittent aphonia and recurrent sore throat; laryngeal fibrosis; laryngeal granulomas and polyps; laryngeal obstruction; laryngeal stenosis; laryngeal ulcers; laryngotracheal membranes and webs; membranous glottic congestion; membranous tracheobronchitis; mild edema of the epiglottis; mucosal sloughing; paresis of the hypoglossal and/or lingual nerves; perforation of esophagus; perforation of the trachea; pneumothorax; replacement of the trachea wall with scar tissue; respiratory obstruction; retrobulbar hemorrhage; retropharyngeal abscess; retropharyngeal dissection, rupture of the trachea; sore throat, dysphagia; stricture of nostril; stridor; subglottic annular cicatricial stenosis; submucosal hemorrhage, submucous puncture of the larynx; superficial epithelial abrasion; swallowed tube; symecia of the vocal cords; teeth trauma; tissue burns; tracheal bleeding; tracheal stenosis; trauma to lips, tongue, pharynx, nose, trachea, glottis, palate, tonsil, etc.; traumatic lesions of the larynx and trachea; ulcerations exposing cartilaginous rings and minor erosions at cuff site; ulceration of the lips, mouth, pharynx; ulcers of the arytenoid; vocal cord congestion; vocal cord paralysis, and vocal cord ulcerations; barotrauma is a potential complication of the use of this catheter for oxygen delivery, especially in patients with compromised lung conditions (ie, previous pneumothoraces, pulmonary fibrosis, COPD, previous lung resection/lobectomies, etc).

Helpful references for more detailed discussions of tracheal tube adverse reactions include the following:

- Baron, S.H., and Kohlmoos, H.W.: Laryngeal sequelae of endotracheal anesthesia. Ann. Otol. 60:767-792, 1951.
- Blanc, V.F., and Tremblay, N.A.G.: The complications of tracheal intubation. Anest. Analg. 53:202-213, 1974.
- Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82:981-986, 1968.
- Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann. Surg. 169:334-348, 1969.
- Gross, C.W., and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
- McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80:556-564, 1971.

Instructions for Use: Use aseptic technique

Intubation and extubation should be performed following currently accepted medical techniques.

Adult Use:

- Expert clinical judgment should be used in choosing the suitable endotracheal tube size and style for each adult patient.
- Remove the sterile AVANOS* MICROCUFF* Endotracheal Tube from its protective package.
- Test the cuff, pilot balloon and valve of each tube by inflation/deflation prior to use. Insert a Luer-tip syringe into the cuff inflation valve housing and inflate and completely deflate the cuff to check appropriate expansion and deflation.
- In general the tube should not be altered by cutting. If precutting of the tracheal tube is desired by the physician, evaluate suitability of the tube for precutting prior to intubation. If the tube is pre-cut, it should be cut at a slight angle to facilitate reinsertion of the 15 mm connector into the tube. **Always ensure the connector is firmly seated in both the tracheal tube and the breathing circuit to prevent disconnection during use.**
- In those situations where it is deemed appropriate to pre-cut the tube, the user is cautioned that anatomical variations, conditions of use or other factors may result in a tracheal tube either too long or too short for a given patient when the standard pre-cut length indicated on the tube is used. Expert clinical judgment should be used in selecting the appropriate tube size and pre-cut length.
- Intubate the patient following currently accepted medical techniques with consideration given to the specific cuff-related WARNINGS and CAUTIONS stated in this product insert.
- Once the patient is intubated, inflate the cuff only with enough gas mixture to provide an effective seal at the desired lung inflation pressure. The use of Minimal Occluding Volume, Minimum Leak techniques and monitoring (measuring) of cuff pressure, can help reduce the occurrence of many of the adverse reactions associated with the use of cuffed tracheal tubes.
- Remove syringe from the valve housing after cuff inflation. Leaving the syringe attached will keep the valve open, permitting the cuff to deflate. Firmly secure and avoid moving the endotracheal tube. Ensure patency of the cuff inflation line and avoid kinking/obstruction of the endotracheal tube.
- Check to verify inflation system is not leaking. Integrity of the system should be verified periodically during the intubation period. Uncorrected failure of the inflation system could result in death. Cuff pressure should be closely monitored and any deviation from the selected seal pressure should be investigated and corrected immediately.
- Prior to extubation, deflate cuff by inserting syringe into valve housing and removing gas mixture until a definite vacuum is noted in the syringe and the pilot balloon is collapsed.
- Extubate the patient following currently accepted medical techniques.
- Discard the tracheal tube.

Pediatric Use: Follow steps 1 – 12 above and pay special attention to the following recommendations for safe and effective use of the pediatric endotracheal tube.

- To assist in the proper selection of the AVANOS* MICROCUFF* Endotracheal Tube, Pediatric Oral/Nasal Magill and Pediatric Oral Curved, refer to the sizing chart below.^{1,2} The chart below may not apply to a child known to have abnormal tracheal anatomy or dimensions.
- During the intubation process pay careful attention to ensure that the black intubation depth mark is placed between the vocal cords. Auscultate for symmetrical breath sounds to avoid endobronchial intubation.
- After intubation, an audible air leak must be present at ≤ 20 cm H₂O airway pressure with the cuff fully deflated. If no air leak is detected, the tube may be too large; consider changing the tube.
- Inflate the cuff to the effective sealing pressure but not higher than 20 cm H₂O cuff pressure. In cases of excessive airway leakage, re-evaluate patient and intubation depth.
- Cuff pressure in pediatric patients should be monitored continuously. Hyperinflation can be minimized by use of an automatic pressure regulator or pressure relief valve² and adjusted if necessary to maintain an appropriate cuff pressure (≤ 20 cm H₂O) throughout the procedure. **Note: Manual compression of the pilot balloon should be avoided.**
- The pediatric cuff should not be fully deflated except at the time of air leak testing and prior to extubation.

Recommended Age or Weight	Tube Interior Diameter
Term ≥ 3 kg up to < 8 months	ID 3.0 mm
8 months to < 2 years	ID 3.5 mm
2 to < 4 years	ID 4.0 mm
4 to < 6 years	ID 4.5 mm
6 to < 8 years	ID 5.0 mm
8 to < 10 years	ID 5.5 mm
10 to < 12 years	ID 6.0 mm
12 to < 14 years	ID 6.5 mm
14 to < 16 years	ID 7.0 mm

¹Dullenkopf, A., A.C. Gerber and M. Weiss, Fit and seal characteristics of a new paediatric tracheal tube with high volume-low pressure polyurethane cuff. Acta Anaesthesiol Scand 2005; 49:232-237.

²Salgo, B., et al. Evaluation of a new recommendation for improved cuffed tracheal tube size selection in infants and small children. Acta Anaesthesiol Scand 2006; 50:557-561.

³Morrow BM, Argent AC, A Comprehensive Review of Pediatric Endotracheal Suctioning: Effects, Indications, and Clinical Practice. Pediatric Critical Care Medicine 2008, 9(5):465-477.

I.D. Inside diameter	O.D. Outside diameter	Single Use Only	Radiopaque	Sterilized using Ethylene Oxide	Do not use if package is damaged	Do not resterilize
Contains PHT DEHP	Rx Only	Caution	Consult instructions for use	Protect from heat and radioactive sources	Batch code	Use by

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1. Baron, S.H., and Kohlmoss, H.W.: Laryngeal sequelae of endotracheal anesthesia. Ann. Otol. 60:767-792, 1951.
2. Blanc, V.F., and Tremblay, N.A.G.: The complications of tracheal intubation. Anest.Anal. 53:202-213, 1974.
3. Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82:981-986, 1968.
4. Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann. Surg. 169:334-348, 1969.
5. Gross, C.W., and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
6. McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80:556-564, 1971.

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1. Baron, S.H., and Kohlmoo, H.W.: Laryngeal sequelae of endotracheal anesthesia. Ann. Otol. 60:767-792, 1951.
2. Blanc, V.F., and Tremblay, N.A.G.: The complications of tracheal intubation. Anest.Analg. 53:202-213, 1974.
3. Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82:981-986, 1968.
4. Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann.Surg. 169:334-348, 1969.
5. Gross, C.W. and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
6. McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80:556-564, 1971.

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³Morrow BM, Argent AC. A Comprehensive Review of Pediatric Endotracheal Suctioning: Effects, Indications, and Clinical Practice. Pediatric Critical Care Medicine 2008; 9(5):465-477.

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1. Baron, S.H., and Kohlmoos, H.W.: Laryngeal sequelae of endotracheal anesthesia. Ann. Otol. 60:767-792, 1951.
2. Blanc, V.F., and Tremblay, N.A.G.: The complications of tracheal intubation. Anesth. Analg. 53:202-213, 1974.
3. Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82:981-986, 1968.
4. Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann. Surg. 169:334-348, 1969.
5. Gross, C.W., and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
6. McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80:556-564, 1971.

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- Baron, S.H., and Kohlmoos, H.W.: Laryngeal sequelae of endotracheal anesthesia. Ann. Otol. 60:767-792, 1951.
- Blanc, V.F., and Tremblay, N.A.G.: The complications of tracheal intubation. Anest. Analg. 53:202-213, 1974.
- Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82:981-986, 1968.
- Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann. Surg. 169:334-348, 1969.
- Gross, C.W. and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
- McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80:556-564, 1971.

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¹Dullenkopf, A., A.C. Gerber and M. Weiss, Fit and seal characteristics of a new paediatric tracheal tube with high volume–low pressure polyurethane cuff. Acta Anaesthesiol Scand 2005;49:232-237.

²Salgob, B., et al. Evaluation of a new recommendation for improved cuffed tracheal tube size selection in infants and small children. Acta Anaesthesiol Scand 2006;50:557-561.

³Morrow BM, Argent AC. A Comprehensive Review of Pediatric Endotracheal Suctioning: Effects, Indications, and Clinical Practice. Pediatric Critical Care Medicine 2008, 9(5):465-477.

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- Baron, S.H., and Kohlmoos, H.W.: Laryngeal sequelae of endotracheal anesthesia. Ann. Otol. 60-767-92, 1951.
- Blanc, V.F., and Tremblay, N.A.G.: The complications of tracheal intubation. Anest.Analg. 53:202-213, 1974.
- Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82-981-986, 1968.
- Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann. Surg. 169-334-348, 1969.
- Gross, C.W. and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
- McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80:556-564, 1971.

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- Baron, S.H., and Kohlmoo, H.W.: Laryngeal sequelae of endotracheal anesthesia. Ann. Otol. 60:767-792, 1951.
- Blanc, V.F., and Tremblay, N.A.G.: The complications of tracheal intubation. Anest.Analg. 53:202-213, 1974.
- Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82:981-986, 1968.
- Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann. Surg. 169:334-348, 1969.
- Gross, C.W. and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
- McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80:556-564, 1971.

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1. Baron, S.H., and Kohlmoos, H.W.: Laryngeal sequelae of endotracheal anesthesia. Ann. Otol. 60-767-792, 1951.
2. Blanc, V.F., and Tremblay, N.A.G.: The complications of tracheal intubation. Anest. Analg. 53:202-213, 1974.
3. Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82-981-986, 1968.
4. Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann. Surg. 169-334-348, 1969.
5. Gross, C.W., and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
6. McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80:556-564, 1971.

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1. Baron, S.H. in Kohlmoss, H.W.: Laryngeal sequelae of endotracheal anesthesia. Ann. Otol. 60:767-792, 1951.
2. Blanc, V.F. in Tremblay, N.A.G.: The complications of tracheal intubation. Anest.Analg. 53:202-213, 1974.
3. Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82:981-986, 1968.
4. Cooper, J.D. in Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann. Surg. 169:334-348, 1969.
5. Gross, C.W. in Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
6. McGovern, F.H., Fitz-Hugh, G.S. in Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80:556-564, 1971.

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1. Baron, S.H., and Kohlmoss, H.W.: Laryngeal sequelae of endotracheal anesthesia. Ann. Otol. 60:767-792, 1951.
2. Blanc, V.F., and Tremblay, N.A.G.: The complications of tracheal intubation. Anest.Analg. 53:202-213, 1974.
3. Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82:981-986, 1968.
4. Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann. Surg. 169:334-348, 1969.
5. Gross, C.W. and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
6. McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80:556-564, 1971.

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- Baron, S.H., and Kohlmoos, H.W.: Laryngeal sequelae of endotracheal anesthesia. Ann. Otol. 60-767-792, 1951.
- Blanc, V.F., and Tremblay, N.A.G.: The complications of tracheal intubation. Anest. Analg. 53:202-213, 1974.
- Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82-981-986, 1968.
- Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann. Surg. 169-334-348, 1969.
- Gross, C.W., and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
- McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80:556-564, 1971.

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- Baron, S.H., and Kohlmoos, H.W.: Laryngeal sequelae of endotracheal anesthesia. Ann. Otol. 60:767-792, 1951.
- Blanc, V.F., and Tremblay, N.A.G.: The complications of tracheal intubation. Anesth.Analg. 53:202-213, 1974.
- Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82:981-986, 1968.
- Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann. Surg. 169:334-348, 1969.
- Gross, C.W., and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
- McGovern, F.H., Fitz-Hugh, G.S., and Edgemon, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80:556-564, 1971.

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¹Dullenkopf, A., A.C. Gerber and M. Weiss, Fit and seal characteristics of a new paediatric tracheal tube with high volume-low pressure polyurethane cuff. Acta Anaesthesiol Scand 2005; 49:232-237.
²Salgo, B., et al. Evaluation of a new recommendation for improved cuffed tracheal tube size selection in infants and small children. Acta Anaesthesiol Scand 2006; 50:557-561.
³Morrow BM, Argent AC, A Comprehensive Review of Pediatric Endotracheal Suctioning: Effects, Indications, and Clinical Practice. Pediatric Critical Care Medicine 2008, 9(5):465-477.

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- Baron, S.H., and Kohlmoos, H.W.: Laryngeal sequelae of endotracheal anesthesia. Ann. Otol. 60:767-792, 1951.
- Blanc, V.F., and Tremblay, N.A.G.: The complications of tracheal intubation. Anesth. Analg. 53:202-213, 1974.
- Campbell, D.: Trauma to larynx and trachea following intubation and tracheostomy. J. Laryng. 82:981-986, 1968.
- Cooper, J.D., and Grillo, H.C.: The evolution of tracheal injury due to ventilatory assistance through cuffed tubes. Ann. Surg. 169:334-348, 1969.
- Gross, C.W., and Gros, J.C.: Rare complications after prolonged translaryngotracheal intubation. Ann. Otol. 80:582-584, 1971.
- McGovern, F.H., Fitz-Hugh, G.S., and Edgeman, L.J.: The hazards of endotracheal intubation. Ann. Otol. 80:556-564, 1971.

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¹Düllenkopf, A., A.C. Gerber and M. Weiss, Fit and seal characteristics of a new paediatric tracheal tube with high volume-low pressure polyurethane cuff. Acta Anaesthesiol Scand 2005; 49:232-237.

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³Morrow BM, Argent AC, A Comprehensive Review of Pediatric Endotracheal Suctioning: Effects, Indications, and Clinical Practice. Pediatric Critical Care Medicine 2008; 9(5):465-477.

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Distributed in the USA by Avanos Medical Sales, LLC, 5405 Windward Parkway, Alpharetta, GA 30004 USA.

In USA, 1-844-4AVANOS (1-844-428-2667). www.avanos.com

Avanos Medical, Inc., 5405 Windward Parkway, Alpharetta, GA 30004 USA.

Avanos Medical Belgium BVBA, Leonardo Da Vinciiaan 1, 1930 Zaventem, Belgium.

Sponsored in Australia by Avanos Medical Australia Pty Ltd, 475 Victoria Avenue, Chatswood, NSW 2067, Australia.

CE 2797

*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.
© 2018 AVNS. All rights reserved. 2019-07-01 15-M1-228-01