Care and Service of Reusable GreenLine® Fiber Optic Laryngoscope Handles with Removable Lamp Cartridge

Not Made with Natural Rubber Latex

R_X ONLY
USA Law Only

CLEAN HANDLE IMMEDIATELY AFTER USE TO PREVENT RESIDUES FROM DRYING

Intended Use: Used for endotracheal tube intubation.

Care and Service Document: These instructions apply to SunMed's reusable GreenLine® fiber optic laryngoscope handles with removable lamp cartridge. SunMed has validated these directions as being capable of preparing its reusable GreenLine® handles for re-use. Fiber optic laryngoscope handles and lamp cartridges must be processed prior to first use, and between each use, using any of the methods presented in this care and service document:

- CLEANING + (A) COLD SOLUTION HIGH-LEVEL DISINFECTION
- CLEANING + (B) GRAVITY AUTOCLAVE STERILIZATION
- CLEANING + (C) PRE-VACUUM STEAM STERILIZATION
- CLEANING + (D) STERIS AMSCO V-PRO STERILIZATION
- CLEANING + (E) STERRAD STERILIZATION

WARNINGS and CAUTIONS

- ⚠ DO NOT USE THE FOLLOWING: Ultrasonic clean, dry heat sterilization or Potassium Hydroxide liquids.
- ⚠ REMOVE LAMP CARTRIDGE AND BATTERIES BEFORE cleaning, high-level disinfecting, or sterilizing laryngoscope handles.
- ⚠ CHECK laryngoscope system functionality, and ensure blade is fully engaged on handle prior to use.
- ⚠ DISCARD any component that shows evidence of damage or deterioration.

NOTE: Batteries, battery cartridge, and the lamp cartridge **MUST** be removed prior to cleaning, disinfection, and/or sterilization of the laryngoscope handle. See instructions below for disinfection of the the lamp cartridge.

HANDLE DISASSEMBLY

Battery Removal To remove batteries, unscrew handle end cap and remove battery(ies) (Fig.1, 2). Batteries in the Stubby Handle (Fig.1) sit in a battery cartridge, which is removeable. Follow with lamp cartridge removal.

IMPORTANT: Lamp Cartridge Removal

Allow the bulb time to cool. The bulb is permanently housed in the lamp cartridge. To remove lamp cartridge, unscrew handle head from the handle body **and unscrew the lamp cartridge by hand. No tooling is required** (Fig. 3). Follow with lamp cartridge disinfection.

LAMP CARTRIDGE DISINFECTION

The lamp cartridge cannot use high-level disinfection or sterilization for preparation for reuse. Instead, it can be disinfected with an intermediate-level disinfection procedure.

First, select a suitable quaternary ammonium and isopropanol alcohol based germicidal wipe that is designed for intermediate level disinfection. Refer to the manufacturer's instructions. To disinfect:

- Following removal of the lamp cartridge (Fig. 3), use a germicidal cloth to wipe the bulb to remove soil.
- With a second cleangermicidal cloth, thoroughly wipe all surfaces of the lamp cartridge, ensuring it remains wet for two minutes. Use additional cloths as needed to keep the cartridge surface wet for the entire contact time.
- Allow the lamp cartridge to air dry completely.
- Set lamp cartridge aside for re-installation following disinfection of laryngoscope handle.

Follow with laryngoscope handle disinfection.


Fig. 1 Stubby Handle
Fig. 2 All Other Handle Sizes


Fig. 3

CLEANING THE HANDLE PRIOR TO DISINFECTION / STERILIZATION

- Remove laryngoscope blade from the handle.
- Disassemble the handle, and remove the lamp cartridge and batteries (Fig. 1, 2, 3). Follow the separate lamp cartridge disinfection procedure.
- Rinse the disassembled handle body, head and end cap under cool running tap water, until all visible soil is removed.
- Immerse the handle body, head and end cap in a presoak enzymatic prepared according to manufacturer's instructions, for a minimum of two minutes.
- After the presoak, scrub the handle body, head and end cap using a soft bristled brush, and rinse with lukewarm running tap water for a minimum of one minute to remove any residues and visible soil.
- Immerse the handle body, head and end cap in an enzymatic detergent prepared according to manufacturer's instructions.
- Scrub the handle body, head and end cap thoroughly by using a soft bristle brush, and thoroughly rinse with cool running tap water to remove all residual detergent.
- Dry the handle body, head and end cap with lint-free clean cloth or filtered pressurized air.
- Reassemble the handle body, head and end cap without batteries or lamp cartridge.

Follow with HIGH-LEVEL DISINFECTION or STERILIZATION of handle.

AUTOMATED CLEANING FOR HANDLES:

- Remove laryngoscope blade from the handle
- Disassemble the handle, and remove the lamp cartridge and batteries (Fig 1, 2, 3). Follow the separate lamp cartridge disinfection procedure.
- Place the handle and bottom cap into an Automated washer per ISO15883 and follow the below recommended cycle and detergent. Equivalent washers, cycles, and detergent can be utilized.
- Reassemble the handle body, head, and end cap without batteries or lamp cartridge.
- Follow with HIGH-LEVEL DISINFECTION or STERILIZATION of handle.

Phase	Recirculation time (Minutes)	Temperature	Detergent Type and Concentration
Pre-Wash	02:00	Cold Tap Water	N/A
Wash	02:00	43°C Tap Water	Valsure Enzymatic 2 mL/L
Rinse	01:00	43°C Tap Water	N/A
Dry Time	07:00	90°C	N/A

(A) COLD SOLUTION HIGH-LEVEL DISINFECTION OF THE HANDLE

- Prepare Cidex® OPA or 2.4% Glutaraldehyde solution according to manufacturer's instructions.
- Immerse the assembled handle body, head and end cap in the disinfectant.
- Thoroughly rinse by immersing the assembled handle in purified water. The handle may need to be disassembled to ensure thorough rinsing.
- Dry with a lint-free clean cloth or filtered pressurized air. The handle may need to be disassembled to ensure thorough drying.
- If necessary, reassemble the handle body, head and end cap.

Follow with testing prior to use.

STERILIZATION OF THE HANDLE

Use one of the following four sterilization methods to sterilize SunMed's reusable GreenLine® fiber optic laryngoscope handles. **IMPORTANT: Prior to sterilization, batteries, battery cartridge and lamp cartridge must be removed.**

Steam Sterilization Insert assembled handle (without batteries, battery cartridge or lamp cartridge) in an appropriate autoclave pouch.

	(B) GRAVITY DISPLACEMENT STEAM	(C) PRE-VACUUM STEAM
Temperature	121° C (250° F)	132° C (270° F)
Cycle Time	30 Minutes	4 Minutes
Dry Time	15 Minutes	20 Minutes

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(E) STERRAD® STERILIZATION SYS.

- STERRAD® 100NX® System (Std and EXPRESS Cycle)
- STERRAD® NX® System (Std Cycle)
- STERRAD® 100S and 200 System (Short Cycle Outside US)
- STERRAD® 50 System

Follow with battery and lamp cartridge installation.

BATTERY & LAMP CARTRIDGE INSTALLATION

HANDLE	BATTERY TYPE (QTY)
Large	1.5 v / D (2)
Stubby	1.5 v / AAA (2)
Medium	1.5 v / C (2)
Penlite	1.5 v / AA (2)
Mini	3 v / Lithium 123 (1)
Micro	3 v / Lithium 123 (1)

Battery Installation

Following handle cleaning, disinfection or sterilization, remove end cap (Fig.1, 2) and insert batteries (and battery cartridge if applicable). See Battery Type chart for accurate battery replacement.

IMPORTANT: Lamp Cartridge Installation

Following handle cleaning, disinfection or sterilization, remove handle head (Fig. 3) **and insert lamp cartridge by hand tightening No tooling is required.**

LAMP TYPE	CARTRIDGE BASE COLOR	RE-ORDER #
High End LED	Yellow	5-0240-63
Eco LED	Orange	5-0240-62
Halogen Xenon	Blue	5-0240-61

See video at: www.Sun-Med.com/resources for assembly method and review Lamp Type chart for accurate cartridge replacement if lamp no longer illuminates.

Follow with testing prior to use.

PRIOR TO USE

Laryngoscope handles must be tested after each reprocessing cycle and prior to use. To test:

- Fully engage compatible blade on laryngoscope handle. If the light guide fails to light, check the batteries and lamp for correct placement. Replace with new batteries or lamp if needed.
- If light guide fails to light following additional blade engagements, do not use laryngoscope.

SunMed, LLC,
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA

Manufactured for SunMed, LLC,
a Subsidiary of SunMed Holdings, LLC
www.Sun-Med.com
Made in Pakistan

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