# Vital Signs™ Pressure Infusors/ Vital Signs Clear Pressure Infusors

<b>REF</b>	Catalogue number		Manufacturer
<b>LOT</b>	Lot Number		Manufacture date
<b>X</b>	Quantity		Do not re-use
	Not made with natural rubber latex	<b>SIZE</b>	Volume
	Consult Instructions for Use vyaire.eifu.com	<b>DATE</b>	Date
<b>EC REP</b>	Authorized Representative in the European Community		Patient
<b>MD</b>	Medical Device		Solution
<b>Rx ONLY</b>	U.S. Federal law restricts this device to sales by or on the order of a physician		Exp. Date
<b>en</b>	<b>English</b>		

## Intended Use:

This pressure infusor is a device consisting of an inflatable cuff which is placed around an I.V. bag. When the device is inflated, it increases the pressure on the I.V. bag to assist the infusion of the fluid.

## Direction For Use:

1. Insert fluid bag between netting or clear window and air bladder. Hang fluid bag on fluid bag hook. Use the loop to hang pressure infusor from the IV pole.
2. Turn stopcock handle marked "OFF" to point toward the open stopcock vent. See diagram 1.
3. Pump inflation bulb until the pressure gauge indicates the desired pressure. Do not exceed 300 mmHg of pressure.
4. **IMPORTANT:** To maintain a leak proof seal, turn the stopcock handle marked "OFF" toward the infusor bag. See diagram 2.
5. To deflate, turn stopcock handle marked "OFF" toward the inflation bulb. See diagram 3.
6. Dispose of the pressure infusor in accordance with hospital policies.

## WARNINGS:

1. Prior to infusion it is recommended that the administration set be properly primed and residual air removed from the fluid bag and tubing to minimize the risk of air embolism.

**Vyair Medical, Inc.**  
26125 North Riverwoods Blvd  
Mettawa, IL 60045 USA  
1-833-327-3284  
customersupport@vyaire.com  
www.vyaire.com

**EC REP** **EMERGO EUROPE**   
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

© 2020 Vyair. Vyair, the Vyair Logo and all other trademarks or registered trademarks are property of Vyair Medical, Inc., or one of its affiliates.

**vyaire™**


2. Make sure pressure gauge can move freely and is clear of all obstructions.
3. Check the infusor pressure periodically during use.
4. Do not connect intravascular fluid systems due to the luer port on the pressure infusor. This may result in an air embolism which could lead to serious patient injury.

Dispose per local and state-environmental regulations.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer at GMB-RS-Complaint-Intake@vyaire.com and the competent authority of the Member State in which the user and/or patient is established.

- This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.

A discontinued use of this product on the same patient, with periods of use separated by a period of non-use, is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.

## Note:

1. Do not use if product is damaged.
2. Do not exceed the pressure infusor capacity.
3. For expiration date, reference product labeling.

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