Gas Sampling Lines, M/M Luer

EN **DEVICE DESCRIPTION:** The device is non-sterile, disposable, single-patient use sampling line that enables sampling of exhaled gases. The sampling line available in different end connectors. **INTENDED PURPOSE:** The Sample Lines are intended to be used where exhaled gas is monitored. **INDICATIONS:** The Sample Lines are intended to be used to sample exhaled gas in intubated or non-intubated patients requiring expired gas monitoring. **ENVIRONMENT:** Hospital, sub-acute, pre-hospital **PATIENT TARGET GROUP:** Infant, pediatric, adult **EXPECTED CLINICAL BENEFITS:**

- Notifies clinicians of changes in the patient's ventilator status
- Provides information on the effectiveness of cardiac resusciation
- Provides verification of placement of ETT in trachea.

INSTRUCTIONS:

- VERIFY universal luer connector on monitor and circuit.
- REMOVE sampling line from package.
- ATTACH line to luer connection of monitor securely.
- ATTACH other end of line to luer connection on circuit.
- TEST for air fl ow and leaks prior to use.

CONTRAINDICATIONS: None known

WARNINGS: Position tubing to avoid strangulation. Do not soak, rinse, wash, gas or steam sterilize these products. To do so may result in harmful, residual gases or erroneous measurements. Do not use in the presence of flammable anesthetics. To reduce the risk of misconnections and patient injury, always trace tubing from gas source to the medical device before connecting. **CAUTIONS:** Consult physician if patient develops infection, skin irritation or material sensitization. Do not place anything on supply tubing that may obstruct flow. Excessive patient secretion or a build-up of liquids in the tubing may occlude the sampling and require more frequent replacement.

INSTRUCTIONS: Contact your country's Competent Authority to report any serious incident. Dispose of device in accordance with local, state or national regulations. Follow instructions for proper care of device to prevent risk of contamination and injury. This product is disposable and is not intended for prolonged use. Check for leaks at connections and for proper placement prior to use. Check for proper position to ensure detection of exhaled gas. Replace sampling line when blockage is indicated by the capnography monitor.

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