# REF 377CBR


## AirLife DuoTherm™ Humidification Chamber

Compatible with AirLife DuoTherm Series  
Heated Humidifier

<b>REF</b>	Catalogue number		Manufacturer
<b>LOT</b>	Lot Number		MR Unsafe
<b>QTY</b>	Quantity		Do not re-use Single Use
	Consult Instructions for Use		Maximum Water Level
	Caution		Not made with natural rubber latex
	U.S. Federal law restricts this device to sales by or on the order of a physician		Expiration Date
	Not made with Phthalates		Keep away from sunlight
	Storage Temperature		Warning
	Humidity limit		

**AU** **REP** **Australian Sponsor:**  
**AirLife Australia Holdings Pty Ltd**  
PO Box 97  
North Ryde BC, NSW, 1670  
Australia

 **Manufactured By**  
**AirLife**  
2710 Northridge Dr. NW, Suite A  
Grand Rapids, MI 49544 USA  
www.myAirLife.com

EN

English

### Intended Use

The AirLife DuoTherm Humidification System is intended to add moisture and warmth to breathing gases administered to patients that require assistance breathing or mucosal humidification. Gases that are available for medical use do not contain sufficient moisture and may damage or irritate the respiratory tract.

**Note:** The AirLife DuoTherm Humidification Chamber is compatible with the 377HTR AirLife DuoTherm Series Heated Humidifier and AirLife DuoTherm 377HTR Series Circuits.

### Indication For Use

The AirLife DuoTherm Humidification Chamber is intended to hold water required to humidify breathing gases delivered to patients ranging from neonates to adults using a heated humidifier. The product is a single use device, non-sterile, and used in professional healthcare environments under the supervision of a licensed healthcare practitioner.

### Warning

- DO NOT use the chamber if water level is above the maximum water level line.
- DO NOT use the chamber if the blue caddy is not intact when received.
- DO NOT operate the chamber at an angle in excess of 10 degrees.
- DO NOT fill the chamber with water in excess of 37°C.
- DO NOT touch the heater plate or chamber base. Temperatures may exceed 85°C.
- DO NOT turn off water level alarms.
- DO NOT soak, rinse, wash, or sterilize.

- Ensure there is a water supply connected to the chamber and that water is present within the chamber.
- Set appropriate ventilator alarms.
- For Single Patient Use Only. Reuse may degrade performance and result in transmission of infectious substances, interruption of treatment, serious harm or death.
- Always position the humidifier with humidification chamber lower than the patient.
- This product is for use only by trained medical professionals under a doctor's supervision.
- Use USP Sterile Water for Inhalation only.
- High Flow Nasal Cannula (HFNC) therapies require the use of a ventilator or pressure limiting device such as a relief valve (AH115-05 or equivalent) to prevent damage to the system in the event of an occlusion.
- DO NOT use the chamber if it has been dropped or damaged.
- DO NOT use the chamber for more than 14 days.
- DO NOT add any attachments or accessories to humidifier that are not listed in the instruction for use of humidifier. Otherwise, the accessory or humidifier may adversely affect device performance, the quality of the therapy or injuring the patient.
- Use above the maximum flow rates and/or operating pressure may affect the quality of therapy and lead to damage resulting in leakage and loss of ventilation pressure.

### Caution

- Chamber may only be used with 377HTR AirLife DuoTherm Series Humidifier and AirLife DuoTherm 377HTR Series Circuits.
- Operation or storage under direct sunlight should be avoided.
- DO NOT administer nebulizer through the heated humidifier, humidification chamber, and circuit. Administering nebulizer may adversely affect device performance.
- The heated humidifier, humidification chamber, and circuit are not intended to be used with inhalers and/

or for administering drugs, medicinal substances or human blood derivatives.

- Do not manually fill the humidification chamber.
- The medical professional is responsible for the compatibility of the humidifier and accessories used to connect to the patient, ventilator and other equipment before use.
- The humidification chamber may have hot fluid temperature that exceeds 85°C.

Specifications Table	
Technical Specifications	
Interface connections: ISO 5356-1 (22 mm Male)	
Compliance	≤15 mL/kPa (1.471 mL/cmH <sub>2</sub> O)
Resistance to flow @ 60 L/min	<0.82 cmH <sub>2</sub> O @ 60 L/min
Maximum Operating Pressure	13.2 kPa (2 psi)
Maximum Peak Flow	200 L/min for 3 seconds
Gas Leakage	≤10 mL/min
Invasive Max Flow Rate	60 L/min
Noninvasive Max Flow Rate	60 L/min
Standard Compliance	ISO 80601-2-74 First Edition
Environment Temperature & Humidity	
When operating	18 to 26°C 0-85% RH, N.C.
During Storage and Transport	-20 to 50°C 0-85% RH, N.C.
Patient	
Population	All patient population
System Components	
Heated Humidifier	377HTR AirLife DuoTherm Series
Humidification Chamber	377CBR
Temperature Probe Set(s)	377PR5 / 377PR6
Heater Wire Adapter(s)	377ASL / 377ADL
Relief Valve	AH115-05

#### Set-up Instructions

1. Unpack the heated humidification chamber and check for visible damage before use.  
DO NOT use a damaged humidification chamber.
2. Insert the bottom edge of the humidification chamber under the front, bottom lip on the AirLife DuoTherm Heated Humidifier and lock the humidification chamber into place under the movable mounting lever. The lever is locked when you hear a click sound (Figure 1).
3. Always position the humidifier with humidification chamber lower than the patient.
4. Remove the blue caddy (Figure 2) from the humidification chamber connection ports and unwind the feedset, leaving the spike in its air-tight protective cap.
5. Connect the circuit and dry line to the 22 mm connection ports on the chamber (Figure 3). **Note:** Ports are not flow direction dependent and either may be used.
6. Perform a ventilator pressure and leak test on the breathing system and check for occlusions before use.
7. When ready for use, remove the spike from the blue caddy and discard. Connect the water spike to the USP Sterile Water for Inhalation source. Ensure that the water source is at least 50 cm (20") above the humidification chamber (Figure 4).
8. Ensure the vent cap (Figure 4) on the water spike is closed when using Flexible Bags and open when using Rigid Containers.

9. Check the water level to ensure the chamber has filled properly and that the fill level remains below the "MAX" water level indicator during use.

Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous materials.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer at [productquality@myAirLife.com](mailto:productquality@myAirLife.com) and the competent authority of the Member State in which the user and/or patient is established.

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