

CLINICAL SUMMARY

AirLife Open™ Test Data Overview – AirLife Open vs. SouthMedic OxyMask™

DATA OVERVIEW

Testing was conducted by our research and development team. It was set up to simulate realistic clinical use situations for two open oxygen delivery masks. The products tested were AirLife™ Open and SouthMedic OxyMask™. The simulation measured the tracheal concentration of FiO₂ over the indicated O₂ flow rate ranges and patient simulations using an anatomical head model and servo-controlled test lung. The outcome shows the following:




Figure 1. AirLife™ Open vs. SouthMedic OxyMask is clinically comparable in efficacy in patient simulation testing


CONCLUSION

- The AirLife Open delivers better average tracheal (O₂) when compared to SouthMedic OxyMask at any given oxygen flow rate
- The AirLife Open oxygen mask products can cover the full 1-15 LPM range of O₂ delivery settings
- During testing masks were not sealed to the models to prevent leaks to simulate more realistic clinical practice




TESTING OVERVIEW

ITEM	DESCRIPTION
AirLife Open	M (AM-1000), L (AM-1100)
SouthMedic Oxymask	M (OM1125-8), Plus (OP1125-8)




WHY CHOOSE AIRLIFE OPEN™ OXYGEN MASK?

- Can be used from 1-15 LPM and flush
- May improve workflow:
 - Can replace multiple O₂ delivery devices with one
 - Staff require training on less devices
- Can create storage and supply room efficiencies
- Helps improve compliance



TEST METHOD

- All tests were conducted using the same end-use simulation apparatus with customized head models and an active servo test lung (ASL500, IngMar)
- Head models are based on anthropometric data and made with soft, pliable silicone, including oral and nasal structures and proportional flow paths
- An Oxigraf analyzer was used to measure O₂ concentrations in the simulated trachea and averaged over the test period for comparison
- Set on assumption that tracheal oxygen concentration correlates most closely to patient oxygenation
- Wide array of breathing patterns were derived with clinical input to cover normal/typical breathing patterns and COPD breathing patterns
- Each product was tested 9 times and placed on the head model using the included strap to simulate realistic clinical use of the masks
- **No sealing materials or tapes were used to seal the mask surfaces to the head model to ensure real case simulation**



1. Data on file at Vyair. CFD Analysis of Open O₂ Mask Rib Structure Efficacy 45 R4680.

For U.S. distribution only