EN AirLife® Oxygen Tubing

Device Description: Crush resistant, 7 ft (2.1m) oxygen tubing with vinyl tip. The tubing is non-sterile, disposable and for single-patient use.

Intended Purpose: Intended for oxygen to be delivered to the patient at a prescribed concentration and flow rate. Oxygen tubing is connected to an oxygen source on one end and a cannula on the other end to deliver oxygen to the patient.

Indications for Use: For patients prescribed supplemental oxygen, used to extend the length of the oxygen tubing on the patient interface.

Environment: Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

MR Safety: MR Safe

Patient Target Group: The intended population includes neonatal, pediatric and adult patients. Spontaneously breathing individuals with a variety of breathing conditions requiring supplemental oxygen.

Expected Clinical Benefits:

- The tubing delivers up to 15 LPM of supplemental oxygen from the gas source to the patient's interface, to treat or prevent hypoxemia in acute and chronic respiratory conditions. It is monitored by improvement in SpO2 values, PaO2 values, vital signs, and/or work of breathing.
- Oxygen tubing lengthens the nasal cannula or oxygen mask supply tubing, increasing patient's mobility and improve health-related quality of life.
- Supplemental oxygen helps reduce COPD complications by stabilizing pulmonary hypertension, reducing secondary polycythemia, or decreasing arrhythmias.
- Supplemental oxygen can decrease anxiety, improve sleep, improve mental alertness, and/or improve health-related quality of life in people with chronic respiratory disease.
- Using oxygen during exercise can improve endurance, heighten performance, and/or decrease the sensation of breathlessness. This helps patients complete physical activities, improve ability to walk and/or perform other physical activities for an extended period.

Contraindications: None known.

Warnings:

- When oxygen is in use, to reduce the risk of fire and burns, do not use near flame or heat source, and do not smoke or vape.
- Position tubing to prevent patient from becoming tangled in tubing or strangulation.
- Patient may become hypoxic if oxygen flow is interrupted.
- Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
- To reduce the risk of misconnections and patient injury, always trace tubing gas and lines back to their point of origin to verify that correct connections are made.
- High levels of oxygen may cause oxygen-induced hypoventilation, oxygen-induced hypercapnia, or oxygen toxicity. In premature infants, it may also cause retrolental fibroplasia.

Cautions:

- Flow rates above 15 LPM may cause increased resistance or back pressure in the tubing.
- To prevent interruption of therapy or injury, do not let children play with the oxygen tubing.
- To prevent excessive moisture accumulation in the tubing, do not use with heated humidification.
- To prevent equipment damage and/or interruption of therapy, do not let pets play with the oxygen tubing.
- Do not place tubing under rugs or other objects as it may obstruct flow.
- To prevent becoming tangled in tubing, position the tubing in front of you before turning to sit or stand.
- Tubing malfunction may occur as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.
- This is a single-patient use device, reuse on more than one patient may cause a risk of cross-contamination.

Residual Risks: Refer to warnings and cautions.

Necessary User Qualifications: This device is restricted to sale by or on the order of a physician. The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings and cautions as stated in the IFU.

Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

Instructions for Use:

- Remove tubing from package.
- Attach one end of the oxygen tubing to oxygen source.
- If needed, insert a tubing connector into the other end of the tubing. Attach the other end of the tubing to an oxygen mask, nasal cannula or another oxygen tubing.
- Adjust the flow control knob to the prescribed liter flow.
- Check connections for leaks and verify flow is coming from the oxygen mask or nasal cannula.
- Discard and replace if the tubing becomes soiled, discolored, or stiff. This product is disposable and is not intended to be cleaned, sterilized, or reused.

Instructions for Use in Combination with Accessories /

Or Other Devices: AirLife oxygen tubing end connectors are designed to work with gas outlets that comply to ISO 5359:2014. This includes, but is not limited to oxygen tanks, oxygen flowmeters, oxygen blender, tubing connectors and bubble humidifiers. The tubing trumpet end connector is compatible with a 5 mm - 7 mm tapered oxygen fitting or ¼ inch (6.35 mm) tapered oxygen fitting (Figure A). The style of the tapered fitting will vary depending on the oxygen gas

EN

source or oxygen device.

- Place the ribbed or trumpet end connector between your thumb and forefinger (Figure B).
- Align the end of the connector to tapered fitting, then push the end connector over the tapered fitting until snug (Figure B).
- Gently tug on the tubing to ensure a secure fit.


Fig. A


Fig. B

Safe Disposal: Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

Incident Reporting: Notice to the user and/or patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. If your physical symptoms worsen or you experience a sudden change in your condition (e.g., increased shortness of breath, fever, dizziness) call your doctor.

Storage Instructions / Conditions:

No special storage conditions.

REF

Catalogue Number

Manufacturer

LOT

Lot Number

Expiration Date

X

Quantity

Authorized Representative in the European Community

Not Made with Natural Rubber Latex

Medical Device

Not Made with DEHP

Cautions

Consult Instructions for Use

Temperature Limit

U.S. Federal Law restricts this device to sales by or on order of a physician

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© 2019 Vyaire. Vyaire, the Vyaire Logo and AirLife are trademarks or registered trademarks of Vyaire Medical, Inc., or of its affiliates.
AirLife Australia Holding Pty. Ltd.
Vyaire Medical, Inc.