

AVANOS*



BALLARD* RESTORE2
ENDOTRACHEAL
TUBE CLEANER

Instructions for Use



Federal (USA) law restricts this device to sale by or on the order of a physician. Single Use Device. Do not resterilize.

Product Description

The BALLAR2® Restore2 Endotracheal Tube Cleaner is a sterile single use device. The Restore2 wiper is a radially expanding squared edged wiper that can be used to clear away mucus, secretions and bacteria from the inside of the endotracheal tube (ETT). The attached Collection Adapter (CA) fits most standard adult endotracheal tubes (ETT).

Indications for Use

The Restore2 is a single use device indicated for the removal of secretions from adult endotracheal tubes sized 7 mm to 8.5 mm.

Warnings

- Do not use excessive force when removing the Restore2. Wiping action should be gentle and smooth. If any hard resistance is encountered, release the pressure on the syringe and carefully remove the device.**
- Do not use the device if it appears damaged in any way.**
- Do not use in ETT smaller than 7 mm or larger than 8.5 mm.**
- Do not cut ETT while Restore2 device is attached.**
- Do not overinflate. (Syringes included are for use with 7.0-7.5 or 8.0-8.5 ETT only.) Use specific syringe for specific ETT size.**
- Do not insert past the appropriate centimeter depth marking.**
- Do not reuse or reprocess this medical device. Reuse or reprocessing may: (a) adversely affect the known biocompatibility of the device, (b) compromise the structural integrity of the device, (c) lead to the device not performing as intended, (d) create a risk of contamination and cause the transmission of infectious diseases resulting in patient injury, illness or death.**
- Do not use on tracheostomy patients.**

Cautions

- Rx only.
- Do not use this device if package is damaged or opened. Do not use this device if it shows signs of damage, e.g., crimps, kinks, or crushed areas.
- This device is for single use only.
- Proper surgical procedures and techniques are necessarily the responsibility of the medical profession.
- These instructions are furnished for informational purposes only. Each user must evaluate the appropriate use of this device, case by case, based on medical training, experience, and the type of procedure employed.
- Dispose of used products in accordance with established hospital protocol for hazardous waste.
- Verify all use connections of this device are secure prior to use.

Set Up

- First perform endotracheal tube (ETT) suctioning with either a closed or open suction catheter according to AARC guidelines.
- Check that the ETT is secure and note centimeter markings at the lip or teeth.
- Using universal precautions and clean gloves, remove the Restore2 with Collection Adapter (CA) from the sterile pouch.
- Select the appropriate syringe (Red Label 7.0-7.5 attached to wiper in package, or Blue Label 8.0-8.5). Connect to the Restore2 catheter as needed and test the wiper by firing the syringe completely and depressing. Withdraw plunger fully and depress fully. Assume that the cleaning member inflates. Note that it may not inflate symmetrically outside of the ETT. This is normal. Withdraw the syringe fully and confirm that the balloon spontaneously returns to its unexpanded position tight against the catheter.
- Rapidly insert the CA and Restore2 between the endotracheal tube and the ventilator, with the ventilator connection superior and the insertion port for the Restore2 device inferior.
- Advance the Restore2 into the ETT taking care to align the centimeter markings on the Restore2 catheter with the visible centimeter markings on the ETT. In most cases on uncut ETTs, the catheter may be fully inserted. **Do not over-insert.**
- Hold the CA adapter and Restore2 sheath with one hand. At the end of patient or ventilator inspiration, completely fire the syringe and pull the Restore2 cleaning member back into the CA over 3 to 5 seconds. **Do not use excessive force to wipe the ETT clear.** Restore2 users should feel a smooth wiping action.
- If the device cannot be pulled back, release the syringe and remove the catheter. Recheck the inflation volume (above) that is recommended and repeat the procedure without fully depressing the syringe. If the device still cannot be pulled back, release the syringe and remove the catheter. Abort the procedure. Remove the CA and reconnect the ventilator. Fiber optic bronchoscopy may be recommended.
- The CA and Restore2 can be removed from the ventilator circuit as a unit and discarded according hospital policy for hazardous waste.
- If factors indicate the need for a second pass the Restore2 can be removed from the CA and withdrawn into its protective sheath. Confirm function of the cleaning member by temporary re-inflation and deflation. Repeat the cleaning procedure in steps 7 through 11 above.
- If endotracheal tube obstruction is still considered, evaluate other means of investigation such as fiber-optic bronchoscopy.

 Single Use Only	 STERILE EO	 Do not use if package is damaged	 Do not resterilize
 Not made with natural rubber latex	Rx Only	 Caution	 Consult instructions for use

Warranty

Avanos Medical warrants the instrument against defects in material and workmanship subject to the conditions that the instrument has been used in accordance with its operating instructions at all times.

The end customer shall bear the cost of transportation and risk of return shipments.

Opening the equipment for performance of any repairs or modifications to the equipment by unauthorized persons shall release Avanos Medical from any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranties expressed or implied.

Representations and warranties made by any person, including but not limited to dealers and representatives of Avanos Medical which are inconsistent or in conflict with the terms of this limited warranty, shall not be binding upon Avanos Medical unless reduced to writing and approved by an expressly authorized officer of Avanos Medical.

This limited warranty is the complete and exclusive statement of warranty which Avanos Medical agrees to provide with respect to this product and it shall supersede all prior and contemporaneous oral or written agreements, understandings, proposals, and communications pertaining to the subject matter hereof.

Distributed in the USA by Avanos Medical Sales, LLC,
5405 Windward Parkway, Alpharetta, GA 30004 USA.
In USA, 1-844-4AVANOS (1-844-428-2667). www.avanos.com

 Avanos Medical, Inc.,
5405 Windward Parkway,
Alpharetta, GA 30004 USA.

 Avanos Medical Belgium BVBA,
 Leonardo da Vincilaan 1,
1930 Zaventem, Belgium.

Sponsored in Australia by Avanos Medical
Australia Pty Ltd, 475 Victoria Avenue,
Chatswood, NSW 2067, Australia.

*Registered Trademark or Trademark of Avanos Medical, Inc.,
or its affiliates. © 2018 AVNS. All rights reserved. 2020-12-10
15-M1-680-01 / 50000311