Care and Service of Conventional Laryngoscope Systems

R_X ONLY Not Made with Natural Rubber Latex
USA Law Only

CLEAN IMMEDIATELY AFTER USE TO PREVENT RESIDUES FROM DRYING

WARNINGS and CAUTIONS

- DO NOT USE THE FOLLOWING: Ultrasonic clean, dry heat sterilization or Potassium Hydroxide liquids.
- REMOVE BATTERIES BEFORE cleaning, high-level disinfecting, or sterilizing of the laryngoscope handles.
- ENSURE blade is fully engaged on handle prior to use.
- CHECK laryngoscope system functionality prior to use.

CLEANING

Laryngoscope Blades: Rinse the laryngoscope blade under cool running tap water, until all visible soil is removed. Immerse laryngoscope blade, in a presoak enzymatic cleaner solution, which was prepared in accordance with enzymatic cleaner manufacturer's recommendations, for a minimum of two minutes. After the presoak, scrub the laryngoscope blade using a soft bristled brush, paying attention to hard to reach areas. Remove laryngoscope blade from enzymatic cleaner solution. Rinse with lukewarm running tap water for a minimum of one minute in order to remove any residues and visible soil. Then, immerse the laryngoscope blade in an enzymatic detergent prepared according with recommendations. Scrub the laryngoscope blade thoroughly by using a soft bristle brush, while ensuring that all hard-to-reach areas are reached and visible soils/residues removed. Remove laryngoscope blade from the enzymatic detergent and thoroughly rinse it with cool running tap water in order to remove all residual detergent. Dry the laryngoscope blade with lint free, clean cloth or filtered pressurized air. Follow with HIGH-LEVEL DISINFECTION or STERILIZATION PROCEDURE.

Laryngoscope Handles: Disassemble the laryngoscope handle by removing the bottom cap and batteries. Rinse the disassembled laryngoscope handle under cool running tap water, until all visible soil is removed. Immerse disassembled laryngoscope handle and caps, in a presoak enzymatic cleaner solution, which was prepared in accordance with enzymatic cleaner manufacturer's recommendations, for a minimum of two minutes. After the soak, scrub the cap and laryngoscope handle using a soft bristled brush paying attention to hard to reach areas. Remove disassembled laryngoscope handle from enzymatic cleaner solution. Rinse with lukewarm running tap water for a minimum of one minute in order to remove any residues and visible soil. Then, immerse disassembled laryngoscope handle and cap in an enzymatic detergent. Scrub all disassembled components thoroughly by using a soft bristle brush, while ensuring that all hard-to-reach areas are reached and visible soils/residues removed. Remove disassembled laryngoscope handle from the enzymatic detergent and thoroughly rinse all parts with cool running tap water in order to remove all residual detergent. Dry all disassembled components (top and bottom cap and handle) with lint free, clean cloth or filtered pressurized air. Reassemble all components, except batteries, and follow with HIGH-LEVEL DISINFECTION or STERILIZATION PROCEDURE.

AUTOMATED CLEANING

Automated cleaning for Blades: Load laryngoscope blade into an Automated washer per ISO15883 or equivalent and follow the below recommended cycle and detergent. Equivalent cycles and detergent can be utilized. Follow with HIGH-LEVEL DISINFECTION or STERILIZATION PROCEDURE

Automated cleaning for Handles: Disassemble the laryngoscope handle by removing the bottom cap and batteries. Place the handle and bottom cap into an Automated washer per ISO15883 and follow the below recommended cycle and detergent. Equivalent washers, cycles, and detergent can be utilized. Follow with HIGH-LEVEL DISINFECTION or STERILIZATION PROCEDURE.

Phase	Recirculation time (Minutes)	Temperature	Detergent Type and Concentration
Pre-Wash	02:00	Cold Tap Water	N/A
Wash	02:00	43°C Tap Water	Valsure Enzymatic 2 mL/L
Rinse	01:00	43°C Tap Water	N/A
Dry Time	07:00	90°C	N/A

DISINFECTION

Laryngoscope Blades: To achieve a high-level disinfection, Cidex® OPA or 2.4% Glutaraldehyde solution may be used according to high-level disinfectant manufacturer's instructions. The laryngoscope blade shall be fully immersed into the disinfectant. After exposure, the laryngoscope blade must be thoroughly rinsed by fully immersing in purified water. Dry with lint free, clean cloth, or filtered pressurized air. If not functional, review battery/lamp testing instructions below.

Laryngoscope Handles: To achieve a high-level disinfection, Cidex® OPA or 2.4% Glutaraldehyde solution may be used according to high-level disinfectant manufacturer's instructions. The laryngoscope handle shall be fully immersed into the disinfectant. Ensure that no air bubbles are present on the exterior of the device. After exposure, the laryngoscope handle must be thoroughly rinsed by fully immersing in purified water. Dry with lint free, clean cloth or filtered pressurized air. Load handle with batteries and test it for proper function prior to use. If not functional, review battery/lamp testing instructions below.

STERILIZATION

Steam Sterilization: Insert laryngoscope handle (without batteries) or laryngoscope blades in appropriate autoclave pouch.

	(A) GRAVITY DISPLACEMENT STEAM
Temperature	121° C (250° F)
Cycle Time	30 Minutes
Dry Time	15 Minutes

STERIS® AMSCO® V-PRO® Conventional Laryngoscope Blades and Handles are compatible with:

- Amsco® V-PRO® 1 Low Temperature Sterilization System
- Amsco® V-PRO® 1 Plus Low Temperature Sterilization System
- Amsco® V-PRO® maX Temperature Sterilization System

STERRAD® Conventional Laryngoscope Blades and Handles are compatible with:


- STERRAD® 100NX® System (Standard and EXPRESS Cycle)
- STERRAD® NX® System (Standard Cycle)
- STERRAD® 100S and 200 System (Short cycle Outside US)
- STERRAD® 50 System

BATTERY / LAMP TESTING

Laryngoscope systems (laryngoscope blades and laryngoscope handles) must be tested after each cleaning/high-disinfection or cleaning/sterilization processing, prior to use, and after battery replacement. Fully engage blade on handle. If the light guide fails to light, check the batteries/lamp for correct placement or replace with new ones.

To Install Batteries: Place two 1.5 Volt batteries, positive end first into Penlight, Medium, or Large handle (Fig 1) and reinstall the end cap. In Stubby handle, insert battery cartridge according to instructions marked on cartridge (Fig 2). If Stubby handle light guide fails to light, batteries may need to be reversed for correct polarity. Attach Conventional blade in operating position. Lamp should illuminate.

To Change Lamp: Prior to lamp removal, allow lamp time to cool. To replace the lamp, unscrew the lamp from the blade, and pull the lamp out. Reinstall a new lamp in the laryngoscope blade, and tighten it. Ensure the lamp is seated properly and verify lamp function prior to use.

SunMed, LLC,
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA

Manufactured for SunMed, LLC,
a Subsidiary of SunMed Holdings, LLC
www.Sun-Med.com
Made in Pakistan

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