FLOCAP Instructions For Use

R_X ONLY

USA Law Only

Not Made with Natural Rubber Latex


NON STERILE

5°C 41°F

37°C 98°F

Single Use

REF R500P21

EN ASSEMBLY INSTRUCTIONS & INSTRUCTIONS FOR USE

1.0 SYSTEM OVERVIEW

1.1 Description: The FLOCAP is a simple, disposable, single-use CO₂, and flow indication device designed for placement between a breathing device and a patient's endotracheal tube or mask for visualization of exhaled CO₂ and to give a visual indication of the presence of expiratory flow. The CO₂ indicator will assist the caregiver in verifying proper ET tube placement. A lack of color change may indicate improper intubation. Exhaled gas passes through the device and indicates a range of end-tidal CO₂. The flow indicator allows the caregiver to visually detect if the patient is still exhaling.

1.2 Indications for use: The FLOCAP is used to provide a semi-quantitative visualization of the CO₂ in the patient airway. It is an adjunct in patient assessment, to be used in conjunction with other methods to determine clinical signs and symptoms by or on the order of a physician.

- The FLOCAP has a visual indicator to visually detect the end of exhalation.
- For use up to 24 hours.
- For patients greater than 15 kg (33 lbs.)
- Environment of use – hospital, sub-acute, pre-hospital, transport

1.3 Contraindications

DO NOT use the FLOCAP for the detection of Hypercapnia/Hypercarbia.
DO NOT use the FLOCAP for the detection of main-stem bronchial intubation.
DO NOT use the FLOCAP during mouth to tube ventilation.
DO NOT use the FLOCAP to detect oropharyngeal tube placement.

- Standard clinical assessment must be used. For patients greater than 15 kg (33 lbs.)

When low pulmonary perfusion coincides with accidental esophageal intubation, colorimetric CO₂ indication cannot be properly interpreted. However, if proper tube placement is ascertained by independent means, then the FLOCAP may be used to help assess the progress of positive pressure ventilation as evidenced by an increase in end-tidal CO₂.

1.4 Warnings

- DO NOT use if you have blue-yellow color blindness.
DO NOT use on patients with body weight less than 15kg (33lbs) due to the potential for rebreathing exhaled CO₂.
DO NOT use the FLOCAP with devices that elevate humidity such as nebulizers or heated humidifiers.
DO NOT use in the presence of the following agents: atropine, infasurf, naloxone, intratracheal epinephrine, trichloroethylene, chloroform.
DO NOT use the FLOCAP for a duration of more than 24 hours.
- For use with air and oxygen only; chemical interactions will affect device accuracy.
 - Read the entire contents of this operating manual before using the FLOCAP.

Periodic color changes reflect the breathing pattern of the patient. However, a permanent purple to purple-beige color indicates a lack of exhaled CO₂, the cause of which requires immediate attention. If the purple-beige color does not appear during a breathing cycle, this may indicate a significant degree of CO₂ rebreathing and should also be of immediate concern.

! 1.5 Cautions

- Single use only, discard after use.
- Interpreting color change before 6 complete breaths may lead to a false result.
- Excessively low cardiac output will result in low CO₂ content in the lungs.
- Excessive CO₂ in the stomach may cause erroneous color change.
- Contamination, liquid water or excessive moisture may cause poor visibility, limit litmus function or inhibit spinner motion. The FLOCAP should be replaced immediately if this occurs.
- Reflux of any kind into the FLOCAP may compromise its accuracy or performance. The FLOCAP should be replaced immediately if this occurs.
- The FLOCAP is an adjunct assessment tool and should not be relied upon as the sole means of verifying proper intubation. Follow your institutional guidelines for verifying proper intubation in addition to use of the FLOCAP.
- During cardiac arrest, CO₂ levels in the lungs may be too low to affect a color change in the FLOCAP. Re-establishment of pulmonary blood flow is required for the FLOCAP to function properly.
- Storage for extended periods at temperatures above 37°C may result in reduced shelf life.
DO NOT use if the package is already unsealed.
DO NOT open the pouch until ready to use, prolonged exposure to ambient air may affect accuracy or performance.
- The bright colored vane of the flow indicator is not intended to be a quantitative flow meter.
- At the end of exhalation, the flow indicator's momentum will cause it to spin briefly after exhalation has ceased.
DO NOT use in an environment with insufficient lighting.
- Inspect the FLOCAP for damage prior to use.
DO NOT use in the presence of acidic liquid or medication.
- Avoid exposure to strong sunlight and other sources of ultraviolet light.
- The FLOCAP does not take the place of traditional end tidal CO₂ monitoring which remains necessary to provide quantitative measurements and patient alarms.

3.0 2.0 INSTRUCTIONS FOR USE

- To open: Open the pouch, remove and inspect the FLOCAP. Ensure the element inside the housing is still purple. If its color looks to be closer to beige the device should be discarded. Some yellowing around the edge is permissible.
- To connect: Attach the FLOCAP securely between the ventilation source and the mask or endotracheal tube, by pushing together and rotating slightly, as illustrated here: (fig.1)

A

Connect to breathing device here.

B

Connect to endotracheal tube here.
- For CO₂ Detection: Ventilate the patient with 6 complete breaths. At the end of the 6 breath cycles, check the color of the FLOCAP indicator.

- If the color of the indicator is yellow, continue to ventilate the patient and monitor clinical cues for adequate ventilation.
 - If the color of the indicator is beige, continue to ventilate with 6 additional breaths and recheck. If the color has not moved towards yellow, this indicates insufficient exhaled CO₂. This requires immediate attention. Employ your institution's protocol to verify proper tube placement and adequate gas exchange.
 - If the color remains purple, or a color similar to purple, the patient may not have a patent airway or the endotracheal tube may not be positioned correctly. This requires immediate attention. Employ your institution's protocol to verify proper tube placement and adequate gas exchange.
 - Continue to monitor the color of the CO₂ indicator throughout the entire time the device is in use.
DO NOT exceed 24 hours of use.
 - During the course of ventilation, if the FLOCAP returns to and remains purple or beige colored, this indicates insufficient exhaled CO₂. This requires immediate attention.
4. Not intended to be a quantitative indicator. Follow your institution's guidelines for breath frequency and tidal volume during resuscitation and for avoiding auto-PEEP during resuscitation. At the end of exhalation, the flow indicator's momentum will cause it to spin briefly after exhalation has ceased.

3.0 3.0 COLOR & FLOW INDICATOR INTERPRETATION

Purple (0%) ... No carbon dioxide detected
Beige (1-2%) ... Exposed to CO₂, approx. 1% to 2%
Yellow (5%) ... Exposed to 5% or greater CO₂
Spinning* ... Patient is exhaling
Not Spinning ... Exhalation flow < 3 LPM

**At the end of exhalation, the flow indicator's momentum will cause it to spin briefly after exhalation has ceased.*

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