Care and Service of Reusable
GreenLine®Fiber Optic
Laryngoscope Blades

REF

Intended Use: A laryngoscope is a device which helps in intubation during the administration of general anesthesia or for mechanical ventilation.

Care and Service Document: Fiber optic laryngoscope blades must be processed prior to first use, and between each use, using any of the methods presented in this care and service document:

- CLEANING + (A) COLD SOLUTION HIGH-LEVEL DISINFECTION
- CLEANING + (B) GRAVITY AUTOCLAVE STERILIZATION
- CLEANING + (C) PRE-VACUUM STEAM STERILIZATION
- CLEANING + (D) STERIS AMSCO V-PRO STERILIZATION
- CLEANING + (E) STERRAD STERILIZATION

WARNINGS and CAUTIONS

- DO NOT USE THE FOLLOWING: Ultrasonic clean, dry heat sterilization or Potassium Hydroxide liquids.
- CHECK laryngoscope system functionality, and ensure blade is fully engaged on handle prior to use.
- DISCARD any component that shows evidence of damage or deterioration.

NOTE: The laryngoscope blade **MUST** be removed from the handle prior to cleaning, disinfection, and/or sterilization. To clean, disinfect or sterilize the handle, consult the Care and Service Instructions for Reusable GreenLine Fiber Optic Laryngoscope Handles with Removable Lamp Cartridge.

MANUAL CLEANING

Immediately after use, the laryngoscope system should be rinsed under cool running tap water until all visible soil is removed. Rinse hard-to reach areas with a cleaning gun for at least one (1) minute. Remove the blade from the handle head (Figure 1).

1. Immerse the blade in a pre-soak enzymatic cleaning solution, which was prepared in accordance to manufacturer's recommendations.
2. Brush the hard-to-reach areas of the blade with a soft brush for at least one (1) minute. Pay attention to hard to reach areas.
3. Make sure all the surfaces are completely moistened with the enzymatic cleaning solution.
4. Remove the blade from the enzymatic cleaning solution.
5. Rinse the blade completely under lukewarm running tap water for a minimum of one (1) minute to remove the cleaning solution completely. Check for cleanliness. If dirt is visible, repeat steps 1-5 until visibly clean.
6. Dry with lint-free, clean cloth or filtered pressurized air. Follow with HIGH-LEVEL DISINFECTION or STERILIZATION.


FIG. 1

AUTOMATED CLEANING

1. Remove laryngoscope blade from the handle head (Fig. 1).
2. Load laryngoscope blade into an automated washer per ISO 15883 or equivalent and follow the below recommended process chart. Equivalent cycles and detergent can be utilized.
3. Follow with HIGH-LEVEL DISINFECTION or STERILIZATION of the blade.

Phase	Recirculation time (Minutes)	Temperature	Detergent Type and Concentration
Pre-Wash	02:00	Cold Tap Water	N/A
Wash	02:00	43°C Tap Water	Valsure Enzymatic 2 mL/L
Rinse	01:00	43°C Tap Water	N/A
Dry Time	07:00	90°C	N/A

5-5230-XX	Miller American Profile Blade, Light on Right	5-5243-XX	Mac English Profile Blades, Left Handed
5-5231-XX	Mac American Profile Blades	5-5247-XX	Polio Mac English Profile Blades
5-5232-XX	Mac English Profile Blades	5-5248-XX	Port Miller American Profile Blades
5-5233-XX	Miller English Profile Blades	5-5249-XX	Snow Miller American Profile Blades
5-5234-XX	Wisconsin Miller American Profile Blades	5-5254-XX	Mac American Profile Blades, Improved View
5-5235-XX	Phillips Miller American Profile Blades	5-5271-XX	Mac English Profile Blades, Reduced Flange
5-5236-XX	Cranwall Miller American Profile Blades	5-5272-XX	Mac German Profile Blades
5-5237-XX	Robertshaw Miller English Profile Blades	5-5273-XX	Miller German Profile Blades
5-5238-XX	Oxford Miller English Profile Blades	5-5282-XX	Sun-Flex™ Mac German Profile Blades

(A) COLD SOLUTION HIGH-LEVEL DISINFECTION

1. Prepare the disinfectant solution CIDEX®-OPA (0.55% ortho-Phthalaldehyde based disinfectant) according to the manufacturer's directions.
2. Completely immerse the blade in the disinfecting solution. Make sure that all surfaces are completely moistened. Brush the hard-to-reach areas during the exposure time for at least one (1) minute. Allow the device to soak for a minimum of twelve (12) minutes.
3. Remove the blade from the disinfecting solution.
4. Rinse by immersing the blade in demineralized water for one (1) minute. Repeat two (2) additional times for a total of three (3) rinse cycles, utilizing fresh demineralized water each time. During the rinsing process, ensure the device is agitated.
5. Dry with a lint-free clean cloth or medical compressed air.

STERILIZATION

Use one (1) of the following four (4) sterilization methods to sterilize SunMed® reusable GreenLine® fiber optic laryngoscope blades.

Insert blade in an appropriate autoclave pouch.

	(B) GRAVITY AUTOCLAVE	(C) PRE-VACUUM STEAM
Temperature	121° C (250° F)	134° C (273° F)
Cycle Time	30 Minutes	3 Minutes
Dry Time	15 Minutes	20 Minutes

D) STERIS® AMSCO® V-PRO® STERILIZATION.

- Amsco® V-PRO® 1 Low Temperature System
- Amsco® V-PRO® 1 Plus Low Temperature System
- Amsco® V-PRO® maX Low Temperature System

(E) STERRAD® STERILIZATION.

- STERRAD® 100NX® System (Std and EXPRESS Cycle)
- STERRAD® NX® System (Std Cycle)
- STERRAD® 100S and 200 System (Short Cycle Outside US)
- STERRAD® 50 System

Fiber Optic blades that are compatible in STERRAD® Systems must be manufactured in or after May 2012 or have a lot number "Me" or greater. Upper case "M" corresponds to the year 2012 and lower case "e" corresponds to the manufacturing month, May. Lot number will be visible on every blade (Fig.1). Lot numbering key:

FIRST LETTER - YEAR							
M = 2012	N = 2013	O = 2014	P = 2015	Q = 2016	R = 2017	S = 2018	T = 2019
U = 2020	V = 2021	W = 2022	X = 2023	Y = 2024	Z = 2025	A = 2026	B = 2027
C = 2028	D = 2029	E = 2030	F = 2031	G = 2032	H = 2033		

SECOND LETTER - MONTH		
a = January	e = May	i = September
b = February	f = June	j = October
c = March	g = July	k = November
d = April	h = August	l = December

Prior to use, follow the BATTERY / LAMP TESTING instructions below:

BATTERY / LAMP TESTING

Laryngoscope systems must be tested after each reprocessing cycle and prior to use. To test:

1. Fully engage blade on compatible laryngoscope handle. If the lamp fails to light, check the batteries and lamp for correct placement. Replace with new batteries or lamp if needed.
2. If lamp fails to light following additional blade engagements, do not use laryngoscope.

SunMed, LLC,
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA

Manufactured for SunMed, LLC,
a Subsidiary of SunMed Holdings, LLC
www.Sun-Med.com
Made in Pakistan

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Fabricado para SunMed, LLC,
una subsidiaria de SunMed Holdings, LLC