
# AirLife®

EN

English

## Salter Labs™ EtCO2 Gas Sampling Lines

### DEVICE DESCRIPTION:

The device is non-sterile, disposable, single-patient multi-use sampling line that enables sampling of exhaled gases. The sampling line is available in different end connectors.

### INTENDED PURPOSE:

EtCO2 Gas Sampling Lines are intended to be used where exhaled gas is monitored.

### INDICATIONS FOR USE:

EtCO2 Gas Sampling Lines are intended to be used to sample exhaled gas in intubated or non-intubated patients requiring expired gas monitoring.

### ENVIRONMENT:

Hospital, sub-acute, pre-hospital.

### PATIENT TARGET GROUP:

Infant, pediatric, adult.

### INTENDED USERS:

Trained medical professionals.

### EXPECTED CLINICAL BENEFITS:

Capnography provides early detection and warning of changes in a patient's ventilatory status and cardiac perfusion. This enables clinicians to adjust medication dosage or apply therapeutic intervention to prevent hypoventilation and apnea. ETCO2 monitoring provides information on endotracheal tube placement, the effectiveness of chest compressions and the return of circulation.

### CONTRAINDICATIONS:

None Known.

### WARNINGS:

- If oxygen is in use, do not use near flame or heat source.
- Position the sampling line to reduce the possibility of patient entanglement or strangulation.
- Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
- To reduce the risk of misconnections and patient injury, always trace tubing from capnography monitor to the patient.
- Do not reuse on more than one patient as this may increase risk of contamination and infection.

### CAUTIONS:

- Consult physician if patient develops infection, skin irritation or material sensitization.
- Do not place anything on supply tubing that may obstruct flow.
- Accuracy of CO2 monitoring can be diminished during the time the patient is mouthbreathing or if the patient has a deviated septum.
- Excessive patient secretion or a build-up of liquids in the tubing may occlude the sampling and require more frequent replacement.
- This is a single-patient multi-use product, reuse on more than one patient may cause a risk of cross-contamination.
- Check for leaks. Loose connections may compromise ventilation or cause an inaccurate EtCO2 measurement.
- Do not sterilize or clean device as it may affect the performance of the device and cause damage to the monitor.
- Check sample line regularly during use to ensure that no kinks are present. Kinked tubing may cause inaccurate EtCO2 readings.
- EtCO2 values, waveforms and duration of use can be affected by certain ambient environment conditions and certain patient conditions.

### RESIDUAL RISKS:

Refer to Warnings and Cautions.

### NECESSARY USER QUALIFICATIONS:

The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings, and cautions as stated in the IFU.

Rx Only: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

### INSTRUCTIONS FOR USE:

Follow instructions to avoid risk of infection and contamination.

- Connect the gas sampling line to the capnography monitor.
- Connect the other end of the gas sampling line to the luer connector of the circuit adaptor or sampling line.
- Check the monitor display for EtCO2 value and/or waveform.

This product is disposable and is not intended for prolonged use greater than 7 days.

### MONITORING AND MAINTENANCE:

- Replace the sampling line when a blockage or an occlusion alarm signals on the monitor.
- If there is no waveform displayed, check to see if the sample line is disconnected or kinked. If the sample line is not disconnected or kinked replace the line.
- For high humidity environment or long-term monitoring, use the EtCO2 sampling line with Nafion tube and the yellow-end connector.
- Check sampling line regularly during use to ensure that no kinks are present. Kinked tubing may cause inaccurate EtCO2 readings.

### CLEANING / REUSE INSTRUCTIONS:

A discontinued use of this product on the same patient, with periods of use separated by a period of non-use, is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.

### SAFE DISPOSAL:

Decontaminate and dispose of all potentially biohazardous material. Dispose of device in accordance with local, state or national regulations.

### INCIDENT REPORTING:

Notice to the user that any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myairlife.com, and the competent authority of the Member State in which the user and/or patient is established.

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