Salter® Reflective EtCO2 Nasal Sampling Line

REF 4MSF1-0, 4MSF2, 4MSF3 and Associated SO Numbers

EN Reflective EtCO2 Nasal Sampling Line

DEVICE DESCRIPTION: The device is a non-sterile, disposable, single patient multi-use cannula that allows sampling of patient's exhaled gases.

INTENDED PURPOSE: The reflective ETCO2 Nasal Sampling Cannula is intended to be used to sample exhaled gas.

INDICATIONS FOR USE: The reflective ETCO2 Nasal Sampling Cannula is intended to be used to sample exhaled gas in spontaneously breathing patients requiring expired gas monitoring.

ENVIRONMENT: Hospital, sub-acute, pre-hospital, surgical centers, dental offices.

PATIENT TARGET GROUP: Infant, pediatric, adult.

EXPECTED CLINICAL BENEFITS: Capnography provides early detection and warning of changes in a patient's ventilatory status and cardiac perfusion. This enables clinicians to adjust medication dosage or apply therapeutic intervention to prevent hypoventilation and apnea. Supplemental oxygen can be delivered to treat hypoxemia measured by changes in SpO2, blood gases and relief of shortness of breath.

CONTRAINDICATIONS: None Known.

WARNINGS:

- Position tubing to avoid strangulation.
- Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
- To reduce the risk of misconnections and patient injury, always trace tubing from gas source to the medical device before connecting.
- Contains Tris(nonylphenyl) phosphite (TNPP) above 0.1% w/w.

CAUTIONS:

- Consult physician if patient develops infection, skin irritation or material sensitization.
- Do not place anything on supply tubing that may obstruct flow.
- Not for use on patients with obstruction of both nostrils.
- Accuracy of CO2 monitoring can be diminished during the time the patient is mouth breathing or if the patient has a deviated septum.
- Excessive patient secretion or a build-up of liquids in the tubing may occlude the sampling and require more frequent replacement.
- Single patient use, reuse may increase risk of contamination and infection

RESIDUAL RISKS: Refer to Warnings and Cautions.

NECESSARY USER QUALIFICATIONS: This device is restricted to sale by or on the order of a physician. There are no additional requirements for special facilities, special training, or particular qualifications for the use of this device.

INSTRUCTIONS: Follow instructions to avoid risk of infection and contamination.

- 1) Connect the gas sampling line to the capnography monitor.
- 2) Insert the prongs into the patient's nares.
- 3) Wrap the headset loop up and over the ears.
- 4) Squeeze the sides of the bolo and glide up under the chin.
- 5) Check the monitor display for EtCO2 value and/or waveform.

INSTRUCTIONS FOR USE IN COMBINATION WITH ACCESSORIES OR OTHER DEVICES: Sampling line connects to capnography monitor.

SAFE DISPOSAL: Dispose of device in accordance with local, state or national regulations.

INCIDENT REPORTING: Notice to the user and/or patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

PERFORMANCE CHARACTERISTICS:

Reference Number	4MSF1, 4MSF3	4MSF2, 4MSF3-PED	4MSF3-INF
Patient Population	Adult	Pediatric	Infant
EtCO2 Sampling Line (I.D.)	0.04" (1.02mm)	0.04" (1.02mm)	0.04" (1.02mm)
EtCO2 Sampling Rate	< 0.065 LPM	< 0.065 LPM	< 0.065 LPM
EtCO2 End Connector	Reflective, Female Luer	Reflective, Female Luer	Reflective, Female Luer
EtCO2 Tubing Length	6' (1.8m), 13' (4.0m)	6' (1.8m), 13' (4.0m)	6' (1.8m), 13' (4.0m)
Prong (I.D.)	0.111"	0.111"	0.078"
Oxygen Flow Rate	≤ 4 LPM	≤ 3 LPM	≤ 3 LPM
Oxygen Tubing Length	7' (2.1m)	7' (2.1m)	7' (2.1m)

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MT Promedt Consulting GmbH
MT Promedt Consulting Ltd.
MedEnvoy UK Limited