

AirLife®

Salter Style®

EN Salter Labs™ EtCO2/O2 Divided Nasal Cannula

DEVICE DESCRIPTION:
The device is a non-sterile, disposable, single patient multi-use cannula that provides supplemental oxygen while sampling patient's exhaled gases. This product is disposable and is not intended for prolonged use greater than 7 days.

INTENDED PURPOSE:
EtCO2/O2 Divided Nasal Cannula is intended to be used where low flow oxygen is administered and exhaled gas is monitored.

INDICATIONS FOR USE:
The EtCO2/O2 Divided Nasal Cannula is intended to be used with spontaneously breathing patients who require supplemental oxygen and exhaled gas monitored.

ENVIRONMENT FOR USE:
Hospital, sub-acute, pre-hospital, surgical center and dental offices.

PATIENT TARGET GROUP:
Adult.

INTENDED USERS:
Trained medical professionals.

NECESSARY USER QUALIFICATIONS:
The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings, and cautions as stated in the IFU.

EXPECTED CLINICAL BENEFITS
Capnography provides early detection and warning of changes in

a patient's ventilatory status and cardiac perfusion. This enables clinicians to adjust medication dosage or apply therapeutic intervention to prevent hypoventilation and apnea. Supplemental oxygen up to 6 LPM can be delivered to treat hypoxemia measured by changes in SpO2, blood gases and relief of shortness of breath.

CONTRAINDICATIONS
None known.

- ⚠ WARNINGS**
- If oxygen is in use, do not use near flame or heat source.
 - Position the tubing to reduce the possibility of patient entanglement or strangulation.
 - Patient may become hypoxic if oxygen flow is interrupted.
 - Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
 - To reduce the risk of misconnections and patient injury, always trace tubing from capnography monitor to the patient.
 - Do not reuse on more than one patient as this may increase risk of contamination and infection.

- ⚠ CAUTIONS**
- Consult physician if patient develops infection, skin irritation or material sensitization.
 - Do not place anything on supply tubing that may obstruct flow.
 - Not for use on patients with obstruction of both nostrils.
 - Accuracy of CO2 monitoring can be diminished during the time the patient is mouth breathing or if the patient has a deviated septum.
 - Excessive patient secretion or a build-up of liquids in the tubing may occlude the sampling line and require more frequent replacement.
 - This is a single-patient multi-use product, reuse on more than one patient may cause a risk of cross-contamination.
 - Check for leaks. Loose connections may compromise ventilation or cause an inaccurate EtCO2 measurement.
 - Do not sterilize or clean device as it may affect the performance of the device and cause damage to the monitor.
 - Check sample line regularly during use to ensure that no kinks are present. Kinked tubing may cause inaccurate EtCO2 readings.
 - EtCO2 values, waveforms and duration of use can be affected by certain ambient environment conditions and certain patient conditions.

RESIDUAL RISKS
Refer to warnings and cautions.

- INSTRUCTIONS FOR USE**
Follow instructions to avoid risk of infection and contamination.
1. Connect the gas sampling line to the capnography monitor.
 2. Connect the oxygen tubing to oxygen source.
 3. Set prescribed oxygen flow rate.
 4. Check for flow from prongs.
 5. Insert the prongs into the patient's nares.
 6. Wrap the headset loop up and over ears.
 7. Squeeze the sides of the bolo and glide up under the chin.
 8. Check the monitor display for EtCO2 value and/or waveform.

- MONITORING AND MAINTENANCE**
- Replace the cannula when a blockage or an occlusion alarm signals on the monitor.
 - If there is no waveform or value displayed, check to see if the sample line is disconnected or kinked. If the sample line is disconnected or kinked replace the cannula.
 - Check sampling line regularly during use to ensure that no kinks are present. Kinked tubing may cause inaccurate EtCO2 readings or affect O2 delivery to patient.

INSTRUCTIONS FOR USE IN COMBINATION WITH ACCESSORIES / OR OTHER DEVICES
The sampling line connects to capnography monitor. The oxygen tubing connects to an oxygen gas source.

AirLife tubing connectors are designed to work with gas outlets that comply to ISO 5359:2014. This includes, but is not limited to, stand-alone and portable oxygen concentrators, oxygen tanks, oxygen flowmeters, and bubble humidifiers. The ribbed style end connector is compatible with a 5 mm - 7 mm tapered oxygen fitting or ¼ inch (6.35 mm) tapered oxygen fitting (Figure A). The style of the tapered fitting will vary depending on the oxygen gas source or oxygen device.

1. Place the ribbed end connector between your thumb and forefinger (Figure B).
2. Align the end of the connector to tapered fitting, then push the end connector over the tapered fitting until snug (Figure B).
3. Gently tug on the tubing to ensure a secure fit.


Fig. A

Fig. B

CLEANING / REUSE INSTRUCTIONS
A discontinued use of this product on the same patient, with periods of use separated by a period of non-use, is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.

SAFE DISPOSAL
Decontaminate and dispose of all potentially biohazardous material. Dispose of device in accordance with local, state or national regulations.

INCIDENT REPORTING
Notice to the user that any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myairlife.com, and the competent authority of the Member State in which the user and/or patient is established.

PERFORMANCE CHARACTERISTICS-SPECIFICATIONS	
Patient Population	Adult
EtCO2 Sampling Line (I.D.)	0.050" (1.27 mm)
EtCO2 Sampling Rate	< 0.5 LPM
EtCO2 End Connector	Male Luer / Female Luer
EtCO2 Tubing Length	7' (2.1m) 10' (3.1m) 14' (4.3m)
O2 Tubing Length	7' (2.1m) 10' (3.1m) 14' (4.3m)

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MT Promedt Consulting Ltd.
MedEnvoy UK Limited