EN

English

Salter Labs™ In-Line Small Volume Nebulizer

DEVICE DESCRIPTION:
The device is non-sterile, disposable and for single-patient, multiple-use. A pneumatic small volume nebulizer that fits into an adult, 22 mm breathing circuit to deliver medicated aerosols to the patient. The anti-drool tee-assembly has a 22mm O.D. X 22mm O.D. X 18mm I.D. The recommended driving gases are air, oxygen, or air-oxygen mixture. The nebulizer may be used with a gas source capable of delivering a flow rate of 4 to 8 LPM at a nominal operating pressure between 25 to 52 psi (172.7-358.5 kPa; 1723.7- 3585.3 hPa.

INTENDED PURPOSE:
The nebulizer is a disposable device which is used where liquids are to be delivered to a patient in an aerosol form. It is for single patient multiple-use.

INDICATIONS FOR USE:
The SVN is indicated for patients requiring aerosolized respiratory medication to treat a variety of breathing conditions. This includes but is not limited to the treatment of asthma, chronic obstructive pulmonary disorder (COPD), cystic fibrosis and other illnesses that impact airway function.

PATIENT TARGET GROUP:
Pediatric and adult patients

INTENDED USERS:
Trained medical professionals.

ENVIRONMENT OF USE:
Hospitals, sub-acute, surgical centers, skilled nursing facilities.

EXPECTED CLINICAL BENEFITS:
Small volume nebulizer aerosolized medication allows for quicker and easier absorption of medication into the lungs. Aerosolized bronchodilators help to relieve bronchospasms, decrease the frequency of exacerbations and improve lung function, therefore, relieving shortness of breath, and improving health status and exercise tolerance. Other aerosolized medications help to reduce sputum viscosity, decrease airway inflammation and improve pulmonary function.

CONTRAINDICATIONS:
None known

- ⚠️ WARNINGS:**
- Prior to use in the home, ensure patient and/or caregiver have received proper training on the use and cleaning of the nebulizer.
 - Position tubing to avoid strangulation.
 - Device contains components that may present choking hazard.
 - Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
 - To reduce the risk of misconnections and patient injury, always trace tubing from gas source to the medical device before connecting
 - The safety and effectiveness of the nebulizer for unapproved use(s) has not been established.
 - The nebulizer is intended to be use with air or oxygen gas source, it may not be compatible with other gases that haven't been evaluated.
 - To reduce the risk of fire and burns, do not use near flame or heat source when oxygen is in use.
 - To reduce the risk of fire and burns, do not smoke or vape when oxygen is in use.

- ⚠️ CAUTIONS:**
- Clinician or caregiver should avoid inhaling aerosolized medication.
 - This is a single-patient multi-use product, reuse on more than one patient may cause a risk of cross-contamination.
 - Do not attempt to sterilize or autoclave this device.
 - The measurement accuracy, system performance, or a malfunction may occur as a result of the product being physically damaged due to disinfection and/or re-sterilization.

- Follow cleaning instructions to avoid risk of infection and contamination.
- Prior to cleaning, during dismantling and reassembling, and after each use, remove the nebulizer from the driving gas source.

RESIDUAL RISKS:
Refer to Warning and Cautions.

NECESSARY USER QUALIFICATIONS:
The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings, and cautions as stated in the IFU.

Rx Only: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

- INSTRUCTIONS:**
Always follow your facility's policy and procedure for the delivery of aerosol treatments and the cleaning of equipment.
1. Using aseptic technique, remove nebulizer from polybag.
 2. Tighten cap securely to the medication cup.
 3. Attach the tubing to the bottom port of the medication cup and the other end to the gas source.
 4. Ensure the green cone is properly seated on the bottom of the medication cup.
 5. Add prescribed medication. Maximum fill volume is 5 mL.
 6. Attach nebulizer to the bottom of the tee-assembly.
 7. Follow your facility's protocol, and insert the tee-assembly into the breathing circuit on the inspiratory side.
 8. Set gas flow rate between 5 L/min and 8 L/min.
 9. Keep nebulizer in upright position. Check for leaks. Verify aerosolization of the medication.
 10. Adjust ventilator settings as needed to compensate for the additional flow.
 11. At the end of the treatment, turn off flow and readjust the ventilator settings.
 12. Remove nebulizer from breathing circuit per your facility's protocol.

CLEANING/REUSE INSTRUCTIONS:
Replace the nebulizer every 7 days, or after 21 cleaning cycles, whichever is first. To Clean: Disconnect nebulizer from gas source prior to cleaning. After each use, thoroughly rinse the nebulizer with sterilized water. DO NOT submerge the supply tubing in water. Disassemble all nebulizer components. Rinse with sterilized water, shake out any excess. Allow the device to dry, re-assemble and place in a dry, dust free location. Or, follow your facilities approved cleaning process for nebulizers. A discontinued use of this product on the same patient, with periods of use separated by a period of non-use, is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.

INSTRUCTIONS FOR USE IN COMBINATION WITH ACCESSORIES / OR OTHER DEVICES:
The tubing end connector is designed to work with gas outlets that comply to ISO 5359:2014. This includes, but is not limited to oxygen tanks, flowmeters, oxygen blenders and air compressors. The tubing trumpet end connector is compatible with a 5 mm - 7 mm tapered oxygen fitting or ¼ inch (6.35 mm) tapered oxygen fitting. The nebulizer may also be used with an aerosol mask or a mouthpiece, not included in the packaging.


SAFE DISPOSAL:
Decontaminate and dispose of all potentially biohazardous material. Dispose of device in accordance with local, state or national regulations.

INCIDENT REPORTING:
Notice to the user that any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myairlife.com, and the competent authority of the Member State in which the user and/or patient is established.

PERFORMANCE CHARACTERISTICS:
Nebulizer performance characteristics stated in the instructions for use are based upon testing that utilizes adult ventilatory patterns and are likely to be different for pediatric or infant populations. The liquid medication for nebulization, in particular a suspension and/or high-viscosity solution, can alter the particle size distribution curve, the mass median aerodynamic diameter (MMAD), aerosol output, and/or aerosol output rate, which can then be different from what has been disclosed. All stated performance characteristics were derived from the nebulization of 4mL of albuterol 0.1% (M/V) concentration in 0.9% sodium chloride solution. The measured aerosol output/output rate are not intended to determine correct dosage and these outputs can differ when the nebulizer is used in-line with a ventilator.

NEBULIZER IN-LINE SVN			
Test Solution: 4mL of albuterol 0.1% (M/V) concentration in 0.9% sodium chloride solution.			
Performance Characteristic	Source Pressure: 50 psi		
	4 LPM	6 LPM	8 LPM
Aerosol Output Rate (ml/min)	0.18	0.27	0.31
Aerosol Output Rate (µg/min)	181	270	311
Aerosol Output (ml)	0.81	0.92	0.97
Aerosol Output (µg)	807	916	967
Percentage of Fill Volume Emitted	20.2%	22.9%	24.2%
Percentage of Fill Volume Emitted per Minute	4.5%	6.7%	7.8%
Nebulization Time (min)	14.0	7.83	6.25
Residual Volume (ml)	1.67	1.68	1.51
Aerosol Output Fraction	0.81	0.92	0.97
Particle Size Data at Source Pressure of 50 psi at 8 LPM			
Mass Median Aerodynamic Diameter (MMAD)			3.73 (µm)
Geometric Standard Deviation (GSD)			2.87

Stage No	Effective Cut-Off Dia. (µm)	Nebulizer in-Line SVN Cumulative % particle mass of albuterol under size			
		Mean	Sample 1	Sample 2	Sample 3
8	0.00	0.0%	0.0%	0.0%	0.0%
7	0.68	8.5%	8.1%	8.8%	9.1%
6	0.90	13.1%	13.5%	11.9%	13.9%
5	1.22	22.9%	24.3%	19.4%	25.0%
4	2.87	40.7%	44.2%	36.9%	40.9%
3	4.85	60.7%	64.0%	57.5%	60.4%
2	7.50	74.0%	76.9%	70.3%	74.8%
1	8.27	79.0%	83.4%	74.2%	79.5%

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