
8900 Series T-Piece Small Volume Nebulizer

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7.

Tubagem de ar-oxigénio

EN 8900 Series T-Piece Small Volume Nebulizer
DEVICE DESCRIPTION: The device is non-sterile, disposable and for single-patient multiple-use. This is a hand-held pneumatic small volume nebulizer (SVN) which connects to a mouthpiece or aerosol mask. The recommended driving gases are air, oxygen, or air-oxygen mixture. The device may be used with an air compressor or a gas source capable of delivering a flow rate of 4 to 8 LPM at a nominal operating pressure between 25 to 52 psi (172.7-358.5 kPa; 1723.7-3585.3 hPa).

INTENDED PURPOSE: The Nebulizer is intended for use where liquid (solutions, suspensions, and emulsions) medications are to be delivered to a patient in aerosol form.

INDICATIONS FOR USE: The SVN is indicated for patients requiring aerosolized respiratory medication to treat a variety of breathing conditions. This includes but is not limited to the treatment of asthma, chronic obstructive pulmonary disorder (COPD), cystic fibrosis and other illnesses that impact airway function.

PATIENT TARGET GROUP: For spontaneously breathing infant, pediatric and adult patients.

ENVIRONMENT OF USE: Hospitals, sub-acute, medical clinics, prehospital, home, surgical centers, skilled nursing facilities.

EXPECTED CLINICAL BENEFITS: Small volume nebulizer aerosolized medication allows for quicker and easier absorption of medication into the lungs. Aerosolized bronchodilators help to relieve bronchospasms, decrease the frequency of exacerbations and improve lung function, therefore, relieving shortness of breath, and improving health status and exercise tolerance. Other aerosolized medications help to reduce sputum viscosity, decrease airway inflammation and improve pulmonary function.

CONTRAINDICATIONS: None known


- WARNING:**
- Prior to use in the home, ensure patient and/or caregiver have received proper training on the use and cleaning of the nebulizer.
 - Position tubing to avoid strangulation.
 - Device contains components that may present choking hazard.
 - Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
 - To reduce the risk of misconnections and patient injury, always trace tubing from gas source to the medical device before connecting.
 - The safety and effectiveness of the nebulizer for unapproved use(s) has not been established.
 - The nebulizer is intended to be used with air or oxygen gas source, it may not be compatible with other gases that haven't been evaluated.
 - To reduce the risk of fire and burns, do not use near flame or heat source when oxygen is in use.
 - To reduce the risk of fire and burns, do not smoke or vape when oxygen is in use.

- CAUTIONS:**
- Clinician or caregiver should avoid inhaling aerosolized medication.
 - This is a single-patient multi-use product, reuse on more than one patient may cause a risk of cross-contamination.
 - Do not attempt to sterilize, autoclave or clean this device. The measurement accuracy, system performance, or a malfunction may occur as a result of the product being physically damaged due to cleaning, disinfection, re-sterilization, or reuse.
 - A discontinued use of this product on the same patient, with periods of use separated by a period of non-use, is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.
 - Follow cleaning instructions to avoid risk of infection and contamination.
 - Prior to cleaning, during dismantling and reassembling, and after each use, remove the nebulizer from the driving gas source.
 - Disconnect nebulizer and tubing from gas source prior to dismantling and reassembling device.

RESIDUAL RISKS: Refer to Warning and Caution statements.
NECESSARY USER QUALIFICATIONS: This device is restricted to sale by or on the order of a physician. The end-user should receive training on the proper use of the device and familiar with the indications for use, warnings and cautions as stated in the instructions for use. There are no additional requirements for special facilities, special training, or particular qualifications for the use of this device. RX-Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

- INSTRUCTIONS:** Always follow your facility's policy and procedure for the delivery of aerosol treatments and the cleaning of equipment.
1. Using aseptic technique, remove nebulizer from polybag.
 2. Tighten cap securely to the medication cup.
 3. Attach the tubing to the bottom port of the medication cup and the other end to the gas source.
 4. Ensure the green cone is properly seated on the bottom of the medication cup.
 5. Add prescribed medication. Maximum fill volume is 5 mL.
 6. **a.)** For use with a mouthpiece. Attach nebulizer to bottom of tee-assembly. Insert mouthpiece into large port of tee-adaptor. If applicable, attach reservoir tube to other side. Instruct patient to keep lips sealed around mouthpiece during treatment

- EN**
7. **b.)** For use with an aerosol mask, insert the mask adaptor into the top the nebulizer cap. The tee-assembly is not required.
 7. Set gas flow rate between 5 L/min and 8 L/min or turn on air compressor.
 8. Check for leaks and verify aerosolization of the medication.
 9. Keep nebulizer at a 45 degree angle or greater for optimal aerosolization.
 10. At the end of the treatment, turn off low and disconnect from gas source.

CLEANING/REUSE INSTRUCTIONS: Replace the nebulizer every 7 days, or after 21 cleaning cycles, whichever is first. To Clean: Disconnect nebulizer from gas source prior to cleaning. After each use, thoroughly rinse the nebulizer with sterilized water. DO NOT submerge the supply tubing in water. Disassemble all nebulizer components. Rinse with sterilized water, shake out any excess. Allow the device to dry, re-assemble and place in a dry, dust free location. Or, follow your facilities approved cleaning process for nebulizers.

INSTRUCTIONS FOR USE IN COMBINATION WITH ACCESSORIES / OR OTHER DEVICES: The tubing end connector is designed to work with gas outlets that comply to ISO 5359:2014. This includes, but is not limited to oxygen tanks, flowmeters, oxygen blenders and air compressors. The tubing trumpet end connector is compatible with a 5 mm - 7 mm tapered oxygen fitting or ¼ inch (6.35 mm) tapered oxygen fitting.

STORAGE INSTRUCTIONS / CONDITIONS: Recommended normal storage temperature -30°C/-22°F to +60°C/+140°F. No other special storage conditions required.

SAFE DISPOSAL: Dispose of device in accordance with local, state or national regulations. Decontaminate and dispose of all potentially biohazardous material.


INCIDENT REPORTING: Notice to the user and/or patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. If the patient's physical symptoms worsen or experiences a sudden change in their condition, contact the physician.

PERFORMANCE CHARACTERISTICS: Nebulizer performance characteristics stated in the instructions for use are based upon testing that utilizes adult ventilatory patterns and are likely to be different for pediatric or infant populations. The liquid medication for nebulization, in particular a suspension and/or high-viscosity solution, can alter the particle size distribution curve, the mass median aerodynamic diameter (MMAD), aerosol output, and/or aerosol output rate, which can then be different from what has been disclosed. All stated performance characteristics were derived from the nebulization of 4mL of albuterol 0.1% (M/V) concentration in 0.9% sodium chloride solution. The measured aerosol output/output rate are not intended to determine correct dosage and these outputs can differ when the nebulizer is used in-line with a ventilator.

Nebulizer In-Line SVN			
Test Solution: 4mL of albuterol 0.1% (M/V) concentration in 0.9% sodium chloride solution.			
Performance Characteristics	Source Pressure: 50 psi		
	4 LPM	6 LPM	8 LPM
Aerosol Output Rate (µg/min)	0.18	0.27	0.31
Aerosol Output Rate (µg/min)	181	270	311
Aerosol Output (ml)	0.81	0.92	0.97

Aerosol Output (µg)	807	916	967
Percentage of Fill Volume Emitted	20.2%	22.9%	24.2%
Percentage of Fill Volume Emitted per Minute	4.5%	6.7%	7.8%
Nebulization Time (min)	14.0	7.83	6.25
Residual Volume (ml)	1.67	1.68	1.51
Aerosol Output Fraction	0.81	0.92	0.97
Particle Size Data at Source Pressure of 50 psi at 8 LPM			
Mass Median Aerodynamic Diameter (MMAD)			3.73 (µm)
Geometric Standard Deviation (GSD)			2.87

Stage No	Effective Cut-Off Dia. (µm)	Nebulizerin-Line SVN Cumulative % particle mass of albuterol under size			
		Mean	Sample 1	Sample 2	Sample 3
8	0.00	0.0%	0.0%	0.0%	0.0%
7	0.68	8.5%	8.1%	8.8%	9.1%
6	0.90	13.1%	13.5%	11.9%	13.9%
5	1.22	22.9%	24.3%	19.4%	25.0%
4	2.87	40.7%	44.2%	36.9%	40.9%
3	4.85	60.7%	64.0%	57.5%	60.4%
2	7.50	74.0%	76.9%	70.3%	74.8%
1	8.27	79.0%	83.4%	74.2%	79.5%

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MT Promedt Consulting GmbH
MT Promedt Consulting Ltd.
MedEnvoy UK Limited