EN

English

Salter Labs™ Oxygen Tubing

DEVICE DESCRIPTION:
The device is non-sterile, disposable and for single-patient multi-use. Oxygen tubing is a conduit to deliver oxygen from the source gas to the patient.

INTENDED PURPOSE:
Oxygen tubing is a conduit to deliver oxygen from the gas source to the patient.

INDICATIONS FOR USE:
For patients prescribed supplemental oxygen, used to extend the length of the oxygen tubing on the patient interface.

ENVIRONMENT:
Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

PATIENT TARGET GROUP:
Infant, pediatric, adult. Spontaneously breathing individuals with a variety of breathing conditions requiring supplemental oxygen.

INTENDED USERS :
Trained medical professionals.

- EXPECTED CLINICAL BENEFITS:**
- The tubing delivers up to 15 LPM of supplemental oxygen from the gas source to the patient's interface, to treat or prevent hypoxemia in acute and chronic respiratory conditions. It is monitored by improvement in SpO₂ values, PaO₂ values, vital signs, and/or work of breathing.
 - The 3-channel tubing helps to reduce kinking and allows for greater than 75% of initial flow to be delivered, should the tubing become kinked.
 - Oxygen tubing lengthens the nasal cannula or oxygen mask supply tubing up to 50 feet, increasing patient's mobility and improve health-related quality of life.
 - The green oxygen tubing is 3 times more visible than clear tubing. This makes it easier to see on the floor and may reduce tripping hazard.
 - Supplemental oxygen helps reduce COPD complications by stabilizing pulmonary hypertension, reducing secondary polycythemia, or decreasing arrhythmias.
 - Supplemental oxygen can decrease anxiety, improve sleep, improve mental alertness, and/or improve health-related quality of life in people with chronic respiratory disease.
 - Using oxygen during exercise can improve endurance, heighten performance, and/or decrease the sensation of breathlessness. This helps patients complete physical activities, improve ability to walk and/or perform other physical activities for an extended period.

CONTRAINDICATIONS:
None Known.

- ⚠️ WARNINGS:**
- If oxygen is in use, do not use near flame or heat source.
 - Position tubing to prevent patient from becoming tangled in tubing or strangulation.
 - Patient may become hypoxic if oxygen flow is interrupted.
 - Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
 - To reduce the risk of misconnections and patient injury, always trace tubing from gas source to the medical device before connecting.
 - High levels of oxygen may cause oxygen-induced hypoventilation, oxygen-induced hypercapnia, or oxygen toxicity. In premature infants, it may also cause retrolental fibroplasia.
 - To reduce the risk of fire during facial surgery, stop supplemental oxygen ≥1 minute before and during the use of electrocautery devices or lasers. Consider using a nasopharyngeal airway to deliver supplemental oxygen or a midfacial seal drape to a achieve a barrier between the excess oxygen and surgical field.

- ⚠️ CAUTIONS:**
- Flow rates above 15 LPM may cause increased resistance or back pressure in the tubing.

- To prevent interruption of therapy or injury, do not let children play with the oxygen tubing.
- To prevent excessive moisture accumulation in the tubing, do not use with heated humidification.
- To prevent equipment damage and/or interruption of therapy, do not let pets play with the oxygen tubing.
- Do not place tubing under rugs or other objects as it may obstruct flow.
- To prevent becoming tangled in tubing, position the tubing in front of you before turning to sit or stand.
- This is a single-patient multi-use product, reuse on more than one patient may cause a risk of cross-contamination.

RESIDUAL RISKS:
Refer to Warnings and Cautions.

NECESSARY USER QUALIFICATIONS:
The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings, and cautions as stated in the IFU.

Rx Only: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

- INSTRUCTIONS:**
Follow instructions to avoid risk of infection and contamination.
- Check for leaks at connections and for proper placement prior to use.
 - This product is disposable and is not intended to be cleaned or for prolonged use (≥ 3 months).
 - Discard and replace if the tubing becomes soiled, discolored, or stiff, or at 3 months.
 - Discard and replace the device if patient has had a respiratory infection.
 - Use water trap if excess condensation is present.

- INSTRUCTIONS FOR USE IN COMBINATION WITH ACCESSORIES / OR OTHER DEVICES:**
Oxygen Tubing is for use within an oxygen delivery system. One end of the Oxygen Tubing connects to the patient interface device (nasal cannula, oxygen mask, nebulizer, or other breathing product) and the other end of the Oxygen Tubing connects to the gas source. Additional devices or accessories may be used within the oxygen delivery system.
1. Connect to respiratory device.
 2. Connect to oxygen source.
 3. Check for proper placement prior to turning on oxygen flow.
 4. Check for leaks at connections.

SAFE DISPOSAL:
Decontaminate and dispose of all potentially biohazardous material. Dispose of device in accordance with local, state or national regulations.

CLEANING / REUSE INSTRUCTIONS:
A discontinued use of this product on the same patient, with periods of use separated by a period of non-use, is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.

INCIDENT REPORTING:
Notice to the user that any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myairlife.com, and the competent authority of the Member State in which the user and/or patient is established. If the patient's physical symptoms worsen or experience a sudden change in their condition (e.g., increased shortness of breath, fever, dizziness) call the physician.

- PERFORMANCE CHARACTERISTICS:**
- Connects to or fits all 5mm to 7mm connectors and accessories.
 - Flow rate: 0-15 LPM.
 - The 3-channel tubing helps reduce kinking and allows for greater than 75% of initial flow to be delivered should the tubing become kinked.

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MedEnvoy UK Limited