Oxygen Cannula- Piggyback Demand

EN Oxygen Cannula- Piggyback Demand

DEVICE DESCRIPTION: Demand Nasal Cannula. The device is non-sterile, disposable and for single patient use. The Demand style has the divide between the nostrils and the Modified Demand style has the divide horizontal so that both nostrils sense and deliver. The cannulas have a dual lumen tubing and standard push on connectors.

INTENDED PURPOSE: A nasal oxygen cannula is a two-pronged device used to administer oxygen to a patient through both nostrils. Depending on the cannula style flow rates will vary from > 0 LPM to ≤ 6 LPM.

INDICATIONS: Used with portable oxygen conserving device to sense breathing patern and delivery oxygen only during the inspiratory phase of the breath.

ENVIRONMENT: Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

PATIENT TARGET GROUP: The intended population includes pediatric and adult patients.

Spontaneously breathing individuals with a variety of breathing conditions requiring supplemental oxygen.

EXPECTED CLINICAL BENEFITS:

- Nasal cannulas are a conduit that deliver supplemental oxygen to a patient's nasopharynx to treat or prevent hypoxemia in acute and chronic respiratory conditions. It is monitored by improvement in SpO₂ values, PaO₂ values, vital signs, and/or work of breathing.
- Supplemental oxygen helps reduce COPD complications by stabilizing pulmonary hypertension, reducing secondary polycythemia, and decreasing arrythmias.
- Supplemental oxygen can decrease anxiety, improve sleep, improve mental alertness, and/or improve health-related quality of life in people with chronic respiratory disease.
- Using oxygen during exercise can improve endurance, heighten performance, and/or decrease the sensation of breathlessness. This helps patients complete physical activities, improve ability to walk and/or perform other physical activities for an extended period.
- The 3-channel tubing helps reduce kinking and allows for greater than 75% of initial flow to be delivered should the tubing become kinked.

CONTRAINDICATIONS: None known.

WARNINGS:

- Patient may become hypoxic of oxygen flow is interrupted.
- If oxygen is in use, do not use near flame or heat source.
- Position tubing to prevent patient from becoming tangled in tubing or strangulation.
- Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
- To reduce the risk of misconnections and patient injury, always trace tubing and lines back to their point of origin to verify that correct connections are made.
- High levels of oxygen may cause oxygen-induced hypoventilation, oxygen-induced hypercapnia, or oxygen toxicity. In premature infants, it may also cause retrolental fibroplasia.
- The facepiece contains Tris (nonylphenol) phosphite (TNPP) above 0.1% w/w.

CAUTIONS:

- Consult physician if the patient develops an infection, skin irritation, or material sensitiza-tion.
- Flow rates above 6 LPM may cause increased resistance or back pressure in the tubing.
- To prevent interruption of therapy or injury, do not let children play with the oxygen tubing.
- To prevent excess moisture accumulation in the tubing, do not use with heated humidifica-tion.
- To prevent equipment damage and/or interruption of therapy, do not let pets play with the oxygen tubing.
- Do not place tubing under rugs or other objects as it may obstruct flow.
- To prevent from becoming tangled in tubing, position the tubing in front of you before turning to sit or stand.

RESIDUAL RISKS: Refer to Warnings and Cautions.

NECESSARY USER QUALIFICATIONS: This device is restricted to sale by or on the order of a physician. There are no additional requirements for special facilities, special training, or particular qualifications for the use of this device.

INSTRUCTIONS:

- Select the correct size cannula for the patient.
- Connect one tubing end to the sensing port of the demand regulator.
- Connect the second tubing end to the oxygen port of the demand regulator.
- Adjust the flow control knob to the prescribed liter flow and check for flow from prongs.
- Position the nasal cannula with the prongs facing upward and curved toward the face. Insert the two prongs into the nostrils.
- Wrap the headset loop up and over the ears and tuck the loop under the chin, or adjust and secure the headset loop behind the head.
- Squeeze the sides of the bolo and glide the bolo up under the chin, or until snug around the head.
- Leave enough space to fit at least the width of two fingers between the bolo and chin.
- Discard and replace cannula if soiled or damaged. Do not sterilize.

SAFE DISPOSAL: Dispose of device in accordance with local, state or national regulations.

INCIDENT REPORTING: Contact your country’s Competent Authority and Manufacturer to report any serious incident.

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INDICATIONS : Used with portable oxygen conserving device to sense breathing patern and delivery oxygen only during the inspiratory phase of the breath.

ENVIRONNEMENT : Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

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MT Promedt Consulting GmbH
MT Promedt Consulting Ltd.
MedEnvoy UK Limited