Adult Salter Style Nasal Cannula840082

ENAdult Salter Style Nasal Cannula

DEVICE DESCRIPTION: The device is non-sterile, disposable and for single-patient use. An oxygen cannula is a device with two prongs used to deliver supplemental oxygen to the patient’s nostrils.

INTENDED PURPOSE: A nasal oxygen cannula is a two-pronged device used to administer oxygen to a patient through both nostrils. The nasal cannula is intended for use with flows of > 0 LPM to ≤6 LPM.

INDICATIONS: Used with portable oxygen conserving device to sense breathing pattern and delivery oxygen only during the inspiratory phase of the breath.

ENVIRONMENT: Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

MR SAFETY: MRI Safe

PATIENT TARGET GROUP: The intended population is adult and pediatric patients. Spontaneously breathing individuals with a variety of breathing conditions requiring supplemental oxygen.

- EXPECTED CLINICAL BENEFITS:**
- Nasal cannulas are a conduit that deliver supplemental oxygen to a patient’s nasopharynx to treat or prevent hypoxemia in acute and chronic respiratory conditions. It is monitored by improvement in SpO₂ values, PaO₂ values, vital signs, and/or work of breathing.
 - Supplemental oxygen helps reduce COPD complications by stabilizing pulmonary hypertension, reducing secondary polycythemia, and decreasing arrhythmias.
 - Supplemental oxygen can decrease anxiety, improve sleep, improve mental alertness, and/or improve health-related quality of life in people with chronic respiratory disease.
 - Using oxygen during exercise can improve endurance, heighten performance, and/or decrease the sensation of breathlessness. This helps patients complete physical activities, improve ability to walk and/or perform other physical activities for an extended period.
 - The 3-channel tubing helps reduce kinking and allows for greater than 75% of initial flow to be delivered should the tubing become kinked.

CONTRAINDICATIONS: None Known.

- WARNINGS:**
- Patient may become hypoxic if oxygen flow is interrupted.
 - When oxygen is in use, to reduce the risk of fire and burns, do not use near flame or heat source, and do not smoke or vape.
 - If oxygen is in use, do not use near flame or heat source.
 - Position tubing to prevent patient from becoming tangled in tubing or strangulation.
 - Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
 - To reduce the risk of misconnections and patient injury, always trace tubing and lines back to their point of origin to verify that correct connections are made.
 - High levels of oxygen may cause oxygen-induced hypoventilation, oxygen-induced hypercapnia, or oxygen toxicity. In premature infants, it may also cause retrolental fibroplasia.
 - Do not reuse on more than one patient as this may increase risk of contamination and infection.
 - Prior to use in the home, ensure patient and/or caregiver is provided instructions on the proper use of oxygen therapy in the home environment.

- CAUTIONS:**
- Consult physician if the patient develops an infection, skin irritation, or material sensitization.
 - Flow rates above 6 LPM may cause increased resistance or back pressure in the tubing.
 - To prevent interruption of therapy or injury, do not let children play with the oxygen tubing.
 - To prevent excess moisture accumulation in the tubing, do not use with heated humidification.
 - To prevent equipment damage and/or interruption of therapy, do not let pets play with the oxygen tubing.
 - Do not place tubing under rugs or other objects as it may obstruct flow.
 - To prevent from becoming tangled in tubing, position the tubing in front of you before turning to sit or stand.

RESIDUAL RISKS: Refer to Warnings and Cautions.

NECESSARY USER QUALIFICATIONS: This device is restricted to sale by or on the order of a physician. The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings, and cautions as stated in the IFU. RX only: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.


- INSTRUCTIONS:**
- Select the correct size cannula for the patient.
 - Connect to the oxygen source.
 - Adjust the flow control knob to the prescribed liter flow. Check for flow from the prongs.
 - Position the nasal cannula with the prongs facing upward and curved toward the face. Insert the two prongs into the nostrils.
 - Wrap the headset loop up and over the ears and tuck the loop under the chin, or adjust and secure the headset loop behind the head.

- EN
- Squeeze the sides of the bolo and glide the bolo up under the chin, or until snug around the head.
 - Leave enough space to fit at least the width of two fingers between the bolo and chin.
 - Discard and replace cannula if soiled or damaged. Do not sterilize.

CLEANING / REUSE INSTRUCTIONS: Single-patient use, do not sterilize.

INSTRUCTIONS FOR USE IN COMBINATION WITH ACCESSORIES / OR OTHER DEVICES: AirLife tubing connectors are designed to work with gas outlets that comply to ISO 5359:2014. This includes, but is not limited to, stand-alone and portable oxygen concentrators, oxygen tanks, oxygen flowmeters, and bubble humidifiers. The ribbed end connector is compatible with a 5 mm - 7 mm tapered oxygen fitting or ¼ inch (6.35 mm) tapered oxygen fitting (Figure A). The style of the tapered fitting will vary depending on the oxygen gas source or oxygen device.

- Place the ribbed end connector between your thumb and forefinger (Figure B).
 - Align the end of the connector to tapered fitting, then push the end connector over the tapered fitting until snug (Figure B).
 - Gently tug on the tubing to ensure a secure fit.
- The female threaded connector is designed to be used with a male DISS 1240 oxygen fitting, 5 mm - 7 mm tapered oxygen fitting, or a ¼ inch (6.35 mm) tapered oxygen fitting, e.g., nipple & nut adapter (Figure C).
- Place the threaded connector between your thumb and forefinger below the threaded section (Figure D).
 - Position the connector’s end with either the threaded male fitting or the tapered fitting.
 - For attachment to male fitting, firmly push the connector onto the oxygen fitting and turn clockwise until a secure fit is achieved. (Figure D)
 - For attachment to a tapered fitting, with slight pressure, slide the connector onto the gas fitting and turn clockwise until a secure fit is achieved.
 - Gently tug on tubing to ensure secure fit.

SAFE DISPOSAL: Dispose of device in accordance with local, state or national regulations. Decontaminate and dispose of all potentially biohazardous material.

STORAGE INSTRUCTIONS / CONDITIONS: There are no special storage conditions or handling requirements.

INCIDENT REPORTING: Notice to the user and/or patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. If your/patient’s physical symptoms worsen or you/ patient’s experience a sudden change in your condition (e.g., increased shortness of breath, fever, dizziness) call your doctor .

Cannula	1602	1602TG	1600, 1601, 1602, 1606	1600TG, 1601TG	1600TLC	1699TLC
Connector	Ribbed	Thread Grip	Ribbed	Thread Grip	Ribbed	Thread Grip
Flow Rate	≤ 6 LPM		≤ 6 LPM			
Prong Size	0.093”		0.114”			
Patient Population	Pediatric		Adult			
Other	N/A		N/A		Ear Wraps	

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MT Promedt Consulting GmbH
MT Promedt Consulting Ltd.
MedEnvoy UK Limited