EN

English

Salter Labs™ Tender Grip Cannula Fixation Device

Fig. 1

Fig. 2

Fig. 3

DEVICE DESCRIPTION:

The device is non-sterile, disposable and for single-patient use. Tender Grip Cannula Fixation Device is an adhesive device used to secure cannula or tubing to patient's skin. Tender Grip Cannula Fixation Device is a circular adhesive tape with a fixation to secure tubing.

INTENDED PURPOSE:

Tender Grip Cannula Fixation Device is used to hold the nasal cannula or tubing in place.

INDICATIONS FOR USE:

Skin fixation accessories to provide a means of affixing cannula or tubing to the face and body.

ENVIRONMENT:

Hospital, out-patient, transport, home, extended care

PATIENT TARGET GROUP:

Newborns, infant, pediatric, adult

INTENDED USERS:

Trained medical professionals.

NECESSARY USER QUALIFICATIONS:

The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings, and cautions as stated in the IFU.

Rx Only: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

EXPECTED CLINICAL BENEFITS:

Tender Grip Cannula Fixation Device helps lift nasal cannula tubing slightly off the skin on a patient's cheek. This prevents rubbing and friction, which then reduces the risk of skin breakdown or potential pressure ulcers.

Tender Grip Cannula Fixation Device securely maintains the cannula position reducing accidental cannula dislodgement and interruption in therapy.

CONTRAINDICATIONS:

- Not for use with endotracheal tubes.
- Not for use over broken skin or on patients with known sensitivity to device materials.

WARNINGS:

- Do not secure to tracheostomy tube with tape or any method other than friction fit.
- Device contains components that may present choking hazard.
- Do not use to secure endotracheal tubes.

CAUTIONS:

- Consult physician if patient develops infection, skin irritation or material sensitization.
- Do not use if known sensitivity to hydrocolloids.
- If skin irritation is noted, discontinue use.
- Do not use on more than one patient to prevent risks of cross contamination.
- Do not apply device over open sores.
- This product is not intended for prolonged or extended use.

- This single-use product is not designed or validated to be reused. Reuse may cause a risk of cross-contamination, affect the measurement accuracy, system performance, or cause a malfunction as a result of the product being physically damaged due to cleaning, disinfection, resterilization, or reuse.

RESIDUAL RISKS:

Refer to warnings and cautions.

INSTRUCTIONS FOR USE:

- Before application, make sure the patient's skin is clean, dry, and free of oil or moisture to ensure proper adhesion.
- Remove Tender Grip Cannula Fixation Device from the protective liner.
- Center the Tender Grip Cannula Fixation Device on the patient's cheek underneath the cannula tubing with the pull tab pointing upward toward the top of the patient's head. (Fig 1)
- Press down on the adhesive disk to secure it firmly to the skin.
- Use the pull tab to peel the clear fixation strap back toward the ear until you feel a resistance. (Fig 2)
- Position the nasal cannula tubing in the center of the Tender Grip Cannula Fixation Device, between the tan adhesive disk and the clear fixation strap.
- Place the clear fixation strap over the cannula tubing and press down to secure it in place. (Fig 3)
- Replace Tender Grip Cannula Fixation Device every 48 hours or sooner if device becomes soiled, becomes loose or starts to lift from skin.
- To remove the Tender Grip Cannula Fixation Device, use a damp cloth with mild soap or baby oil to slowly peel back the adhesive disk. This will loosen the adhesive and minimize discomfort. Replace Tender Grip Cannula Fixation Device.

SAFE DISPOSAL:

Decontaminate and dispose of all potentially biohazardous material. Dispose of device in accordance with local, state or national regulations.

INCIDENT REPORTING:

Notice to the user that any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myairlife.com, and the competent authority of the Member State in which the user and/or patient is established. If the patient's physical symptoms worsen or experience a sudden change in their condition (e.g., increased shortness of breath, fever, dizziness) call the physician.

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