EN

English

Salter Labs™ In-Line Water Trap

DEVICE DESCRIPTION:
The device is non-sterile, disposable and for single-patient multi-use. A device used to collect excess condensation that may have collected in oxygen tubing.

INTENDED PURPOSE:
In-Line Water Trap is intended to collect excess condensation that may have collected in oxygen tubing.

INDICATIONS FOR USE:
To capture excess condensation within oxygen delivery system.

ENVIRONMENT:
Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

PATIENT TARGET GROUP:
Infant, pediatric, adult. Spontaneously breathing individuals with a variety of breathing conditions requiring supplemental oxygen.

INTENDED USERS
Trained medical professionals.

EXPECTED CLINICAL BENEFITS:
(Indirect) Water trap collects moisture that may accumulate in the tubing, which prevents airway aspiration of the excess moisture.

CONTRAINDICATIONS:
None Known.

- **WARNINGS:**
- If oxygen is in use, do not use near flame or heat source.
 - Device contains components that may present choking hazard.
 - To reduce the risk of misconnections and patient injury, always trace tubing from gas source to the medical device before connecting.
 - Position tubing to prevent patient from becoming tangled in tubing or strangulation.
 - Patient may become hypoxic if oxygen flow is interrupted.
 - Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
 - To reduce the risk of misconnections and patient injury, always trace tubing gas and lines back to their point of origin to verify that correct connections are made.
 - High levels of oxygen may cause oxygen-induced hypoventilation, oxygen-induced hypercapnia, or oxygen toxicity. In premature infants, it may also cause retrolental fibroblastoplasia.
 - To reduce the risk of fire during facial surgery, stop supplemental oxygen ≥1 minute before and during the use of electrocautery devices or lasers.
 - Consider using a nasopharyngeal airway to deliver supplemental oxygen or a midfacial seal drape to achieve a barrier between the excess oxygen and surgical field.

- **CAUTIONS:**
- This is a single-patient multi-use product, reuse on more than one patient may cause a risk of cross-contamination.

RESIDUAL RISKS:
Refer to Warnings and Cautions.

NECESSARY USER QUALIFICATIONS:
The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings, and cautions as stated in the IFU.

Rx Only: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

- INSTRUCTIONS:**
Follow instructions to avoid risk of infection and contamination.
- Check for leaks at connections and for proper placement prior to use.
 - This product is disposable and not intended to be cleaned or for prolonged use (≥3 months).

- Discard and replace if the device becomes soiled, discolored, or cracked, or discard at 3 months.
- Empty the In-line Water Trap at least once a day, or more frequently if the water level reaches end cap tubes. The water level should always be kept below the openings on the end cap tubes.
- To empty the In-line Water Trap, remove one of the end caps and drain water into the sink or container.

INSTRUCTIONS FOR USE IN COMBINATION WITH ACCESSORIES / OR OTHER DEVICES:
The In-Line Water Trap is for use within an oxygen delivery system, between the cannula tubing and any supplemental oxygen tubing. Additional devices or accessories may be used within the oxygen delivery system. Connect one end of the water trap to the supply tubing and the other end to the nasal cannula.

CLEANING / REUSE INSTRUCTIONS:
A discontinued use of this product on the same patient, with periods of use separated by a period of non-use, is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.

SAFE DISPOSAL:
Decontaminate and dispose of all potentially biohazardous material. Dispose of device in accordance with local, state or national regulations.

INCIDENT REPORTING:
Notice to the user that any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myairlife.com, and the competent authority of the Member State in which the user and/or patient is established. If the patient's physical symptoms worsen or experience a sudden change in their condition (e.g., increased shortness of breath, fever, dizziness) call the physician.

PERFORMANCE CHARACTERISTICS:

- Fits Oxygen Tubing with standard 5 mm to 7 mm end connectors.
- Flow rate: 0-15 LPM.

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MT Promedt Consulting Ltd.
MedEnvoy UK Limited