


Whisper Cannula

Nasal Cannula delivers up to 8 LPM of oxygen to the nostrils

Contains 1 adult nasal cannula with 7 ft extension



- Quieter
- Larger nose piece
- Same amount of oxygen

REF OTC -103-7



Who to Contact

- For questions about your oxygen equipment, contact your local home care provider.
- If your physical symptoms worsen or you experience a sudden change in your condition (e.g., increased shortness of breath, fever, dizziness) call your doctor .
- If you experience severe physical problems (e.g., chest pains, increasing shortness of breath), call 911.
- For questions or comments about your Whisper C annula, contact Salter Labs customer service at 800-421-0024, Mon-Fri 8 AM to 5 PM CST; or email customer care@salterlabs.com. You may also write to us at Salter Labs, Customer Service, 30 Spur Drive, El Paso, TX 79906.

Material

This product is manufactured without the use of natural latex rubber or phthalates, such as DEHP.

Storage Temperature:
-4 °F to 122°F

Operating Temperature:
41 °F to 104 °F



Whisper Cannula

REF OTC -103-7



LOT#

US Patent 5,573,280; D369,409
Made in Mexico
El Paso, TX 79906 USA
30 Spur Drive
Salter Labs



Read all the information in the package insert before using this product. If you do not understand the instructions, call Salter Labs customer service at 800-421-0024, Mon-Fri 8 am to 5 pm CST; or email customer care@salterlabs.com.

PRECAUTIONS FOR USE OF NASAL CANNULA WITH OXYGEN: Follow your homecare provider's instructions for safe use and operation for your portable and home oxygen systems. Do not smoke or allow anyone to smoke around you. Do not use petroleum or oil based creams and lotions in and around your nose. Avoid using electrical equipment that may cause a spark. Do not let children or pets play with your oxygen equipment. Do not kink, bend or tie your nasal cannula tubing. Replace the nasal cannula if it becomes damaged. Have an extra nasal cannula for backup. Use oxygen as prescribed by your doctor. Contact your doctor if you experience worsening or sudden change in shortness of breath or chest tightness.

Non-sterile, disposable, single-use: This product is intended for use by one person and should not be shared. Recommended replacing the Whisper Cannula at least every 7 days.

Do not use if seal is broken.

Reimbursement


The direct purchase of this product may not be eligible for reimbursement through Medicare or your health care insurance company.

SAFETY PRECAUTIONS

<ul style="list-style-type: none"> • Oxygen is a nonflammable gas, but does support combustion. Follow your homecare provider's instructions for the care and safe operation of your oxygen delivery system (e.g., oxygen cylinder, oxygen concentrator, liquid oxygen). 	
<ul style="list-style-type: none"> • Do not smoke or allow anyone to smoke around you. 	
<ul style="list-style-type: none"> • Keep oxygen equipment at least 6 feet away from flames or any heat source. 	
<ul style="list-style-type: none"> • Do not use petroleum based products, or oil-based creams and lotions in or around your nose. 	
<ul style="list-style-type: none"> • Do not use flammable products such as aerosol sprays, cleaning products near your oxygen. 	
<ul style="list-style-type: none"> • Avoid using electrical equipment that may cause a spark (i.e., electric razor). 	
<ul style="list-style-type: none"> • Use oxygen as prescribed by your doctor. 	
<ul style="list-style-type: none"> • The total length of your nasal cannula and oxygen tubing should not exceed 57 feet to ensure there is enough pressure to deliver prescribed oxygen. 	
<ul style="list-style-type: none"> • Do not kink, bend or tie your oxygen tubing, 	
<ul style="list-style-type: none"> • Do not set anything on your tubing that may obstruct flow. 	
<ul style="list-style-type: none"> • Keep excess tubing loosely coiled and out of the way to prevent tripping on oxygen tubing. 	
<ul style="list-style-type: none"> • Do not place oxygen tubing or nasal cannula under blankets, bedsheets, rugs, etc. 	
<ul style="list-style-type: none"> • Use caution so your oxygen tubing does not become entangled in your furniture. 	
<ul style="list-style-type: none"> • Keep an extra nasal cannula and other oxygen supplies available for use. 	
<ul style="list-style-type: none"> • Do not let children or pets play with your oxygen equipment. 	
<ul style="list-style-type: none"> • Recommend use of swivel adapter to reduce the possibility of twisting for nasal cannula and supply tubing longer than 14 feet. 	



Product Description



The Whisper Cannula is available with a 4-feet (OTC-103-4) or 7-feet (OTC-103-7) oxygen extension tubing and standard end connector.


Intended Use

The Whisper Cannula is used to deliver up to 8 LPM of supplemental oxygen to adults who have a prescription for home oxygen therapy and have received instruction from a trained healthcare professional. The Nasal Cannula can be used in the home, outpatient, extended care, transport and hospital environments.


Instructions for Use

Please read and follow the Instruction for Use prior to using the Whisper Cannula. If you do not understand the instructions, contact Salter Labs.

1. Wash hands. Remove nasal cannula from package.
2. Attach ribbed end connector (A) to oxygen source, e.g., oxygen extension tubing (B) or oxygen flow control outlet (C).




3. Adjust flow control knob to your prescribed liter flow. Check for gas flow from nasal prongs.




4. Position the nasal cannula with the nasal prongs (F) facing upward and curved toward you.


5. Insert the two nasal prongs (F) into your nostrils.




6. Wrap the headset loop (E) up and over both ears.



7. Squeeze the sides of the bolo and glide the bolo (D) up under your chin.



8. Leave enough space to fit at least two fingers between the bolo (D) and chin.



9. Alternative placement: The headset loop secured behind your head.



Nasal Cannula Troubleshooting Guide

Problem	Possible Cause	Corrective Action
No oxygen flow from nasal prongs	<ol style="list-style-type: none"> 1. Cannot feel the airflow in your nostrils. 2. Flow control valve is not turned on. 3. Oxygen system is not functioning properly or oxygen container is empty. 4. The nasal cannula is disconnected from oxygen device or supply tubing. 5. Nasal cannula or oxygen tubing kinked or blocked. 	<ol style="list-style-type: none"> 1. Place nasal prongs into a small container of clean water. Bubbles will appear if there is oxygen flow. 2. Set flow control to prescribed setting. 3. Switch to backup oxygen source and contact your homecare provider. 4. Reconnect oxygen tubing. Ensure all tubing connections are tight and secure. 5. Inspect cannula and oxygen tubing for kinks or damage.
Water in nasal cannula or oxygen supply tubing	<ol style="list-style-type: none"> 1. Humidifier bottle overfilled. 	<ol style="list-style-type: none"> 1. Pour out the excess water. Consider adding a water trap to your oxygen supply tubing.
Nasal dryness or irritation	<ol style="list-style-type: none"> 1. Gas flow is dry. 2. No humidifier is being used. 	<ol style="list-style-type: none"> 1. Use normal saline spray or water soluble ointment, (i.e., AYRSaline Nasal Gel) to moisten the inside your nostrils. If condition worsens, contact your doctor. 2. Contact your doctor or homecare provider to request humidification.
Soreness or irritation around ears	<ol style="list-style-type: none"> 1. Headset tubing too tight. 2. Tubing pressing against skin. 	<ol style="list-style-type: none"> 1. Loosen headset tubing. 2. Place a cotton padding or cushion (i.e., E-Z Wrap) between the tubing and the sore area.
Skin rash and/or sores caused by nasal cannula	<ol style="list-style-type: none"> 1. Sensitivity or reaction to nasal cannula material. 	<ol style="list-style-type: none"> 1. Contact your health care provider and/or doctor.

Care and Cleaning Instructions

1. Once a day, use a clean damp cloth to wipe off the nasal prongs and headset tubing.
2. Do not use strong detergents, alcohol based products or lotion soaps.
3. Do not sterilize your nasal cannula.
4. Discard and replace your nasal cannula every 7 days, or more frequently if your cannula becomes damaged.