EN

English

Salter Labs™ Oxygen Tubing Connectors

DEVICE DESCRIPTION:

- The Nipple Nut connector is a tapered and ribbed, one piece, nipple nut assembly for threaded DISS gas fittings.
 - The Swivel/Static connector has ribbed body design and male conical tip.
 - The Y-connector has an open ribbed male connector, and the other two ends are bonded to green, 3-channel tubing.
- These devices are non-sterile, disposable, and for single-patient multi-use.

INTENDED PURPOSE:

- Nipple Nut connector: Intended to connect oxygen tubing to a gas source.
- Swivel/Static connector: Intended to join two lengths of oxygen supply tubing.
- Static Connectors are intended to join two lengths of oxygen supply tubing.
- Y-connector: Intended to connect two oxygen sources to one supply tubing.

INDICATIONS FOR USE:

- Nipple Nut connector: Connect oxygen supply tubing to a DISS 1240 fitting.
- Swivel/Static connector: Join two oxygen supply tubings.
- Y-connector: Join two oxygen sources to deliver higher flow rates to patient interface.

ENVIRONMENT:

Hospitals, sub-acute, medical clinics, pre-hospital, home, surgical centers, skilled nursing facilities.

PATIENT TARGET GROUP:

Infant, pediatric, adult. Spontaneously breathing individuals with a variety of breathing conditions requiring supplemental oxygen.

INTENDED USERS :

Trained medical professionals.

EXPECTED CLINICAL BENEFITS:

(Indirect)

- The tubing connectors are part of the oxygen setup to deliver supplemental oxygen from gas outlet to patient interface.
- The Swivel connector has a 360 degree rotation that helps to reduce tubing torque and twisting preventing interruption of flow.
 - The connectors join oxygen tubing to increase patient's mobility and improve health related quality of life.
 - The Y-Connector joins two oxygen sources allowing for higher flow rates to relieve hypoxemia.

CONTRAINDICATIONS:

None Known.

WARNINGS:

- If oxygen is in use, do not use near flame or heat source.
- Device contains components that may present choking hazard.
- To reduce the risk of misconnections and patient injury, always trace tubing from gas source to the medical device before connecting.
- Position tubing to prevent patient from becoming tangled in tubing or strangulation.
- Patient may become hypoxic if oxygen flow is interrupted.
- Keep excess tubing loosely coiled and out of the way to prevent kinking and tripping hazard.
- To reduce the risk of misconnections and patient injury, always trace tubing gas and lines back to their point of origin to verify that correct connections are made.
- High levels of oxygen may cause oxygen-induced hypoventilation, oxygen-induced hypercapnia, or oxygen toxicity. In premature infants, it may also cause retrolental fibroblastasia.
- To reduce the risk of fire during facial surgery, stop supplemental oxygen ≥1 minute before and during the use of electrocautery devices or lasers. Consider using a nasopharyngeal airway to deliver supplemental oxygen or a midfacial seal drape to achieve a barrier between the excess oxygen and surgical field.

CAUTIONS:

- This is a single-patient multi-use product, reuse on more than one patient may cause a risk of cross-contamination.

RESIDUAL RISKS:

Refer to Warnings and Cautions.

NECESSARY USER QUALIFICATIONS:

The end-user should receive training on the proper use of the device and be familiar with the indications for use, warnings, and cautions as stated in the IFU.

Rx Only: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.

INSTRUCTIONS:

- Follow instructions to avoid risk of infection and contamination.
- Check for leaks at connections and for proper placement prior to use.
 - This product is disposable and is not intended to be cleaned or for prolonged use (≥ 3 months).
 - Discard and replace if the tubing becomes soiled, discolored, or stiff, or at 3 months.
 - Discard and replace the device if patient has had a respiratory infection.
 - Use water trap if excess condensation is present.

INSTRUCTIONS FOR USE IN COMBINATION WITH ACCESSORIES / OR OTHER DEVICES:

- Nipple Nut Connector is for use within an oxygen delivery system and connects tubing to a DISS 1240 fitting.
 - Attach the nipple adaptor to oxygen gas source.
 - Turn wing nut until tight.
 - Attach oxygen supply tubing to barbed nipple end.
- Oxygen Tubing Connectors (Swivel/Static) are for use within an oxygen delivery system and connects tubing to tubing.
 - Connect one end of the tubing to the oxygen source.
 - Insert, push, and twist the tapered end of the tubing connector into oxygen tubing end connector.
 - Repeat to attach the other side of the connector to the nasal cannula or oxygen mask tubing.
 - Gently tug on connection to check for a secure fit.
 - To remove, hold the center of the tubing connector, twist and pull on the tubing end connector.
- Y-Connector is used to connect two oxygen sources to one supply tubing.
 - Connect the tubing on the Y-connector to the two oxygen sources.
 - Insert the barb end of the Y-connector into the end of nasal cannula or oxygen mask tubing.
 - Gently pull on the connection to ensure secure fit.
 - Check connections for leaks.
 - To remove, hold the center of the Y-connector while pulling on the nasal cannula or mask tubing.
 - Disconnect the tubing from oxygen sources.

CLEANING / REUSE INSTRUCTIONS:

A discontinued use of this product on the same patient, with periods of use separated by a period of non-use, is acceptable. The user must ensure that the product is not damaged or contaminated between usages on the same patient.

SAFE DISPOSAL:

Decontaminate and dispose of all potentially biohazardous material. Dispose of device in accordance with local, state or national regulations.

INCIDENT REPORTING:

Notice to the user that any serious incident that has occurred in relation to the device should be reported to the manufacturer at productquality@myairlife.com, and the competent authority of the Member State in which the user and/or patient is established. If the patient's physical symptoms worsen or experience a sudden change in their condition (e.g., increased shortness of breath, fever, dizziness) call the physician.

PERFORMANCE CHARACTERISTICS:

- Fits Oxygen Tubing with standard 5mm to 7mm end connectors.
- Flow rate: 0-15 LPM.

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