

SUPERNO₂VA[®] ET MASK AND SYSTEM



The SuperNO₂VA[®] Et mask and system and mask are designed to:


Maintain upper airway patency, decrease incidence of hypoxemia^{1,2}, and consistently and accurately capture end tidal CO₂ from either the mouth or nose.

The SuperNO₂VA[®] Et system includes the mask, 2 L hyperinflation bag, and other accessories including a gas sampling line, delivering positive pressure while providing access to the oral cavity.




DECREASED INCIDENCE OF HYPOXEMIA


Incidence of Hypoxemia¹




Incidence of Hypoxemia²



1% CO₂ Delivered to Subject



5% CO₂ Delivered to Subject



1. Zeping Xu. Comparison of oxygenation and ventilation in patients undergoing colonoscopy during anesthesia using the SuperNO₂VA[™] nasal PAP ventilation device vs. routine care: A prospective randomized trial.
2. Francesca Dimou. Nasal positive pressure with the SuperNO₂VA device decreases sedation-related hypoxemia during pre-bariatric surgery EGD. Surgical Endoscopy 2019;

Nasal cannula and
oxygen face mask

SuperNO₂VA® Et

Confidently maintain airway patency and monitor ventilation throughout the perioperative experience.

SAMPLING ET_{CO}₂ WITH SUPERNO₂VA® ET:

- Continuously captures expired CO₂
- Simultaneously switches from nasal or oral sampling
- Accurately detects hyper/hypoventilation, apnea, and other pathologies

SuperNO₂VA Et accurately and consistently measures EtCO₂ from the mouth and the nose:

- Virtually all CO₂ that was delivered was captured on exhalation
- 0.04% absolute error in capture over a wide range of oxygen flow rates, respirator rates, and tidal volumes

SuperNO₂VA Et tested at various O₂ flow rates, respiratory rates, and tidal volumes
mean ± standard deviation results

INPUT CO ₂	CAPTURE CO ₂	CAPTURED CO ₂ DEVIATION MEAN	STANDARD DEVIATION
1%	0.96 – 1.01%	0.01	± 0.05
5%	4.86 – 5.00%	0.08	± 0.07

The SuperNO₂VA Et mask and system treat upper airway obstruction and optimize both oxygenation and ventilation in sedated and general anesthesia patients. The SuperNO₂VA Et mask and system can sample expired gases from the patient's oral and/or nasal cavities.

SuperNO₂VA® Et Competitive Advantage




	NASAL CANNULA	HIGH-FLOW NASAL CANNULA	ANESTHESIA MASK	NASAL CPAP	SUPERNO ₂ VA ET
Oral Access	●	●	○	●	●
Passive Oxygenation	●	●	●	●	●
High Flow O ₂	○	●	○	○	●
Rescue Ventilation/RC	○	○	●	○	●
Nasal and/or Oral CO ₂ Capture	●	●	○	○	●
Additional Capital Equipment	Not Required	●	Not Required	●	Not Required

PRE-OP


INTRA-OP

POST-OP


- EGD/TEE
- Colonoscopy
- Bronchoscopy
- Fiberoptic Intubation
- Laryngoscopy
- MRI



Used to maintain
O₂ saturation



Deep sedation
and intra-oral



Provides positive
pressure in the PACU

Benefits

- Pre-oxygenation
- Relieve upper airway obstruction
- Assist and monitor ventilation
- Rescue ventilation
- Use during oral procedures
- Use in transport