

Validation of SuperNO₂VA[®] Et Nasal Performance CO₂ Measurement and Pressure Test

INTRODUCTION

Moderate and deep sedation has long been associated with high rates of respiratory complications such as hypoxemia and hypoventilation.¹⁻³ These complications have been shown to be caused by inadequate monitoring and sedation medications which result in upper airway obstruction and central respiratory depression.¹⁻³ In an attempt to mitigate respiratory complications, the measurement of exhaled carbon dioxide has recently become accepted as an improvement over pulse oximetry for the monitoring of ventilation.⁴ In addition, the American Society of Anesthesiologists (ASA) has endorsed end-tidal capnography (EtCO₂) as a standard of care for moderate and deep procedural sedation.⁵ A key reason for capnography being adopted is because it measures ventilation by capturing the concentration of exhaled carbon dioxide in a process known as end tidal CO₂ while pulse oximetry measures only oxygenation. Another reason why pulse oximetry is not the preferred method for monitoring ventilation is because of its inability to detect hypoventilation or apnea.⁶ This is especially true in patients undergoing procedural sedation receiving supplemental oxygen.⁶⁻⁷

Although capnography has been shown to effectively detect hypoventilation and apnea much faster than pulse oximetry, accurate and consistent measurements of the EtCO₂ during minimally invasive procedures under deep sedation have historically been challenging.⁸

One reason for this is because it is open to air, causing atmospheric gases to be entrained and sampled.⁹ A second reason is due to patients receiving supplemental oxygen, which at flows > 5 LPM “wash-out” or dilute the sample of exhaled carbon dioxide resulting in either a falsely low reading or no reading at all.¹⁰

Recent prospective randomized controlled trials (RCTs) reported up to 54% of all patients experience severe hypoxemia secondary to sedation-related upper airway obstruction (UAO) and respiratory depression.¹¹ Although passive oxygenating devices have the ability to provide higher concentrations of oxygen, they are incapable of generating positive pressure to maintain airway patency. Continuous Positive Airway Pressure (CPAP) has been shown to relieve upper airway obstruction by creating a pneumatic stent.¹² However, its utility is limited by the fact that the machines are very large and expensive, and the high flows required to maintain pressure also dilute the EtCO₂ sample.¹³⁻¹⁴ Another promising therapy is the SuperNO₂VA[®], which is a completely sealed nasal positive airway pressure (PAP) device that provides positive pressure to maintain upper airway patency without the use of capital equipment. A recent RCT comparing the SuperNO₂VA vs nasal cannula with capnography during deep sedation showed a significantly higher Minute Ventilation and reduction in the incidence of severe hypoxemia in the SuperNO₂VA group compared to the nasal cannula with capnography.¹⁵ However, a disadvantage is its inability to capture EtCO₂, especially when the patient exhales from their mouth resulting in false apnea alarms.




Figure 1: SuperNO₂VA Et Nasal Mask features a EtCO₂ Hood and EtCO₂ nasal sampling port

The SuperNO₂VA Et (figure 1) is a novel completely sealed nasal positive airway pressure (PAP) device designed to capture EtCO₂ exhaled from either the patient’s mouth or nose and provide positive pressure to maintain upper airway patency. Combining capnography with positive pressure may prove to be a methodology to further improve outcomes in deep sedation as opposed to passive oxygenation techniques with capnography.

The objectives of this study are to validate that the SuperNO₂VA Et can capture EtCO₂ exhaled from the nose and the mouth, provide ≤ 20 cm H₂O of positive pressure, quantify leak rates, and summarize the performance testing of the SuperNO₂VA Et compared to predicate devices.

METHODS

A simulated patient setup was utilized to compare the accuracy of CO₂ measurements within the SuperNO₂VA® Et and predicate device (Smart CapnoLine Plus Adult/Intermediate CO₂, “Oral-Nasal Set”). Specifically, the Device Under Test (DUT) was placed on a face surrogate and a ventilator (Harvard Apparatus Model 613) simulated breathing while flow controllers set O₂ flow rate and CO₂ concentration. Eight combinations of Input CO₂ (1 ± 0.25%; 5 ± 0.5%), breath rate and tidal volume (12 BPM/500 mL; 20 BPM/300 mL), and O₂ flow rates (1;5 LPM) were tested. A CO₂ monitor (Ohmeda 5250) recorded 16s of data and trials were repeated three times for each condition. Absolute and relative errors between CO₂ Max (i.e., maximum CO₂ over 16s trial) and Input CO₂ were quantified for each DUT: Absolute Error=CO₂ Max-DUT-InputCO₂; Relative Error=(CO₂ MaxDUT-InputCO₂)/InputCO₂*100%. CO₂ Max errors were compared to the ISO80601-2-55:2018 accuracy specification (± 0.43%vol + 8% of gas level). Unpaired t-tests were used to compare CO₂ Max errors between the two devices for tests with Input CO₂ of 1% and 5%. Leak rate and ability to hold a positive pressure were also tested for three SuperNO₂VA Et and Full Face Anesthesia Mask (Ventlab VR5100) samples. The DUT was placed on a surrogate face and sealed with 10 lb of force. Leak rate was determined as the minimum O₂ flow rate to maintain 20 cm H₂O.

RESULTS

The SuperNO₂VA Et had lower CO₂ Max error than the Oral-Nasal Set for all eight conditions (Figure 2). For 1% Input CO₂, CO₂ Max errors were -0.12 ± 0.03%vol (-12.2 ± 3.3%, mean ± SD) for the Oral-Nasal Set compared to -0.01 ± 0.02%vol (-1.3 ± 2.2%) for the SuperNO₂VA Et (p=0.0005). For 5% Input CO₂, the Oral-Nasal Set dramatically underestimated CO₂ (error = -0.93 ± 0.16%vol (-18.6 ± 3.2%)) while the SuperNO₂VA Et had error of -0.08±0.06%vol (-1.5 ± 1.2%) (p<0.0001). At 5% InputCO₂, eight out of 12 trials for the Oral-Nasal Set were outside the ISO error bound while all SuperNO₂VA Et trials were within this bound.

Both the SuperNO₂VA Et and anesthesia mask successfully held a pressure of 20 cm H₂O for the 5-min trials. The SuperNO₂VA Et had a leak rate of 2.0 LPM for all three samples while the anesthesia mask had a mean leak rate of 2.7 (range: 2.5-3.0 LPM) (table 1).

CONCLUSIONS

The testing described in this report show that measurements of CO₂ within the SuperNO₂VA Et mask are accurate over a range of respiratory rates, tidal volumes, O₂ flows, and CO₂ concentrations and are within the error bounds specified by the International Organization for Standardization. Measurements of CO₂ from within the SuperNO₂VA Et were significantly more accurate than measurements taken within the predicate device, an oral-nasal sampling set. Furthermore, the SuperNO₂VA Et can maintain a positive pressure of 20 cm H₂O within the mask with a low leak rate of 2.0 LPM.

SuperNO₂VA Et is differentiated from other methods of airway management as it can combat upper airway obstruction without sacrificing end tidal measurements. Additionally, it can provide positive pressure to force the airways open, maintaining airway patency. This may help prevent patients from becoming hypoxemic and improve overall outcomes in this patient population.

1. Yilmaz M, Aydin A, Karasu Z, Et al. Risk factors associated with changes in oxygenation and pulse rate during colonoscopy. *Turk J Gastroenterol* 2002;13(4):203-8. [PMID: 16378306]
2. Qadeer MA, Lopez AR, Dumot JA, Et al. Risk factors for hypoxemia during ambulatory gastrointestinal endoscopy in ASA I-II patients. *Dig Dis Sci* 2009;54:1035-40. [PMID: 19003534 DOI: 10.1007/s10620-008-0452-2]
3. Qadeer MA, Lopez AR, Dumot JA, Et al. Hypoxemia during moderate sedation for gastrointestinal endoscopy: causes and associations. *Digestion* 2011;84:37-45. [PMID: 21304242 DOI: 10.1159/000321621]
4. Lam, Thach, Et al. "Continuous pulse oximetry and capnography monitoring for postoperative respiratory depression and adverse events: a systematic review and meta-analysis." *Anesthesia & Analgesia* 125.6 (2017): 2019-2029
5. Weaver, Joel. "The latest ASA mandate: CO₂ monitoring for moderate and deep sedation." *Anesthesia progress* 58.3 (2011): 111.
6. Burton, John H., Et al. "Does end-tidal carbon dioxide monitoring detect respiratory events prior to current sedation monitoring practices?." *Academic Emergency Medicine* 13.5 (2006): 500-504.
7. Arakawa H, Kaise M, Sumiyama K, Et al. Does pulse oximetry accurately monitor a patient's ventilation during sedated endoscopy under oxygen supplementation. *Singapore Med J* 2013;54(4):212-5. [PMID: 23624448 DOI: 10.11622/smedj.2013075]
8. Frasca, D., Et al. "Comparison of acoustic and impedance methods with mask capnometry to assess respiration rate in obese patients recovering from general Anesthesia." *Anesthesia* 70.1 (2015): 26-31.
9. Morimoto, Keiko, Et al. "Respiratory rate is an inadequate parameter of ventilation in non-intubated sedation." *Odontology* 107.2 (2019): 219-222.
10. Teng, Wei-Nung, Et al. "Oral capnography is more effective than nasal capnography during sedative upper gastrointestinal endoscopy." *Journal of clinical monitoring and computing* 32.2 (2018): 321-326.
11. Mehta, PP Et al. Capnographic monitoring in routine EGD and colonoscopy with moderate sedation: a prospective, randomized, controlled trial. *The American journal of Gastroenterology* 2016;111(3):395-404.
12. Mathru, Mali Et al. Magnetic Resonance Imaging of the Upper Airway Effects of Propofol Anesthesia and Nasal Continuous Positive Airway Pressure in Humans. *The Journal of the American Society of Anesthesiologists* 1996;84(2):273-79
13. Chung F, Nagappa M, Singh M, Et al. CPAP in the perioperative setting: evidence of support. *Chest* 2016;149(2):586-97.
14. Schonhofer, B. and S. Sortor-Leger. Equipment needs for noninvasive mechanical ventilation. *Eur Respir J*. 2002. 20(4): p. 1029-36.
15. Zeping Xu. Comparison of oxygenation and ventilation in patients undergoing colonoscopy during anesthesia using the SuperNO₂VA nasal PAP ventilation device vs. routine care: A prospective randomized trial.




Figure 2: CO₂ Max Error (in % vol) for the eight performance tests for Oral-Nasal Set (orange) and SuperNO₂VA Et (blue). Horizontal green dashed line corresponds to the ISO 80601-2-55:2018 error limit (0.51% and 0.83% for 1% and 5% input CO₂ respectively). Filled circles are individual trials and bars represent mean error across the three trials for each test.

FLOW LEAK RATE (LPM)

SAMPLE NUMBER	FULL FACE ANESTHESIA MASK	SUPERNO ₂ VA ET
1	3.0	2.0
2	2.5	2.0
3	2.5	2.0
Average	2.7	2.0

Table 1: Flow leak rate results for Full Face Anesthesia Mask and SuperNO₂VA Et.