
REF **FN2446**

AirLife™

en **Misty Fast™ Small Volume Nebulizer**

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en **AirLife™ Misty Fast™ Small Volume Nebulizer**

With 7 foot (2.1 m) Female U/Connect-It Crush Resistant Oxygen Tubing, Baffled Mouthpiece, and 6 inch (15 cm) Flextube

Symbols Glossary				
	Catalogue number		U.S. Federal law restricts this device to sales by or on the order of a physician	
	Lot Number		Clean, Ready to Use	
	Consult instructions for use		Quantity	
	Not made with natural rubber latex		Expiration Date	
	Manufacturer		Caution	
	Contains Phthalates		Temperature Limit	
	Do not re-use			

Indications for Use

This device is intended to be used to aerosolize liquid medication into gases that are delivered directly to the patient for breathing. Its use is indicated when a licensed healthcare professional prescribes or administers medical aerosols to a patient using a small volume nebulizer. The patient population includes Adults, Pediatrics, and Infants that are spontaneously breathing. The product is a single patient, non-sterile, disposable, prescriptive device intended to be used in a hospital setting.

Directions for Use

- Unscrew nebulizer bottle from cap and insert medication into bottle, as prescribed by physician. Reattach bottle to cap. Medication may be added to nebulizer from the top of the nebulizer cap without unscrewing the two parts. Twist to tighten the cap to bottle.
- Using a push-twist motion, insert the nebulizer into the base of the baffled mouthpiece.
- Attach flextube to opposite end of mouthpiece.
- Using a push-twist motion, attach the blue rigid connector on the tubing to the port of the nebulizer bottle. Connect the other end, female U/Connect-It to the oxygen or air source.
- Set gas flow as prescribed by a health care professional and confirm aerosol mist is flowing.
- If using oxygen mixtures in excess of 23% as a gas source, please refer to the safe handling instructions provided by your oxygen supplier.
- Disconnect from pneumatic power source after use.

Recommended Nebulizer Set Flow Rate 8 ± 1 L/min				
Drug Description	Aerosol Characteristic ¹	7 L/min	8 L/min	9 L/min
Albuterol Sulfate Dose 3 mg/3 mL (0.1% by volume) or 3000 µg of drug substance	MMAD (µm) ²	3.8 ± 0.1 µ	3.2 ± 0.1 µ	3.0 ± 0.1 µ
	GSD ³	2.1 ± 0.06	2.1 ± 0.04	2.2 ± 0.04
	Particles (<2 µm) ⁴	22 ± 1%	26 ± 2%	30 ± 1%
	Particles (2 to 5 µm) ⁴	38 ± 1%	41 ± 1%	39 ± 1%
	Particles (>5 µm) ⁴	40 ± 1%	33 ± 2%	31 ± 2%
	Respirable Particles (% <5 µm) ⁴	60 ± 1%	67 ± 2%	69 ± 2%
	Aerosol Output Rate	160 ± 5 µg/min	172 ± 6 µg/min	169 ± 5 µg/min
	Total Delivered Mass of Aerosol ⁵	533 ± 12 µg	493 ± 11 µg	453 ± 9 µg
	% Fill Volume Emitted per minute	5.3 ± 0.2%	5.7 ± 0.2%	5.6 ± 0.2%
	% Fill Volume Emitted ⁶	21.1 ± 0.5%	19.7 ± 0.7%	18.2 ± 0.6%
	Residual Volume ⁷	1.4 ± 0.03 mL	1.4 ± 0.03 mL	1.4 ± 0.02 mL
	Driving Gas Pressure	11 psig	26 psig	35 psig

¹ Performance data using Cascade Impactor (NGI) with an extraction flow of 15 L/min and a representative adult breathing patient for 0.1% Albuterol Sulfate solution measured according to BS ISO 27247:2013. Treatment time was defined as onset of audible sputtering plus one minute. Values represent the Mean and 95% confidence interval of the mean.

² MMAD is Mass Median Aerodynamic Diameter.

³ GSD is the Geometric Standard Deviation.

⁴ Values are the mass of drug substance represented as a percentage of the total delivered mass of aerosol.

Recommended Flow Rate: 8 ± 1 liters per minute.

8. **For Use without Flow meter Nipple Adapter**
Attach female U/Connect-It adapter to threads of flow meter by pressing directly onto threads and twisting to seal.
- Nebulizer may be used up to an angle of 60 degrees from vertical.
 - Using a solution, suspension or emulsion different from that recommended by the manufacturer can alter the particle size distribution curve, the mass median aerodynamic diameter (MMAD), aerosol output, and/or aerosol output rate, which can then be different from those disclosed by the manufacturer.
 - The disclosures listed for nebulizer performance are based upon testing that utilizes adult ventilator patterns and are likely to be different from those stated for the pediatric or infant populations.

Warning

- Do not use in-line with breathing circuits.
- O₂ dependent patients must be monitored during use as air entrainment can reduce FiO₂ levels.

Caution

- Device to be used in a hospital setting only.
- Device is not intended to be used as an oxygen delivery device.

Cleaning Instructions: After each use: Unscrew the nebulizer cap and bottle. Remove the one-piece jet by pulling and twisting the jet off the jet stem. Wash all components in warm soapy water and rinse well. Sterile water may also be used in place of soapy water. Shake out/off excess water or sterile water and air dry for 30 minutes. Alternatively, shake out/off excess water or sterile water and wipe the device with a clean lint free cloth until dry. To aid in the drying process, you may also use the compressed air supply to clear the orifice. Reseat the jet by placing the jet over the jet stem and snap into place. Reattach the nebulizer cap and bottle.

Cumulative duration of use not to exceed 24 hours or 50 cleaning cycles.

Storage

Recommended normal storage temperature: 0°C/+32°F to +50°C/+122°F.

Contains Phthalates (DEHP).

Risks and Precautionary Measures related to Phthalates: This instruction pertains to the phthalate symbol marked on the device or its packaging. If this device is used for the treatment of children, or treatment of pregnant or nursing women; please note that the following types of procedures may increase the risk of exposure to phthalates: Exchange transfusion in neonates, total parenteral nutrition in neonates, multiple procedures in sick neonates, haemodialysis in peripuberal males, male foetus and male infant of pregnant women, and lactating women; and massive blood infusion into trauma patients. Although these procedures have the potential for increased risk of exposure, conclusive evidence of human health risks has not been established. As a precautionary measure, to reduce the potential for unnecessary exposures to phthalates, the product must be used in accordance with the instructions for use, and practitioners should refrain from using this product beyond the period of time the product is medically necessary or needed.

Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

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en **Aerosol Characterization using Cascade Impaction (NGI) for three common aerosolized medications¹**

Recommended Nebulizer Set Flow Rate 8 ± 1 L/min						
Drug Description	Aerosol Characteristic	7 L/min		8 L/min		9 L/min
Albuterol Sulfate Typical dose 2.5 mg/3 mL (0.083% by volume) or 2500 µg of drug substance.	MMAD (µm) ²	3.5 ± 0.1 µ		3.3 ± 0.1 µ		3.0 ± 0.1 µ
	GSD ³	2.1 ± 0.02		2.2 ± 0.04		2.2 ± 0.02
	Extra-fine particles (< 1 µm) ⁴	8 ± 1%	100 ± 4 µg	10 ± 1%	105 ± 8 µg	10 ± 1% 92 ± 5 µg
	Fine particles (< 5 µm) ⁴	64 ± 1%	757 ± 23 µg	65 ± 1%	696 ± 23 µg	68 ± 1% 596 ± 17 µg
	Coarse particles (> 5 µm) ⁴	36 ± 1%	430 ± 22 µg	35 ± 1%	378 ± 17 µg	32 ± 1% 281 ± 14 µg
	Respirable particles (1-5 µm) ⁴	55 ± 1%	658 ± 20 µg	55 ± 1%	591 ± 18 µg	58 ± 1% 506 ± 15 µg
	Total mass of drug (in nebulizer) ⁵	2500 µg		2500 µg		2500 µg
	Total Delivered Mass of Aerosol ⁶	48 ± 2%	1188 ± 41 µg	43 ± 1%	1074 ± 30 µg	35 ± 1% 877 ± 18 µg
	Average Treatment Time (minutes) ⁷	4.3 ± 0.1		3.5 ± 0.1		3.1 ± 0.1

Budesonide Suspension Typical dose 0.50 mg/2 mL or 500 µg of drug substance.	MMAD (µm) ²	4.8 ± 0.1 µ		4.5 ± 0.1 µ		4.3 ± 0.1 µ
	GSD ³	1.9 ± 0.03		2.0 ± 0.05		1.9 ± 0.04
	Extra-fine particles (< 1 µm) ⁴	6 ± 1%	11 ± 1 µg	9 ± 1%	17 ± 2 µg	9 ± 1% 13 ± 2 µg
	Fine particles (< 5 µm) ⁴	49 ± 1%	89 ± 4 µg	51 ± 1%	96 ± 4 µg	54 ± 1% 78 ± 3 µg
	Coarse particles (> 5 µm) ⁴	51 ± 1%	93 ± 7 µg	49 ± 1%	91 ± 4 µg	46 ± 2% 67 ± 4 µg
	Respirable particles (1-5 µm) ⁴	43 ± 1%	77 ± 4 µg	42 ± 1%	79 ± 3 µg	45 ± 1% 65 ± 2 µg
	Total mass of drug (in nebulizer) ⁵	500 µg		500 µg		500 µg
	Total Delivered Mass of Aerosol ⁶	36 ± 2%	182 ± 11 µg	37 ± 1%	187 ± 7 µg	29 ± 1% 145 ± 7 µg
	Average Treatment Time (minutes) ⁷	2.2 ± 0.1		2.0 ± 0.1		1.9 ± 0.1

Cromolyn Sodium Typical dose 20 mg/2 mL or 20000 µg of drug substance.	MMAD (µm) ²	3.7 ± 0.1 µ		3.3 ± 0.1 µ		3.1 ± 0.1 µ
	GSD ³	2.1 ± 0.01		2.1 ± 0.01		2.1 ± 0.02
	Extra-fine particles (< 1 µm) ⁴	12 ± 1%	800 ± 29 µg	14 ± 1%	747 ± 24 µg	14 ± 1% 743 ± 48 µg
	Fine particles (< 5 µm) ⁴	61 ± 1%	4020 ± 120 µg	65 ± 1%	3390 ± 130 µg	68 ± 3% 3560 ± 150 µg
	Coarse particles (> 5 µm) ⁴	39 ± 1%	2610 ± 170 µg	35 ± 1%	1860 ± 130 µg	32 ± 3% 1700 ± 170 µg
	Respirable particles (1-5 µm) ⁴	49 ± 1%	3220 ± 110 µg	50 ± 1%	2640 ± 110 µg	54 ± 2% 2820 ± 110 µg
	Total mass of drug (in nebulizer) ⁵	20000 µg		20000 µg		20000 µg
	Total Delivered Mass of Aerosol ⁶	33 ± 1%	6640 ± 270 µg	26± 1%	5250 ± 220 µg	26 ± 1% 5270 ± 200 µg
	Average Treatment Time (minutes) ⁷	2.3 ± 0.1		1.9 ± 0.1		1.8 ± 0.1

¹ Performance data using Cascade Impactor (NGI) per USP 34 <1601> chilled to 5 °C with an extraction flow of 15 L/min for three common aerosolized medications. Treatment time was defined as onset of audible sputtering plus one minute. Values represent the Mean and 95% confidence interval of the mean.

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⁴ Values are the mass of drug substance delivered in the specified size range of aerosol particles and the mass represented as a percentage of the Total Delivered Mass of Aerosol.

⁵ Total mass of drug in nebulizer is the mass of drug substance placed in the nebulizer cup for typical dose listed.

⁶ The total delivered mass of aerosol is the mass of drug substance recovered from the cascade impactor after treatment time. The mass is stated in micrograms of drug substance and as a percentage of the total mass of drug substance placed in the nebulizer.

⁷ Treatment time, defined as onset of sputtering plus one minute, to deliver the aerosol masses starting with the initial mass of drug listed.

Specifications
- Recommended Driving Gas: Oxygen or Compressed Air - Recommended Driving Gas Flow: 8 L/min - Maximum Recommended Fill Volume: 10 mL - Maximum A-weighted Sound Pressure Level = 52.0 dB - Empty Mass of Nebulizer = 0.019 kg

- (A)** Flextube
- (B)** Mouthpiece
- (C)** Nebulizer Cap
- (D)** Entrainment Port (Do not block)
- (E)** One-Piece Jet
- (F)** Nebulizer Bottle
- (G)** Jet Stem
- (H)** Crush Resistant Tubing

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