AirLife™

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en **AirLife Misty Finity Medium Volume High Flow Continuous Nebulizer**

With 2.1 m (7 foot) Crush Resistant Oxygen Tubing and Adult Aerosol Mask

Indications for Use

This device is intended to be used to aerosolize liquid medication into gases that are delivered to the patient for breathing. It's use is indicated whenever a physician or healthcare professional administers or prescribes medical aerosol products to a patient using the AirLife Misty Finity Large Volume Continuous Nebulizer. This product is a single patient use, non-sterile prescriptive device and is designed to be used in either a hospital, nursing homes, extended care facilities or outpatient clinics.

Directions for Use

1. Pour prescribed medication and saline (diluent) into nebulizer through the aerosol port atop the nebulizer.
Note: Nebulizer cannot be disassembled.
2. Using a push-twist motion, attach the blue rigid connector on the tubing to the side port of the nebulizer. Connect the other end of the tubing to the flow meter of the compressed air or oxygen source.
3. Insert nebulizer into aerosol mask.
4. Place mask over patient's face, pinch nose clip to tighten around nose, and adjust strap above the ears until comfortable for patient.
5. Adjust flow meter to achieve flow rate of 8 l/min. Output is 20 ml/hr.
Note: The listed flow rates and outputs are only intended as a general guideline (± 20%). This is not a substitute for the clinical judgement of qualified personnel.
6. Adjust oxygen concentration per physician order if using oxygen blender.
7. Check liquid consumption at 60 minutes and 120 minutes to determine output of setup, and adjust accordingly per desired output.
8. After use, discard remaining solution and rinse with sterile water.
9. Replace Misty Finity Nebulizer every 48 hours.

Note: The output of the nebulizer can be affected by several factors: It is important for users to follow the “user correctable” approach when utilizing continuous nebulization therapy.

Using a solution, suspension, or emulsion different from that recommended by the manufacturer can alter the particle size distribution curve, the mass median aerodynamic diameter (MMAD), aerosol output, and/or aerosol output rate, which can then be different from those disclosed by the manufacturer

The disclosures listed for nebulizer performance are based upon testing that utilizes adult ventilator patterns and are likely to be different from those stated for the pediatric or infant populations.

Note: See dosing guidelines.

Specifications	
Recommended Driving Gas: Oxygen or Compressed Air	Maximum A-Weighted Sound Pressure Level: 56.5 dB
Recommended Driving Gas Flow: 8 L/min	Empty Mass of Nebulizer = 0.025kg
Maximum Recommended Fill Volume: 30mL	

USA	8 L/min Flow, Regulated at 50 psi**	
	MMAD	1.74 µm
	GSD	2.76
	%<4.7 m	80%

** Measured utilizing an Anderson eight-stage cascade impactor (ACI) with USP inlet and assayed with a spectrophotometer. Particle size characterization was performed using albuterol/saline solution.

Dosing Guidelines

Use the guidelines below (± 20%) for **one hour** of nebulization. These guidelines are no substitute for the good clinical judgement of the user.

First pick your treatment type, then add the specified medication & diluent to the nebulizer reservoir to achieve a total volume of 20 mL. Upon completion, set the flow rate to 8 liters per minute and begin the treatment. The treatment will take approximately 1 hour.

If using Unit Dose Albuterol, use the following:	
Low dose Albuterol Unit dose only (7.5 mg Albuterol/Tx)	3 unit doses Albuterol + 13 ml NaCl 1 hr TX
High dose Albuterol Unit dose only (15 mg Albuterol/Tx)	6 unit doses Albuterol + 3 ml NaCl 1 hr TX

If using multidose Albuterol, use the following:	
High dose bulk Albuterol (7.5 mg Albuterol/Tx)	1½ mL bulk Albuterol + 18 ml NaCl 1 hr TX
High dose bulk Albuterol (15 mg Albuterol/Tx)	3 mL bulk Albuterol + 18 ml NaCl 1 hr TX

Risks and Precautionary Measures related to Phthalates: This instruction pertains to the phthalate symbol marked on the device or its packaging. If this device is used for the treatment of children, or treatment of pregnant or nursing women; please note that the following types of procedures may increase the risk of exposure to phthalates: Exchange transfusion in neonates, total parenteral nutrition in neonates, multiple procedures in sick neonates, haemodialysis in peripuberal males, male foetus and male infant of pregnant women, and lactating women; and massive blood infusion into trauma patients. Although these procedures have the potential for increased risk of exposure, conclusive evidence of human health risks has not been established. As a precautionary measure, to reduce the potential for unnecessary exposures to phthalates, the product must be used in accordance with the instructions for use, and practitioners should refrain from using this product beyond the period of time the product is medically necessary or needed.

Storage

Recommended normal storage temperature 0°C/+32°F to +50°C/+122°F.

Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

en **Misty Finity Medium Volume High Flow Continuous Nebulizer**

Recommended Nebulizer Set Flow Rates 8 L/min		
Drug Description	Aerosol Characteristic ¹	8 L/min
Albuterol SulfateDose 8 mg/8 mL (0.1% by volume) or 8000 µg of drug substance	MMAD [µm] ²	3.4 ± 0.2
	GSD ³	2.6 ± 0.1
	% Particles (< 2µm) ⁴	28.9 ± 1.5
	% Particles (2 to 5µm) ⁴	33.0 ± 1.3
	% Particles (> 5µm) ⁴	38.1 ± 2.3
	% Respirable Particles (< 5µm) ⁴	61.9 ± 2.3
	Aerosol Output Rate [µg/min]	133 ± 4
	Total Delivered Mass of Aerosol [µg] ⁵	2,002 ± 17
	% Fill Volume Emitted per minute	2.0 ± 0.1
	% Fill Volume Emitted ⁶	73.9 ± 0.0
	Residual Volume at Sputter + 1 minute [mL] ⁷	2.1 ± 0.1
	Driving Gas Pressure [psig]	38

Symbols Glossary			
REF	Catalogue number	Rx ONLY	U.S. Federal law restricts this device to sales by or on the order of a physician
LOT	Lot Number	Clean, Ready to Use	Clean, Ready to Use
	Consult instructions for use	QTY	Quantity
	Product and packaging are not made with natural rubber latex		Expiration Date
	Manufacturer	 Single Use	Do not re-use
PHT DEHP	Contains Phthalates	 STORAGE TEMP	Temperature limit

-  Aerosol Mask
-  Flow Meter
-  Oxygen Blender (Optional)
-  Oxygen Supply Tubing
-  Nebulizer

¹ Performance data using Cascade Impactor (NGI) with an extraction flow of 15 L/min and a representative adult breathing patient for 0.1% Albuterol Sulfate solution measured according to BS ISO 27247:2013. Treatment time was defined as onset of audible sputtering plus one minute. Values represent the Mean and 95% confidence interval of the mean.

² MMAD is Mass Median Aerodynamic Diameter.

³ GSD is the Geometric Standard Deviation.

⁴ Values are the mass of drug substance represented as a percentage of the total delivered mass of aerosol.

⁵ The total delivered mass of aerosol is the mass or aerosol emitted by the nebulizing system using a representative adult breathing patient for the given fill volume (i.e. 8mL) that is emitted by the nebulizer for a treatment time of sputter plus one minute

⁶ The percentage of fill volume emitted is the total delivered mass of aerosol expressed as a percentage of the fill volume (i.e. 8mL) that is emitted by the nebulizer for a treatment time of sputter plus one minute.

⁷ Residual Volume is the volume of liquid remaining in the nebulizer system for a treatment time of sputter plus one minute.

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AirLife Australia Holdings Pty Ltd