REF 001227

AirLife™

Adult Tracheostomy Mask (Large)

For Single Patient Use

REF	Catalogue number	R _X ONLY	U.S. Federal law restricts this device to sales by or on the order of a physician
LOT	Lot Number	Clean, Ready to Use	Clean, Ready to Use
QTY	Quantity	Manufacturer	Manufacturer
LATEX	Product and packaging are not made with natural rubber latex	Date of manufacture	Date of manufacture
PHT DEHP	Contains Phthalates	Consult instructions for use	Consult instructions for use

EN

English

Intended Use

Aerosol Masks are intended to provide the patient with the means of receiving aerosol. Drainage bags are devices intended to trap and drain water that collects in tubing during respiratory therapy, thereby preventing an increase in breathing resistance.

Directions for Use

1. Prior to use, confirm that there is adequate clearance between the mask and the tracheostomy tube.
2. A large bore tube 22 mm ID should be connected to outlet arm of ball swivel on top of mask. The other end of tubing is connected to aerosol generator.
3. Aerosol mist to be delivered as specified by physician and in accordance with manufacturer's instructions.
4. Place mask onto neck region covering tracheostomized area.
5. Extend strap behind patient's neck and slide plastic connector on strap into plastic retainer on mask until it snaps into locking position.
6. Adjust strap to comfort of patient.
7. To remove mask, depress plastic connector on strap and pull out.

Note

If suctioning is required, remove mask.

AirLife,
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com

AirLife®

Contains Phthalates (DEHP). Risks and Precautionary

Measures related to Phthalates: This instruction pertains to the phthalate symbol marked on the device or its packaging. If this device is used for the treatment of children, or treatment of pregnant or nursing women; please note that the following types of procedures may increase the risk of exposure to phthalates: Exchange transfusion in neonates, total parenteral nutrition in neonates, multiple procedures in sick neonates, haemodialysis in peripubertal males, male foetus and male infant of pregnant women, and lactating women; and massive blood infusion into trauma patients. Although these procedures have the potential for increased risk of exposure, conclusive evidence of human health risks has not been established. As a precautionary measure, to reduce the potential for unnecessary exposures to phthalates, the product must be used in accordance with the instructions for use, and practitioners should refrain from using this product beyond the period of time the product is medically necessary or needed.

Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

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