
AirLife™

- EN Adult Mouth-to-Mask Resuscitator**
Oxygen Port, Filter, 6 inch (15 cm) Tubing, Mouthpiece
- FR Insufflateur bouche-masque pour adulte**
Orifice d'oxygène, filtre, tuyau de 15 cm (6 po), embout buccal
- DE Mund-zu-maske-Beatmungsgerät für Erwachsene**
Sauerstoffanschluss, Filter, 15 cm (6") Schlauch, Mundstück
- IT Rianimatore con maschera per bocca per adulti**
Porta di ossigeno, filtro, tubo di 15 cm (6 poll.) e boccaglio
- ES Resucitador de boca a máscara para adultos**
Puerto para oxígeno, filtro, tubo de 15 cm (6 pulgadas), boquilla
- PT Ressuscitador Boca-a-Máscara para Adultos**
Porta de Oxigénio, Filtro, Tubulação de 15 cm (6 polegadas), Bocal
- NL Mond-maskerresuscitator voor volwassenen**
Zuurstofpoort, filter, slang van 15 cm (6 inch), mondstuk
- SV Mun-mot-mask-återupplivningsanordning för vuxna**
Syrgasport, filter, 15 cm (6 tum) slangar, munstycke
- DA Mund-til-maske-genoplivningsudstyr til voksne**
Iltindtag, filter, 15 cm (6 tommer) slange, mundstykke
- NO Munn-til-maske gjenopplivningsenhet for voksne**
Oksygenport, filter, 15 cm (6 tommer) slange, munnstykke
- FI Aikuisten elvytyssuoja**
Happiliitin, suodatin, 15 cm:n (6 tuumaa) letku, suukappale
- PL Maska resuscytacyjna dla dorosłych**
Port tlenu, filtr, przewody o długości 15 cm (6 cali), ustnik

LOT

AirLife®

Rx ONLY

Clean, Ready to Use

AirLife
2710 Northridge Dr. NW, Suite A
Grand Rapids, MI 49544 USA
www.myAirLife.com

Made in Mexico

2K8012.INIFU

Rev.1 2025/03

EN Adult Mouth-to-Mask Resuscitator

Oxygen Port, Filter, 6 inch (15 cm) Tubing, Mouthpiece

Intended Use

The resuscitator is a single use resuscitator intended for pulmonary resuscitation

Cautions

- This product must be used by personnel thoroughly trained in the techniques of pulmonary resuscitation.
- For single use only. Re-use may degrade the performance of the product or contribute to cross contamination.

Warning

- Permanent distortion will occur if the resuscitator bag is stored in a compressed state. This may result in reduced ventilation efficiency.

Precautions

Do

- When not using supplemental oxygen, ensure cap covers oxygen port "C".
- Clear patient's airway before using resuscitator.
- Always check for proper function of resuscitator:
 - Verify proper valve action.
 - Verify that the valve is free of obstruction.
- Verify patient is being ventilated by observing alternate rise and fall of the patient's chest and color of lips and face during resuscitation.

Do Not

- Do not use in contaminated atmosphere (e.g., poisonous gases, smoke, etc.).
- If oxygen is used, do not smoke or use in presence of sparking equipment or open flame.
- Do not autoclave, gas or chemically sterilize the Resuscitator.

Instructions For Use

Set up—Resuscitator

- Remove resuscitator and mask from package and firmly seat the resuscitator housing "B" into the mask port "A".
- Prior to using the resuscitator, visually verify proper valve action.
- If resuscitating with supplemental oxygen, attach oxygen port "C" to proper oxygen source.
- Set oxygen flow not to exceed recommendation of a physician.

Operation—Resuscitator

- Clear patient's airway, if obstructed.
- Tilt patient's head back and pull chin up.
- After establishing this position, place mask firmly over the nose and mouth and hold in place with thumb and index finger, grip chin with other fingers.
- Using mouthpiece "D", inflate lungs using currently accepted mouth-to-mask breathing technique.
- Verify that patient's chest rises and falls during resuscitation. If movement is absent during resuscitation, check patient's airway.
- Time resuscitation with any spontaneous breathing to prevent blockage of exhalation.
- If foreign matter, such as blood or vomitus appears in the mask, remove the resuscitator from patient and clear patient's airway. Clear any obstructions from the resuscitator by blowing into the mouthpiece "D" and

shaking any remaining matter out of the mask "A" and valve housing "B".

8. Discard resuscitator after use.

- Dispose of all materials in accordance with local, state, and federal regulations. Decontaminate and dispose of all potentially biohazardous material.

Risks and Precautionary Measures related to

Phthalates: This instruction pertains to the phthalate symbol marked on the device or its packaging. If this device is used for the treatment of children, or treatment of pregnant or nursing women; please note that the following types of procedures may increase the risk of exposure to phthalates: Exchange transfusion in neonates, total parenteral nutrition in neonates, multiple procedures in sick neonates, haemodialysis in peripuberal males, male foetus and male infant of pregnant women, and lactating women; and massive blood infusion into trauma patients. Although these procedures have the potential for increased risk of exposure, conclusive evidence of human health risks has not been established. As a precautionary measure, to reduce the potential for unnecessary exposures to phthalates, the product must be used in accordance with the instructions for use, and practitioners should refrain from using this product beyond the period of time the product is medically necessary or needed.

Symbols Glossary			
	Catalogue number		Manufacturer
	Lot Number		Quantity
	Clean, Ready to Use		Date of Manufacture
	Consult instructions for use		U.S. Federal law restricts this device to sales by or on the order of a physician
	Not made with natural rubber latex		Caution
	Contains Phthalates		For Single use only

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