AirLife™ Closed Suction Catheter System for use with Neo-Verso™ Airway Access Adapters

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







Product and packaging are not made with natural rubber latex or DEHP.

WARNINGS:

1. Fully withdraw catheter before cutting endotracheal tube, otherwise the catheter may also be cut.
2. If using the Neo-Verso 'Y' adapter, be certain to use the correct size tip; otherwise the 'Y' adapter may loosen from the endotracheal tube and cause loss of ventilating pressure.

CAUTIONS:

1. Select the appropriate French-sized AirLife catheter. Suggested catheter sizes are:

ET Tube I.D	2.0 mm	2.5 mm	3.0 mm	3.5 mm	4.0 mm
Maximum Size Catheter	5 French	5 French	6 French	6 French	8 French

2. Inspect product packaging before opening. Do not use product if packaging has been compromised.
3. For single use only. Re-use may degrade the performance of the product or contribute to cross contamination.
4. Do not sterilize or re-sterilize.
5. After suctioning the patient, always pull the catheter back until the black bar is visible within the sleeve and the catheter is fully retracted in order to avoid increased airway resistance.
6. The Neo-Verso Adapters should be replaced with circuit change-outs, following 7 days of use, or earlier if soiled.
7. The AirLife catheter is approved for use up to 72 hours. Depending on frequency of use and amount of secretions, change the catheter every 24-72 hours or as needed to maximize product performance.
8. Use appropriate suction technique and regulated vacuum levels.
9. During and immediately after a suctioning procedure, monitor the patient for any signs of suction intolerance. Adjustments to the ventilator or oxygen settings may need to be made.
10. Always reinsert the adapter's access port plug when no instrument or device is attached to the adapter.

Product set up AirLife Catheter Systems

1. Select the appropriate sized AirLife closed suction catheter and attach thumb control valve of catheter to suction tubing.
2. Depress and hold thumb valve and simultaneously adjust vacuum regulator to desired level (Fig 1).
3. Release thumb control valve.

Product set up for the Neo-Verso Adapter

1. Place the Neo-Verso between the patient endotracheal tube and the ventilator circuit (Fig 2).
2. Stabilize the Neo-Verso, remove port plug and attach the AirLife closed suction catheter by pushing catheter connector fitting snugly into the Neo-Verso access port (Fig 2).

Product set up for the Neo-Verso 'Y' Adapter

1. Select the matching 'Y' adapter tip to the endotracheal tube size and firmly press it on to the Neo-Verso 'Y' body. Assure the connection (Fig 3.). Once the tip has been attached to the Neo-Verso 'Y', it cannot be removed.
2. Remove the original fitting from the endotracheal tube and firmly attach Neo-Verso 'Y' adapter directly to endotracheal tube and the ventilator circuit (Fig 4).
3. Stabilize the adapter, remove access port plug and attach the AirLife closed suction catheter by pushing catheter connector fitting snugly into the Neo-Verso 'Y' access port (Fig 2).

Suggested Suction Procedure:

1. Stabilize the Neo-Verso with one hand and push the catheter into the endotracheal tube with the opposite hand.
2. Advance catheter to desired depth following instructions below (Fig 5).
3. Depress and hold thumb control valve to apply suction, then gently withdraw catheter. Stop withdrawal when black bar marking is visible inside the sleeve (Fig 6).
4. Release thumb control valve.
5. Repeat steps 1-4 above as necessary.
6. Note: For open suctioning, follow the same above procedure using a sterile open suction catheter while being sure to completely remove the catheter from the adapter after use. Be sure to attach the Neo-Verso adapter's access port plug when finished.

7. Maintain superior suctioning performance by utilizing the valved irrigation port on the AirLife catheter to flush the catheter with saline after each use (see instructions below).

Directions for Bronchoscopy Use (seals sold separately):

1. Remove the access port plug (or closed suction catheter) from the Neo-Verso access adapter port.
2. Select the appropriate sized AirLife Adapter Seal based on the bronchoscope diameter.
3. Firmly press and rotate the AirLife Adapter Seal over the lip of the access port assuring the lip is covered by the seal all around.
4. Remove and discard the adapter seal and reinsert the adapter's access port plug when the procedure is completed.

Using Depth Markings:

1. Align any printed depth number on the catheter with the similar number printed on the endotracheal tube. When these numbers match, the catheter tip has reached the end of the endotracheal tube.
2. Alternatively, note the length on the ET tube where it meets the tip of the plastic adapter.

Thumb Control Valve Operation:

1. The thumb control valve can be locked to prevent inadvertent or accidental suction. Lift the Thumb Control on the valve and rotate it 180 degrees so that its arrow aligns with the valve body lock symbol (Fig 7).
2. To unlock, repeat this action and align the arrow with the valve body unlocked symbol.

Day Sticker Usage:

1. Apply the appropriate day sticker to the thumb control valve or its cap tether. Example: If the AirLife closed suction catheter is to be changed at 48 hours and opened on Monday, apply the "Change on Wednesday" sticker.

Flushing Closed Suction Catheters:

1. Insure that the catheter is pulled back completely. You should be able to visualize the black bar marking on the catheter through the catheter sleeve.
2. Open the white cap on the flush port of the AirLife Closed Suction Catheter. A built-in check valve maintains a closed ventilator circuit.
3. Attach a single-use vial containing water or saline to the flush port.
4. Depress and hold the thumb control valve on the suction catheter.
5. If fluid is not automatically suctioned through the catheter, slowly squeeze the solution into the Neo-Verso adapter while applying suction pressure.
6. After releasing the thumb control valve, remove the vial from the flush port.
7. Close the cap on the flush port.

Patient Lavage / Saline Instillation:

1. Attach a syringe or single-use vial containing normal saline to the flush port.
2. To avoid introducing saline into the circuit or HME, rotate the Neo-Verso Adapter so that the circuit is higher than the endotracheal tube.
3. Push/Squeeze the saline into the flush port. Saline will flow into the endotracheal tube.
4. Advance the suction catheter to the desired depth and suction as previously described.

Changing Closed Suction Catheters Without Disconnecting the Circuit:

1. Stabilize Neo-Verso adapter with one hand and rotate the suction catheter connector fitting with the opposite hand to loosen (Fig 8).
2. Remove suction catheter.
3. Check that new suction catheter is the correct size for the patient's artificial airway.
4. Attach new AirLife closed suction catheter to the access port of the Neo-Verso adapter.

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AirLife Australia Holdings Pty Ltd