AirLife™ Closed Suction Catheter System for use with Verso™ Airway Access Adapters

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






USA Rx Only • Clean, Ready to Use

WARNINGS:

1. Fully withdraw catheter before cutting endotracheal tube, otherwise the catheter may also be cut.

CAUTIONS:

1. Select the appropriate French-sized AirLife catheter. Clinical experts typically select catheters that will occupy no more than ½ of the internal diameter of the artificial airway.
2. Inspect product packaging before opening. Do not use product if packaging has been compromised.
3. For single use only. Re-use may degrade the performance of the product or contribute to cross contamination.
4. Do not sterilize or re-sterilize.
5. After suctioning the patient, always pull the catheter back until the black bar is visible within the sleeve and the catheter is fully retracted in order to avoid increased airway resistance.
6. The Verso Adapter should be replaced with circuit change-outs, following 7 days of use, or earlier if soiled.
7. The AirLife catheter is approved for use up to 72 hours. Depending on frequency of use and amount of secretions, change the catheter every 24-72 hours or as needed to maximize product performance.
8. If introducing fluid into the Verso Adapter, insure that fluid does not come in contact with HMEs (Heat & Moisture Exchangers) as this may increase gas flow resistance.
9. During and immediately after a suctioning procedure, monitor the patient for any signs of suction intolerance. Adjustments to the ventilator or oxygen settings may need to be made.
10. For tracheostomy patients, use tracheostomy-length AirLife closed suction catheters to reduce mucosal damage.
11. Always reinsert the Verso Adapter's access port plug when no instrument or device is attached to the adapter.

Product set up

1. Place the Verso Adapter between the patient endotracheal or tracheostomy tube and the ventilator circuit (Fig 1).
2. Select the appropriate sized AirLife closed suction catheter and attach thumb control valve of catheter to suction tubing.
3. Depress and hold thumb valve and simultaneously adjust vacuum regulator to desired level. Release thumb control valve (Fig 6).
4. Stabilize the Verso Adapter, remove port plug and attach the closed suction catheter by pushing catheter connector fitting snugly into the Verso access port (Fig 1).

Suggested Suction Procedure:

1. Stabilize the Verso Adapter with one hand and push the catheter into the endotracheal tube with the opposite hand (Fig 2).
2. Advance catheter to desired depth.
3. Depress and hold thumb control valve to apply suction, then gently withdraw catheter. Stop withdrawal when black bar marking is visible inside the sleeve (Fig 3).
4. Release thumb control valve.
5. Repeat steps 1-4 above as necessary.
6. Note: For open suctioning, follow the same above procedure using a sterile open suction catheter while being sure to completely remove the catheter from the adapter after use. Be sure to attach the Verso adapter's airway access port plug when finished.
7. Maintain superior suctioning performance by utilizing the valved irrigation port on the Verso Adapter to flush the catheter with saline after each use (see instructions below).

Directions for Bronchoscopy Use (seals sold separately):

1. Remove the access port plug (or closed suction catheter) from the Verso Adapter port.
2. Select the appropriate sized AirLife Adapter Seal based on the bronchoscope diameter.
3. Firmly press and rotate the AirLife Adapter Seal over the lip of the access port assuring the lip is covered by the seal all around. Lubricate bronchoscope.
4. At the completion of the procedure, withdraw the bronchoscope to the inside lip of the adapter seal and flush the chamber with saline while aspirating through the bronchoscope.
5. Remove and discard the adapter seal and reinsert the adapter's access port plug when the procedure is completed.

Directional Tip Procedure:

1. For optimal catheter directional control, tracheostomy or endotracheal tube should be positioned 4-5 cm above carina.
2. The tip of the catheter is bent away from the printed depth numbers on the catheter and indicates the direction the catheter tip will follow.
3. Direct the catheter into the desired lung by keeping the printed numbers on the catheter oriented away from the desired side.

Using Depth Markings:

1. Align any printed depth number on the catheter with the similar number printed on the endotracheal tube. When these numbers match, the catheter tip has reached the end of the endotracheal tube.

Thumb Control Valve Operation:

1. The thumb control valve can be locked to prevent inadvertent or accidental suction. Lift the Thumb Control on the valve and rotate it 180 degrees so that its arrow aligns with the valve body lock symbol (Fig 4).
2. To unlock, repeat this action and align the arrow with the valve body unlocked symbol.

Day Sticker Usage:

1. Apply the appropriate day sticker to the thumb control valve or its cap tether. Example: If the AirLife closed suction catheter is to be changed at 72 hours and opened on Monday, apply the "Change on Thursday" sticker.

Flushing Closed Suction Catheters or Bronchoscopes:

1. Insure that the catheter is pulled back completely. You should be able to visualize the black bar marking on the catheter through the catheter sleeve.
2. Open the white cap on the side flush port of the Verso Adapter. A built-in check valve maintains a closed ventilator circuit.
3. Attach a single-use vial containing water or saline to the flush port.
4. Depress and hold the thumb control valve on the suction catheter or bronchoscope.
5. If the fluid is not automatically suctioned through the catheter, slowly squeeze the solution into the Verso Adapter while applying suction pressure.
6. After releasing the thumb control valve, remove the vial from the flush port.
7. Close the cap on the flush port.

Patient Lavage / Saline Instillation:

1. Attach a syringe or single-use vial containing normal saline to the flush port.
2. To avoid introducing saline into the circuit or HME, rotate the swivel on the Verso Adapter so that the circuit is higher than the endotracheal tube.
3. Push the suction catheter far enough into the Verso Adapter so that it opens the check valve and enters the artificial airway.
4. Push/Squeeze the saline into the flush port. Saline will flow into the endotracheal tube.
5. Advance the suction catheter to the desired depth and suction as previously described.

Changing Closed Suction Catheters Without Disconnecting the Circuit:

1. Stabilize the Verso Adapter with one hand and rotate the suction catheter connector fitting with the opposite hand to loosen (Fig 5).
2. Remove suction catheter.
3. Check that new suction catheter is the correct size for the patient's artificial airway.
4. Attach new AirLife closed suction catheter to the access port of the Verso Adapter.

CSC110T, CSC112T, CSC114T, CSC116T, CSC114D: Sterile Devices (EO), do not use if package is damaged.

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AirLife Australia Holdings Pty Ltd